

SOMATOM go.Up - Exam Protocols *syngo* CT VA40

The dose values stated in the exam protocol list are estimations and provide orientation and overview. For individual exam protocols, the valid dose values are displayed at the CT scanner.

The descriptions below refer to the corresponding column headings of the exam protocol tables on the following pages.

FAST kV IQ level

'FAST kV IQ level' specifies the relative image quality that the system would achieve under specific reference conditions with regards to the patient's size and anatomic shape and standard tube voltage of 120 kV, independent of the specific scanner conditions.

Eff. mAs

'Eff mAs' indicates the values for operation with deactivated CARE Dose4D. The dose values are defined as a typical adult (70 - 80 kg or 155 - 180 lbs) for adult protocols, or as a typical child (5 years, appr. 20 kg or 45 lbs) for child protocols.

CTDIvol [mGy]

Measured in a 16 cm phantom for head modes and in a 32 cm phantom for body scans.

CTDIvol [mGy]; DLP [mGy*cm]

The indicated dose values 'CTDIvol', and 'DLP' refer to the indicated settings for operation with deactivated CARE Dose4D, regardless if CARE Dose4D is activated by default or not.

Shaped Filter

Displays the setting of the filter in collimator, standard or narrow.

Absolute Start Delay [s]

The delay time of all ranges of an auto range always refers to 0.

Iterative Recon (SAFIRE or ADMIRE)

It is an iterative reconstruction and the Strength can be adjusted to the customer's preference.

For each Siemens Healthineers exam protocol [factory], the Decision Trees and the Recon Strategies are listed. The scan strategies can be combined with patient characteristics and categories or patient questions into a decision tree in the myExam Cockpit workspace in the Exam Designer.

For each Siemens Healthineers exam protocol [factory], the Contrast Protocols are listed. The contrast protocols are displayed in the Contrast Protocols workspace in the Exam Designer, where you can set up the injection phases and contrast agent information for CARE Contrast CT examinations.

Protocol Name	Protocol Description	Range Name	Recon Name	Organ Characteristic	Scan Mode	Decision Trees	Acquisition [mm]	CARE Dose4D	FAST kV IQ level	FAST kV optimized for	kV	Eff. mAs	mA	CTDIvol [mGy]	DLP [mGy*cm]	Shaped Filter	Feed [mm]	Pitch	Rot. Time [s]	Scan Time [s]	Absolute Start Delay [s]	Contrast Protocol	Slice Thickness [mm]	Increment [mm]	Recon Strategies			Inline Results	Multi Recons	Image Type	Kernel	Iterative Recon	IR Strength	Recon Type	Recon Phase Step	Recon Phase(s)	Recon Phase Unit				
Brain [factory] (Adult)	>Indications: stroke, hemorrhage, cerebral atrophy, hydrocephalus, inflammation, w/o CM >Scan modes: RoutineSpiralAdultHead >Decision Tree: *Brain contrast medium necessary [factory]	Topogram			ScanplanningTopoAdult		6 x 0.70				130		30	SPR (16 cm) 0.25	6.23	Standard					3.0																				
			Topogram 0.70 Tr20 sag																														Tr20								
		Brain		Brain	RoutineSpiralAdultHead	*Brain contrast medium necessary [factory]	32 x 0.70	Off			130	230		(16 cm) 51.0	716	Standard		0.55	1.00		3.0																				
			Brain 5.00 Hr40 I2																				5.00	5.00		Brain contrast medium necessary [factory],w/o CM; Brain contrast medium necessary [factory],w CM	None	Impression Multi Recons (2)			Hr40	On	2								
			Brain 5.00 Hr60 I2																												Hr60	On	2								
Brain Multi Phase [factory] (Adult)	>Indications: stroke, tumors, cerebral atrophy, hydrocephalus, inflammation,w and w/o CM >Scan modes: RoutineSpiralAdultHead >Contrast Protocol: Brain [factory]	Topogram			ScanplanningTopoAdult		6 x 0.70				130		30	SPR (16 cm) 0.25	6.23	Standard					3.0												Tr20								
			Topogram 0.70 Tr20 sag																																						
		Non Contrast Brain		Brain	RoutineSpiralAdultHead		32 x 0.70	Off			130	230		(16 cm) 51.0	716	Standard		0.55	1.00		3.0																				
			Non Contrast Brain 5.00 Hr40 I2																				5.00	5.00			None	Impression Multi Recons (2)			Hr40	On	2								
			Non Contrast Brain 5.00 Hr60 I2																												Hr60	On	2								
		Venous Phase Brain		Brain	RoutineSpiralAdultHead		32 x 0.70	Off			130	251		(16 cm) 55.7	782	Standard		0.55	1.00		180.0	Brain [factory]																			
			Venous Phase Brain 5.00 Hr40 I2																				5.00	5.00			None	None			Hr40	On	2								
Brain Perfusion [factory] (Adult)	>Indications: stroke, hemorrhage, ischemia, malperfusion, aneurysm, w and w/o CM >Possible application: CT Neuro Perfusion >Scan modes: RoutineSpiralAdultHead, DynamicSequenceAdultHead, PremonitoringSequenceAdultI2d, MonitoringSequenceAdultI2d, RoutineSpiralAdultAngioNeck >Contrast Protocol: Brain Perfusion [factory], CTA [factory]	Topogram			ScanplanningTopoAdult		6 x 0.70				130		30	SPR (16 cm) 0.25	6.23	Standard					3.0											Tr20									
			Topogram 0.70 Tr20 sag																																						
		Non Contrast Brain		Brain	RoutineSpiralAdultHead		32 x 0.70	Off			130	230		(16 cm) 51.0	716	Standard		0.55	1.00		3.0																				
			Non Contrast Brain 5.00 Hr40 I2																				5.00	5.00			None	Impression Multi Recons (2)			Hr40	On	2								
			Non Contrast Brain 5.00 Hr60 I2																												Hr60	On	2								
			Non Contrast Brain 5.00 Hr40 I2 ASPECT																				5.00	5.00			ASPECT Score	None			Hr40	On	2								
			Non Contrast Brain 4.00 Hr40 I3 BrainHem																				4.00	4.00			Brain Hemorrhage	None				On	3								
		Perfusion Brain		Brain	DynamicSequenceAdultHead		32 x 0.70	Off			80	96		(16 cm) 171	382	Standard	0.0		0.80	Full 0.80	3.0	Brain Perfusion [factory]		5.00	3.00							Hr36	Off	2							
			Perfusion Brain 5.00 Hr36 3.80s - 48.60s																																						
		Topogram			ScanplanningTopoAdult		6 x 0.70				110	15		SPR (32 cm) 0.03	1.60	Standard					7.0												Tr20								
			Topogram 0.70 Tr20 cor																																						
		Premonitoring		Vascular Body			1 x 10.00	On	25	Soft Tissue with Contrast	110	22		(32 cm) 1.41	1.41	Standard	0.0		0.80	Full 0.80	3.0																				
			Premonitoring 10.00 Br36																				10.00					None			Br36										
		Monitoring		Vascular Body			1 x 10.00	On	25	Soft Tissue with Contrast	110	22		(32 cm) 42.3	42.3	Standard	0.0		0.80	Full 0.80	10.0	CTA [factory]																			
			Monitoring 10.00 Br36																				10.00					None			Br36										
		CTA Brain and Neck		Neck	RoutineSpiralAdultAngioNeck		32 x 0.70	On	145	Auto (Vascular)	110	111		(32 cm) 7.78	341	Standard		1.20	0.80			CTA [factory]																			
			CTA Brain and Neck 3.00 Hv36 I2 sag CTA Brain and Neck 3.00 Hv36 I2 cor																				3.00	3.00			None	Orientation Multi Recons (3)	MIP Thin	Hv36	On	2									
			CTA Brain and Neck 3.00 Hv36 I2 ax																																						

Protocol Name	Protocol Description	Range Name	Recon Name	Organ Characteristic	Scan Mode	Decision Trees	Acquisition [mm]	CARE Dose4D	FAST KV IQ level	FAST KV optimized for	kV	Eff. mAs	CTDIvol [mGy]	DLP [mGy*cm]	Shaped Filter	Feed [mm]	Pitch	Rot. Time [s]	Scan Time [s]	Absolute Start Delay [s]	Contrast Protocol	Slice Thickness [mm]	Increment [mm]	Recon Strategies			Inline Results	Multi Recons	Image Type	Kernel	Iterative Recon	IR Strength	Recon Type	Recon Phase Step	Recon Phase(s)	Recon Phase Unit		
Brain Perfusion Prolonged [factory] (Adult)	>Indications: stroke, hemorrhage, ischemia, malperfusion, aneurysm, w and wo CM >Possible application: CT Neuro Perfusion >Scan modes: RoutineSpiralAdultHead, DynamicSequenceAdultHead, PremonitoringSequenceAdult2d, MonitoringSequenceAdult2d, RoutineSpiralAdultAngioNeck >Contrast Protocol: Brain Perfusion [factory], CTA [factory]	Topogram			ScanplanningTopoAdult		6 x 0.70				130	30	SPR (16 cm) 0.25	6.23	Standard						3.0																	
		Topogram 0.70 Tr20 sag																													Tr20							
		Non Contrast Brain	Brain	RoutineSpiralAdultHead			32 x 0.70	Off			130	230	(16 cm) 51.0	716	Standard		0.55	1.00			3.0																	
		Non Contrast Brain 5.00 Hr40 I2																				5.00	5.00		None	Impression Multi Recons (2)		Hr40	On	2								
		Non Contrast Brain 5.00 Hr60 I2																									Hr60	On	2									
		Non Contrast Brain 5.00 Hr40 I2 ASPECT																				5.00	5.00		ASPECT Score	None		Hr40	On	2								
		Non Contrast Brain 4.00 Hr40 I3 BrainHem																				4.00	4.00		Brain Hemorrhage	None			On	3								
		Perfusion Brain	Brain	DynamicSequenceAdultHead			32 x 0.70	Off			80	96	(16 cm) 171	382	Standard	0.0		0.80	Full 0.80		3.0	Brain Perfusion [factory]																
		Perfusion Brain 5.00 Hr36 3.80s - 61.40s																				5.00	3.00							Hr36	Off	2						
		Topogram			ScanplanningTopoAdult		6 x 0.70					110	15	SPR (32 cm) 0.03	1.60	Standard						7.0																
		Topogram 0.70 Tr20 cor																													Tr20							
		Premonitoring	Vascular Body				1 x 10.00	On	25	Soft Tissue with Contrast	110	22	(32 cm) 1.41	1.41	Standard	0.0		0.80	Full 0.80		3.0																	
		Premonitoring 10.00 Br36																					10.00						None		Br36							
		Monitoring	Vascular Body				1 x 10.00	On	25	Soft Tissue with Contrast	110	22	(32 cm) 42.3	42.3	Standard	0.0		0.80	Full 0.80		10.0	CTA [factory]																
		Monitoring 10.00 Br36																					10.00						None		Br36							
		CTA Brain and Neck	Neck	RoutineSpiralAdultAngioNeck			32 x 0.70	On	145	Auto (Vascular)	110	111	(32 cm) 7.78	341	Standard		1.20	0.80				CTA [factory]																
		CTA Brain and Neck 3.00 Hv36 I2 sag																					3.00	3.00		None	Orientation Multi Recons (3)	MIP Thin	Hv36	On	2							
		CTA Brain and Neck 3.00 Hv36 I2 cor																																				
		CTA Brain and Neck 3.00 Hv36 I2 ax																																				
		Inner Ear [factory] (Adult)	>Indications: inflammation, tumorous processes of pyramids, acoustic neuroma, pathologies of the mastoid process, cerebellopontine angle tumors, cochlea implants, trauma, wo CM >Scan modes: RoutineSpiralAdultHeadSn	Topogram			ScanplanningTopoAdult		6 x 0.70				130	30	SPR (16 cm) 0.25	6.23	Standard						3.0															
Topogram 0.70 Tr20 sag																															Tr20							
Inner Ear	Brain			RoutineSpiralAdultHeadSn			32 x 0.70	On	110	Auto (Bone/Calcium)	Sn 130	316	(16 cm) 16.9	119	Standard		0.55	1.00			3.0																	
Inner Ear 1.50 Hr60 I3 ax																						1.50	1.50		None	Impression Multi Recons (2)	MPR	Hr60	On	3								
Inner Ear 1.50 Hr40 I3 ax																											Hr40	On	3									
Inner Ear 0.80 Hr60 I3 cor R																							0.80	0.80		Parallel Ranges	None	MPR	Hr60	On	3							
Inner Ear 0.80 Hr60 I3 cor L																					0.80	0.80		Parallel Ranges	None	MPR	Hr60	On	3									

Protocol Name	Protocol Description	Range Name	Recon Name	Organ Characteristic	Scan Mode	Decision Trees	Acquisition [mm]	CARE Dose4D	FAST KV IQ level	FAST KV optimized for	KV	Eff. mAs	mA	CTDIvol [mGy]	DLP [mGy*cm]	Shaped Filter	Feed [mm]	Pitch	Rot. Time [s]	Scan Time [s]	Absolute Start Delay [s]	Contrast Protocol	Slice Thickness [mm]	Increment [mm]	Recon Strategies			Inline Results	Multi Recons	Image Type	Kernel	Iterative Recon	IR Strength	Recon Type	Recon Phase Step	Recon Phase(s)	Recon Phase Unit
Sinus [factory] (Adult)	>Indications: sinusitis, mucocoele, pneumatization, polyposis, tumor, corrections, trauma, wo CM >Scan modes: RoutineSpiralAdultHeadsSnXcare >Decision Tree: Sinus sinusitis or soft tissue [factory]	Topogram			ScanplanningTopoAdult		6 x 0.70				130		30	SPR (16 cm) 0.25	6.23	Standard					3.0																
		Topogram 0.70 Tr20 sag																														Tr20					
		Sinus		Brain	RoutineSpiralAdultHeadsSnXcare	Sinus sinusitis or soft tissue [factory]	32 x 0.70	On	50	Bone/Calcium	Sn 110	236	(16 cm) 6.43	64.6	Standard		0.55	1.00		3.0																	
		Sinus 3.00 Hr60 I3 ax																					3.00	3.00		Sinus sinusitis or soft tissue [factory].low dose	None	Impression Multi Recons (2)	MPR	Hr60	On	3					
		Sinus 3.00 Hr40 I3 ax																								Sinus sinusitis or soft tissue [factory].low dose	Parallel Ranges	None	MPR	Hr60	On	3					
		Sinus 3.00 Hr60 I3 cor																					3.00	3.00		Sinus sinusitis or soft tissue [factory].low dose	None	Impression Multi Recons (2)	MPR	Hr40	On	3					
		Sinus 3.00 Hr40 I3 ax																								Sinus sinusitis or soft tissue [factory].standard dose	None	Impression Multi Recons (2)	MPR	Hr60	On	3					
		Sinus 3.00 Hr60 I3 ax																											Hr60	On	3						
		Sinus 3.00 Hr40 I3 cor																					3.00	3.00		Sinus sinusitis or soft tissue [factory].standard dose	Parallel Ranges	None	MPR	Hr40	On	3					
Orbit [factory] (Adult)	>Indications: fracture, tumor, trauma, wo CM >Scan modes: RoutineSpiralAdultHeadsSnXcare	Topogram			ScanplanningTopoAdult		6 x 0.70				130		30	SPR (16 cm) 0.25	6.23	Standard					3.0																
		Topogram 0.70 Tr20 sag																														Tr20					
		Orbit		Brain	RoutineSpiralAdultHeadsSnXcare		32 x 0.70	On	110	Auto (Bone/Calcium)	Sn 110	475	(16 cm) 12.9	91.1	Standard		0.55	1.00		3.0																	
		Orbit 2.00 Hr60 I3																					2.00	2.00			None	Impression Multi Recons (2)		Hr60	On	3					
		Orbit 2.00 Hr40 I3																											Hr40	On	3						
		Orbit 2.00 Hr60 I3 cor R																					2.00	2.00			Parallel Ranges	None	MPR	Hr60	On	3					
		Orbit 2.00 Hr60 I3 cor L																					2.00	2.00			Parallel Ranges	None	MPR	Hr60	On	3					
Orbit 2.00 Hr40 I3 cor Radial																											Parallel Ranges	None	VRT	Hr40	On	3					
Dental [factory] (Adult)	>Indications: display and measurement of the bone structures of the upper and lower jaw as the basis for planning in oral surgery, wo CM >Possible application: CT Dental >Scan modes: RoutineSpiralAdultHeadsSnXcare	Topogram			ScanplanningTopoAdult		6 x 0.70				130		30	SPR (16 cm) 0.25	6.23	Standard					3.0																
		Topogram 0.70 Tr20 sag																														Tr20					
		Dental		Brain	RoutineSpiralAdultHeadsSnXcare		32 x 0.70	On	65	Bone/Calcium	130	44	(16 cm) 9.76	68.8	Standard		0.55	1.00		3.0																	
		Dental 1.00 Hr60 I3																					1.00	0.70			None	None		Hr60	On	3					
	Dental 5.00 Hr60 I3 ax																						5.00	5.00			None	None	MPR	Hr60	On	3					
Neck [factory] (Adult)	>Indications: masses, tumor, metastases, lymphoma, lymph nodes in the larynx or pharynx region, w/wo CM >Scan modes: RoutineSpiralAdultHeadsSnXcare >Decision Tree: Neck contrast medium necessary [factory]	Topogram			ScanplanningTopoAdult		6 x 0.70				130		30	SPR (32 cm) 0.12	3.01	Standard					7.0																
		Topogram 0.70 Tr20 sag																														Tr20					
		Neck		Neck	RoutineSpiralAdultHeadsSnXcare	Neck contrast medium necessary [factory]	32 x 0.70	On	155	Non-Contrast	130	102	(32 cm) 10.9	229	Standard		0.80	1.00		7.0																	
		Neck 5.00 Br40 I3 sag																					5.00	5.00		Neck contrast medium necessary [factory].wo CM; Neck contrast medium necessary [factory].w CM	None	Orientation Multi Recons (3)	MPR	Br40	On	3					
		Neck 5.00 Br40 I3 cor																																			
	Neck 5.00 Br40 I3 ax																																				

Protocol Name	Protocol Description	Range Name	Recon Name	Organ Characteristic	Scan Mode	Decision Trees	Acquisition [mm]	CARE Dose4D	FAST kV IQ level	FAST kV optimized for	kV	Eff. mAs	mA	CTDIvol [mGy]	DLP [mGy*cm]	Shaped Filter	Feed [mm]	Pitch	Rot. Time [s]	Scan Time [s]	Absolute Start Delay [s]	Contrast Protocol	Slice Thickness [mm]	Increment [mm]	Recon Strategies			Inline Results	Multi Recons	Image Type	Kernel	Iterative Recon	IR Strength	Recon Type	Recon Phase Step	Recon Phase(s)	Recon Phase Unit		
Shoulder [factory] (Adult)	>Indications: evaluation of joint cavities, trauma, dislocations, orthopedic indications, shoulder joint, w/o CM >Scan modes: RoutineSpiralAdultShoulder >Decision Tree: Side [factory]	Topogram			ScanplanningTopoAdult		6 x 0.70				110		15	SPR (32 cm) 0.04	0.98	Standard						3.0																	
		Topogram 0.70 Tr20 cor																													Tr20								
		Shoulder		Shoulder	RoutineSpiralAdultShoulder	Side [factory]	32 x 0.70	On	130	Auto (Bone/Calcium)	130	84		(32 cm) 9.01	161	Standard		0.80	1.00		3.0																		
		Shoulder 5.00 Br60 I3 ax																					5.00	5.00		Side [factory].recon both; Side [factory].recon right; Side [factory].recon left	None	None	MPR	Br60	On	3							
		Shoulder 3.00 Br60 I3 sag R																					3.00	3.00		Side [factory].recon both; Side [factory].recon right	Parallel Ranges	Orientation Multi Recons (3)	MPR	Br60	On	3							
		Shoulder 3.00 Br60 I3 cor R																																					
		Shoulder 3.00 Br60 I3 ax R																																					
		Shoulder 3.00 Br40 I3 sag R																																					
		Shoulder 3.00 Br40 I3 cor R																																					
		Shoulder 3.00 Br40 I3 ax R																																					
		Shoulder 3.00 Br60 I3 sag L																					3.00	3.00		Side [factory].recon both; Side [factory].recon left	Parallel Ranges	Orientation Multi Recons (3)	MPR	Br60	On	3							
		Shoulder 3.00 Br60 I3 cor L																																					
		Shoulder 3.00 Br60 I3 ax L																																					
		Shoulder 3.00 Br40 I3 sag L																																					
		Shoulder 3.00 Br40 I3 cor L																																					
		Shoulder 3.00 Br40 I3 ax L																																					
Thorax [factory] (Adult)	>Indications: masses, tumor, metastases, lymphoma, lymph nodes, tuberculosis, empyema, infection, early visualization of pulmonary nodules, screening of the lung, w/o CM >Scan modes: RoutineSpiralAdultThorax >Decision Tree: Thorax FoV [factory]	Topogram			ScanplanningTopoAdult		6 x 0.70				110		15	SPR (32 cm) 0.03	1.60	Standard						6.0																	
		Topogram 0.70 Tr20 cor																													Tr20								
		Thorax		Thorax	RoutineSpiralAdultThorax	Thorax FoV [factory]	32 x 0.70	On	80	Non-Contrast	130	51		(32 cm) 5.47	180	Standard		0.80	0.80		6.0																		
		Thorax 5.00 Br40 I3 cor																					5.00	5.00		Thorax FoV [factory].standard acquisition; Thorax FoV [factory].fast acquisition	None	Orientation Multi Recons (2)	MPR	Br40	On	3							
		Thorax 5.00 Br40 I3 ax																																					
		Thorax 5.00 BI60 I3 sag																					5.00	5.00		Thorax FoV [factory].standard acquisition; Thorax FoV [factory].fast acquisition	None	Orientation Multi Recons (3)	MIP Thin	BI60	On	3							
		Thorax 5.00 BI60 I3 cor																																					
		Thorax 5.00 BI60 I3 ax																																					
		Thorax 10.00 BI60 I3 cor																					10.00	10.00		Thorax FoV [factory].standard acquisition; Thorax FoV [factory].fast acquisition	None	None	MIP Thin	BI60	On	3							
		Thorax 1.00 BI60 I3																					1.00	0.70		Thorax FoV [factory].standard acquisition; Thorax FoV [factory].fast acquisition	Lung CAD	None		BI60	On	3							

Protocol Name	Protocol Description	Range Name	Recon Name	Organ Characteristic	Scan Mode	Decision Trees	Acquisition [mm]	CARE Dose4D	FAST kV IQ level	FAST kV optimized for	kV	Eff. mAs	mA	CTDIvol [mGy]	DLP [mGy*cm]	Shaped Filter	Feed [mm]	Pitch	Rot. Time [s]	Scan Time [s]	Absolute Start Delay [s]	Contrast Protocol	Slice Thickness [mm]	Increment [mm]	Recon Strategies				Inline Results	Multi Recons	Image Type	Kernel	Iterative Recon	IR Strength	Recon Type	Recon Phase Step	Recon Phase(s)	Recon Phase Unit			
Thorax Dynamic [factory] (Adult)	>Indications: abnormalities, upper airway obstruction, trachea bronchomalacia, pulmonary air trapping, abnormal diaphragmatic motion, tumor invasion, routine scan in inspiration w/ CM, dynamic scan during full breathing w/ CM >Scan modes: PremonitoringSequenceAdult2d, MonitoringSequenceAdult2d, RoutineSpiralAdultThorax, DynamicSequenceAdultBody >Decision Tree: Thorax FoV [factory] >Contrast Protocol: Thorax [factory]	Topogram			ScanplanningTopoAdult		6 x 0.70				110		15	SPR (32 cm) 0.03	1.60	Standard					6.0																				
			Topogram 0.70 Tr20 cor																														Tr20								
		Premonitoring		Vascular Body			1 x 10.00	On	25	Soft Tissue with Contrast	130	17			(32 cm) 1.67	1.67	Standard	0.0		0.80	Full 0.80	3.0																			
			Premonitoring 10.00 Br36																					10.00						None		Br36									
		Monitoring		Vascular Body			1 x 10.00	On	25	Soft Tissue with Contrast	130	17			(32 cm) 50.2	50.2	Standard	0.0		0.80	Full 0.80	10.0	Thorax [factory]																		
			Monitoring 10.00 Br36																					10.00						None		Br36									
		Thorax		Thorax	RoutineSpiralAdultThorax	Thorax FoV [factory]	32 x 0.70	On	80	Auto (Soft Tissue with Contrast)	130	56			(32 cm) 6.01	198	Standard		0.80	0.80			Thorax [factory]																		
			Thorax 5.00 Br40 I3 cor																					5.00	5.00	Thorax FoV [factory].standard acquisition; Thorax FoV [factory].fast acquisition				None	Orientation Multi Recons (2)	MPR	Br40	On	3						
			Thorax 5.00 Br40 I3 ax																					5.00	5.00	Thorax FoV [factory].standard acquisition; Thorax FoV [factory].fast acquisition				None	Orientation Multi Recons (3)	MIP Thin	BI60	On	3						
			Thorax 5.00 BI60 I3 sag																																						
			Thorax 5.00 BI60 I3 cor																																						
			Thorax 5.00 BI60 I3 ax																																						
			Thorax 10.00 BI60 I3 cor																					10.00	10.00	Thorax FoV [factory].standard acquisition; Thorax FoV [factory].fast acquisition				None	None	MIP Thin	BI60	On	3						
			Thorax 1.00 BI60 I3																					1.00	0.70	Thorax FoV [factory].standard acquisition; Thorax FoV [factory].fast acquisition				Lung CAD	None		BI60	On	3						
		Dynamic in free breathing		Abdomen	DynamicSequenceAdultBody		32 x 0.70	Off				130	13			(32 cm) 9.77	21.9	Standard	0.0		0.80	Full 0.80	3.0																		
			Dynamic in free breathing 5.00 BI60 I3 3.60s - 10.80s																						5.00	5.00								BI60	On	3					
		Dynamic in free breathing		Abdomen	DynamicSequenceAdultBody		32 x 0.70	Off				130	13			(32 cm) 9.77	21.9	Standard	0.0		0.80	Full 0.80	3.0																		
			Dynamic in free breathing 5.00 BI60 I3 3.60s - 10.80s																						5.00	5.00									BI60	On	3				
		Dynamic in free breathing		Abdomen	DynamicSequenceAdultBody		32 x 0.70	Off				130	13			(32 cm) 9.77	21.9	Standard	0.0		0.80	Full 0.80	3.0																		
			Dynamic in free breathing 5.00 BI60 I3 3.60s - 10.80s																						5.00	5.00									BI60	On	3				
		Dynamic in free breathing		Abdomen	DynamicSequenceAdultBody		32 x 0.70	Off				130	13			(32 cm) 9.77	21.9	Standard	0.0		0.80	Full 0.80	3.0																		
			Dynamic in free breathing 5.00 BI60 I3 3.60s - 10.80s																						5.00	5.00									BI60	On	3				

Protocol Name	Protocol Description	Range Name	Recon Name	Organ Characteristic	Scan Mode	Decision Trees	Acquisition [mm]	CARE Dose4D	FAST kV IQ level	FAST kV optimized for	kV	Eff. mAs	CTDIvol [mGy]	DLP [mGy*cm]	Shaped Filter	Feed [mm]	Pitch	Rot. Time [s]	Scan Time [s]	Absolute Start Delay [s]	Contrast Protocol	Slice Thickness [mm]	Increment [mm]	Recon Strategies			Inline Results	Multi Recons	Image Type	Kernel	Iterative Recon	IR Strength	Recon Type	Recon Phase Step	Recon Phase(s)	Recon Phase Unit	
Thorax and Abdomen [factory] (Adult)	>Indications: masses, tumor, lymphoma, lymph nodes of the lung, liver, pancreas, gall bladder, spleen, kidneys, w CM >Scan modes: RoutineSpiralAdultAbdomenFlexdose >Decision Tree: Thorax and Abdomen contrast medium necessary [factory]	Topogram			ScanplanningTopoAdult		6 x 0.70				110		15	SPR (32 cm) 0.03	2.83	Standard					6.0																
			Topogram 0.70 Tr20 cor																																		
		Thorax and Abdomen		Abdomen	RoutineSpiralAdultAbdomenFlexdose	Thorax and Abdomen contrast medium necessary [factory]	32 x 0.70	On	145	Non-Contrast	130	74	(32 cm) 7.97	501	Standard		0.80	0.80		6.0																	
			Thorax and Abdomen 5.00 Br40 I3 cor																			5.00	5.00		Thorax and Abdomen contrast medium necessary [factory].wo CM; Thorax and Abdomen contrast medium necessary [factory].w CM	None	Orientation Multi Recons (2)	MPR	Br40	On	3						
			Thorax and Abdomen 5.00 Br40 I3 ax																																		
			Thorax and Abdomen 5.00 BI60 I3 sag Lung																			5.00	5.00		Thorax and Abdomen contrast medium necessary [factory].wo CM; Thorax and Abdomen contrast medium necessary [factory].w CM	None	Orientation Multi Recons (3)	MIP Thin	BI60	On	3						
			Thorax and Abdomen 5.00 BI60 I3 cor Lung																																		
			Thorax and Abdomen 5.00 BI60 I3 ax Lung																																		
			Thorax and Abdomen 1.00 BI60 I3																			1.00	0.70		Thorax and Abdomen contrast medium necessary [factory].wo CM; Thorax and Abdomen contrast medium necessary [factory].w CM	Lung CAD	None			BI60	On	3					
		Thorax and Abdomen Multi Phase [factory] (Adult)	>Indications: masses, tumor, lymphoma, lymph nodes of the lung, liver, pancreas, gall bladder, spleen, kidneys, w CM, 2 phases preconfigured >Scan modes: PremonitoringSequenceAdult2d, MonitoringSequenceAdult2d, RoutineSpiralAdultThorax, RoutineSpiralAdultAbdomen >Contrast Protocol: Thorax [factory]	Topogram			ScanplanningTopoAdult		6 x 0.70				110		15	SPR (32 cm) 0.03	2.83	Standard					6.0														
	Topogram 0.70 Tr20 cor																																				
Premonitoring				Vascular Body			1 x 10.00	On	25	Soft Tissue with Contrast	110	22	(32 cm) 1.41	1.41	Standard	0.0		0.80	Full 0.80	3.0																	
Premonitoring	PremonitoringSequenceAdult2d, MonitoringSequenceAdult2d, RoutineSpiralAdultThorax, RoutineSpiralAdultAbdomen						1 x 10.00	On	25	Soft Tissue with Contrast	110	22	(32 cm) 42.3	42.3	Standard	0.0		0.80	Full 0.80	10.0	Thorax [factory]																
Monitoring	Monitoring 10.00 Br36																					10.00					None			Br36							
Monitoring	Monitoring 10.00 Br36																					10.00					None				Br36						
Arterial Phase Thorax				Thorax	RoutineSpiralAdultThorax		32 x 0.70	On	80	Auto (Soft Tissue with Contrast)	130	56	(32 cm) 6.01	204	Standard		1.20	0.80			Thorax [factory]																
	Arterial Phase Thorax 5.00 Br40 I3 cor																					5.00	5.00				None	Orientation Multi Recons (2)	MPR	Br40	On	3					
	Arterial Phase Thorax 5.00 Br40 I3 ax																																				
	Arterial Phase Thorax 5.00 BI60 I3 sag																					5.00	5.00				None	Orientation Multi Recons (3)	MIP Thin	BI60	On	3					
	Arterial Phase Thorax 5.00 BI60 I3 cor																																				
	Arterial Phase Thorax 5.00 BI60 I3 ax																																				
	Arterial Phase Thorax 1.00 BI60 I3																					1.00	0.70				Lung CAD	None			BI60	On	3				
Venous Phase Abdomen				Abdomen	RoutineSpiralAdultAbdomen		32 x 0.70	On	145	Auto (Soft Tissue with Contrast)	130	100	(32 cm) 10.7	246	Standard		0.80	0.80			Thorax [factory]																
	Venous Phase Abdomen 5.00 Br40 I3 cor																					5.00	5.00				None	Orientation Multi Recons (2)	MPR	Br40	On	3					
	Venous Phase Abdomen 5.00 Br40 I3 ax																																				

Protocol Name	Protocol Description	Range Name	Recon Name	Organ Characteristic	Scan Mode	Decision Trees	Acquisition [mm]	CARE Dose4D	FAST KV IQ level	FAST KV optimized for	kV	Eff. mAs	CTDIvol [mGy]	DLP [mGy*cm]	Shaped Filter	Feed [mm]	Pitch	Rot. Time [s]	Scan Time [s]	Absolute Start Delay [s]	Contrast Protocol	Slice Thickness [mm]	Increment [mm]	Recon Strategies			Inline Results	Multi Recons	Image Type	Kernel	Iterative Recon	IR Strength	Recon Type	Recon Phase Step	Recon Phase(s)	Recon Phase Unit
Lung Screening [factory] (Adult)	>Indications: masses, tumor, metastases, lymphoma, lymph nodes, tuberculosis, empyema, infection, early visualization of pulmonary nodules, screening of the lung, w/wo CM >Scan modes: RoutineSpiralAdultThoraxSn >Decision Tree: Lung screening [factory]	Topogram			ScanplanningTopoAdult		6 x 0.70				110	15	SPR (32 cm) 0.03	1.60	Standard					6.0																
		Topogram 0.70 Tr20 cor																													Tr20					
		Lung		Thorax	RoutineSpiralAdultThoraxSn	Lung screening [factory]	32 x 0.70	On	15	Non-Contrast	Sn 110	62	(32 cm) 0.88	29.1	Standard		0.80	0.80		6.0																
		Lung 5.00 Br40 I3 cor																				5.00	5.00		Lung screening [factory].low dose wo X-CARE; Lung screening [factory].low dose w X-CARE; Lung screening [factory].high resolution wo X-CARE; Lung screening [factory].low dose w X-CARE	None	Orientation Multi Recons (2)	MPR	Br40	On	3					
		Lung 5.00 Br40 I3 ax																							Lung screening [factory].low dose wo X-CARE; Lung screening [factory].high resolution wo X-CARE; Lung screening [factory].low dose w X-CARE	None	Orientation Multi Recons (3)	MIP Thin	BI60	On	3					
		Lung 5.00 BI60 I3 sag																				5.00	5.00													
		Lung 5.00 BI60 I3 cor																																		
		Lung 5.00 BI60 I3 ax																																		
		Lung 2.00 BI64 I3 sag																				2.00	2.00		Lung screening [factory].high resolution wo X-CARE; Lung screening [factory].high resolution w X-CARE	None	Orientation Multi Recons (3)	MIP Thin	BI64	On	3					
		Lung 2.00 BI64 I3 cor																																		
		Lung 2.00 BI64 I3 ax																																		
Lung 10.00 BI60 I3 cor																					10.00	10.00		Lung screening [factory].low dose wo X-CARE; Lung screening [factory].low dose w X-CARE; Lung screening [factory].high resolution wo X-CARE; Lung screening [factory].high resolution w X-CARE	None	None	MIP Thin	BI60	On	3						
Lung 1.00 BI60 I3																					1.00	0.70		Lung screening [factory].low dose wo X-CARE; Lung screening [factory].low dose w X-CARE; Lung screening [factory].high resolution wo X-CARE; Lung screening [factory].high resolution w X-CARE	Lung CAD	None		BI60	On	3						
Lung Low Dose [factory] (Adult)	>Indications: early visualization of pulmonary nodules, screening of the lung, wo CM >Scan modes: RoutineSpiralAdultThoraxSn >Decision Tree: Lung Low Dose breathhold and X-CARE [factory]	Topogram			ScanplanningTopoAdult		6 x 0.70				110	15	SPR (32 cm) 0.03	1.60	Standard					6.0																
		Topogram 0.70 Tr20 cor																													Tr20					
		Lung		Thorax	RoutineSpiralAdultThoraxSn		32 x 0.70	On	15	Auto (Non-Contrast)	Sn 110	62	(32 cm) 0.88	30.0	Standard		1.20	0.80		6.0																
		Lung 5.00 Br40 I3 sag																				5.00	5.00			None	Orientation Multi Recons (3)	MPR	Br40	On	3					
		Lung 5.00 Br40 I3 cor																																		
		Lung 5.00 Br40 I3 ax																																		
		Lung 5.00 BI60 I3 sag																				5.00	5.00			None	Orientation Multi Recons (3)	MIP Thin	BI60	On	3					
		Lung 5.00 BI60 I3 cor																																		
Lung 5.00 BI60 I3 ax																																				
Lung 10.00 BI60 I3 cor																					10.00	10.00				None	None	MIP Thin	BI60	On	3					
Lung 1.00 BI60 I3																					1.00	0.70			Lung CAD	None		BI60	On	3						
Abdomen [factory] (Adult)	>Indications: masses, tumor, lymphoma, lymph nodes of the liver, pancreas, gall bladder, spleen, kidneys, intestine, w/wo CM >Scan modes: RoutineSpiralAdultAbdomen >Decision Tree: Abdomen contrast medium necessary [factory]	Topogram			ScanplanningTopoAdult		6 x 0.70				110	15	SPR (32 cm) 0.03	1.60	Standard					6.0																
		Topogram 0.70 Tr20 cor																												Tr20						
		Abdomen		Abdomen	RoutineSpiralAdultAbdomen	Abdomen contrast medium necessary [factory]	32 x 0.70	On	145	Non-Contrast	130	92	(32 cm) 9.87	325	Standard		0.80	0.80		6.0																
		Abdomen 5.00 Br40 I3 cor																								Abdomen contrast medium necessary [factory].wo CM; Abdomen contrast medium necessary [factory].w CM	None	Orientation Multi Recons (2)	MPR	Br40	On	3				

Protocol Name	Protocol Description	Range Name	Recon Name	Organ Characteristic	Scan Mode	Decision Trees	Acquisition [mm]	CARE Dose4D	FAST KV IQ level	FAST kv optimized for	kV	Eff. mAs	CTDIvol [mGy]	DLP [mGy*cm]	Shaped Filter	Feed [mm]	Pitch	Rot. Time [s]	Scan Time [s]	Absolute Start Delay [s]	Contrast Protocol	Slice Thickness [mm]	Increment [mm]	Recon Strategies			Inline Results	Multi Recons	Image Type	Kernel	Iterative Recon	IR Strength	Recon Type	Recon Phase Step	Recon Phase(s)	Recon Phase Unit		
Abdomen Multi Phase [factory] (Adult)	>Indications: masses, tumor, lymphoma, lymph nodes of the liver, pancreas, gall bladder, spleen, kidneys, intestine, w and wo CM, 3 phases preconfigured >Scan modes: RoutineSpiralAdultAbdomen, PremonitoringSequenceAdult2d, MonitoringSequenceAdult2d, RoutineSpiralAdultAbdomen >Contrast Protocol: Abdomen [factory]	Topogram			ScanplanningTopoAdult		6 x 0.70				110	15	SPR (32 cm) 0.03	1.60	Standard					6.0																		
			Topogram 0.70 Tr20 cor																											Tr20								
		Non Contrast Liver		Abdomen	RoutineSpiralAdultAbdomen		32 x 0.70	On	145	Auto (Non-Contrast)	130	92	(32 cm) 9.87	226	Standard		0.80	0.80		6.0				5.00	5.00		None	Orientation Multi Recons (2)	MPR	Br40	On	3						
			Non Contrast Liver 5.00 Br40 I3 cor																																			
			Non Contrast Liver 5.00 Br40 I3 ax																																			
		Premonitoring		Vascular Body			1 x 10.00	On	25	Soft Tissue with Contrast	110	22	(32 cm) 1.41	1.41	Standard	0.0		0.80	Full 0.80	3.0																		
			Premonitoring 10.00 Br36																				10.00					None			Br36							
		Monitoring		Vascular Body			1 x 10.00	On	25	Soft Tissue with Contrast	110	22	(32 cm) 42.3	42.3	Standard	0.0		0.80	Full 0.80	10.0	Abdomen [factory]																	
			Monitoring 10.00 Br36																					10.00					None			Br36						
		Arterial Phase Abdomen		Abdomen	RoutineSpiralAdultAbdomen		32 x 0.70	On	145	Auto (Soft Tissue with Contrast)	130	100	(32 cm) 10.7	353	Standard		0.80	0.80			Abdomen [factory]																	
			Arterial Phase Abdomen 5.00 Br40 I3 cor																					5.00	5.00		None	Orientation Multi Recons (2)	MPR	Br40	On	3						
			Arterial Phase Abdomen 5.00 Br40 I3 ax																																			
		Venous Phase Abdomen		Abdomen	RoutineSpiralAdultAbdomen		32 x 0.70	On	145	Auto (Soft Tissue with Contrast)	130	100	(32 cm) 10.7	246	Standard		0.80	0.80		84.6	Abdomen [factory]																	
			Venous Phase Abdomen 5.00 Br40 I3 cor																					5.00	5.00		None	Orientation Multi Recons (2)	MPR	Br40	On	3						
			Venous Phase Abdomen 5.00 Br40 I3 ax																																			
Kidney DE Stones [factory] (Adult)	>Indications: query and analyze kidney stones, wo CM >Scan modes: RoutineSpiralAdultAbdomenTwinSpiral (max. 100 cm) >Decision Tree: Kidney TwinSpiral DE possible [factory]	Topogram			ScanplanningTopoAdult		6 x 0.70				110	15	SPR (32 cm) 0.03	1.60	Standard					6.0																		
			Topogram 0.70 Tr20 cor																											Tr20								
		Kidney Stones		Abdomen	RoutineSpiralAdultAbdomenTwinSpiral (max. 100 cm)	Kidney TwinSpiral DE possible [factory]	32 x 0.70	On	145	Auto (Non-Contrast)			(32 cm) 9.22	303	Standard		0.80	0.80		11.0																		
			Kidney Stones 5.00 Br40 I3 cor																				5.00	5.00		Kidney TwinSpiral DE possible [factory].TwinSpiral DE; Kidney TwinSpiral DE possible [factory].Single Energy	None	Orientation Multi Recons (2)	MPR	Br40	On	3						
			Kidney Stones 5.00 Br40 I3 ax																																			
			Kidney Stones 1.50 Qr40 I3 sag																				1.50	1.00		Kidney TwinSpiral DE possible [factory].TwinSpiral DE	DE Parallel Ranges	Orientation Multi Recons (3)	MPR	Qr40	On	3						
			Kidney Stones 1.50 Qr40 I3 cor																																			
			Kidney Stones 1.50 Qr40 I3 ax																																			
			Kidney Stones 5.00 Qr40 I3 cor																				5.00	5.00		Kidney TwinSpiral DE possible [factory].TwinSpiral DE	None	None	MPR	Qr40	On	3						

Protocol Name	Protocol Description	Range Name	Recon Name	Organ Characteristic	Scan Mode	Decision Trees	Acquisition [mm]	CARE Dose4D	FAST KV IQ level	FAST KV optimized for	kV	Eff. mAs	mA	CTDIvol [mGy]	DLP [mGy*cm]	Shaped Filter	Feed [mm]	Pitch	Rot. Time [s]	Scan Time [s]	Absolute Start Delay [s]	Contrast Protocol	Slice Thickness [mm]	Increment [mm]	Recon Strategies				Inline Results	Multi Recons	Image Type	Kernel	Iterative Recon	IR Strength	Recon Type	Recon Phase Step	Recon Phase(s)	Recon Phase Unit			
Colon Low Dose [factory] (Adult)	>Indications: early visualization of diverticles, polyps, nodules, screening of the colon, w/o CM >Scan modes: RoutineSpiralAdultAbdomenSn	Topogram			ScanplanningTopoAdult		6 x 0.70				110		15	SPR (32 cm) 0.03	1.60	Standard						6.0																			
		Topogram 0.70 Tr20 cor																																Tr20							
		Supine Colon		Abdomen	RoutineSpiralAdultAbdomenSn		32 x 0.70	On	39	Non-Contrast	Sn 110	152		(32 cm) 2.17	93.0	Standard		0.80	0.80		6.0																				
		Supine Colon 5.00 Br40 I3 cor																						5.00	5.00			None	Orientation Multi Recons (2)	MPR	Br40	On	3								
			Supine Colon 5.00 Br40 I3 ax																																						
		Supine Colon 1.50 Br40 I3																							1.50	1.00			None	None		Br40	On	3							
		Topogram			ScanplanningTopoAdult		6 x 0.70					110		15	SPR (32 cm) 0.03	1.60	Standard						6.0																		
		Topogram 0.70 Tr20 cor																																			Tr20				
		Prone Colon		Abdomen	RoutineSpiralAdultAbdomenSn		32 x 0.70	On	20	Non-Contrast	Sn 110	81		(32 cm) 1.16	49.5	Standard			0.80	0.80		6.0																			
		Prone Colon 5.00 Br40 I3 cor																							5.00	5.00			None	Orientation Multi Recons (2)	MPR	Br40	On	3							
Prone Colon 5.00 Br40 I3 ax																																									
Prone Colon 1.50 Br40 I3																							1.50	1.00			None	None		Br40	On	3									
Whole Body [factory] (Adult)	>Indications: degenerative changes, trauma, osteolysis, masses, tumor, lymphoma, lymph nodes of the torso region, w/o CM >Scan modes: RoutineSpiralAdultAbdomen >Decision Tree: Whole Body soft tissue or bone [factory]	Topogram			ScanplanningTopoAdult		6 x 0.70				110		15	SPR (32 cm) 0.03	2.83	Standard						6.0																			
		Topogram 0.70 Tr20 cor																																Tr20							
		Whole Body		Abdomen	RoutineSpiralAdultAbdomen	Whole Body soft tissue or bone [factory]	32 x 0.70	On	145	Non-Contrast	130 92			(32 cm) 9.87	621	Standard		0.80	0.80		6.0																				
		Whole Body 5.00 Br40 I3 cor																						5.00	5.00			Whole Body soft tissue or bone [factory].standard dose	None	Orientation Multi Recons (2)	MPR	Br40	On	3							
			Whole Body 5.00 Br40 I3 ax																																						
		Whole Body 5.00 BI60 I3 sag Lung																						5.00	5.00			Whole Body soft tissue or bone [factory].standard dose	None	Orientation Multi Recons (3)	MIP Thin	BI60	On	3							
		Whole Body 5.00 BI60 I3 cor Lung																																							
		Whole Body 5.00 BI60 I3 ax Lung																																							
		Whole Body 1.00 BI60 I3																						1.00	0.70			Whole Body soft tissue or bone [factory].standard dose	Lung CAD	None		BI60	On	3							
		Whole Body 5.00 Br60 I3 sag																							5.00	5.00			Whole Body soft tissue or bone [factory].low dose	None	Orientation Multi Recons (3)	MPR	Br60	On	3						
			Whole Body 5.00 Br60 I3 cor																																						
		Whole Body 5.00 Br60 I3 ax																																							
		Whole Body 5.00 Br60 I3																																							
		Whole Body 5.00 Br60 I3 ax																																							
Whole Body 5.00 Br40 I3 ax																																									

Protocol Name	Protocol Description	Range Name	Recon Name	Organ Characteristic	Scan Mode	Decision Trees	Acquisition [mm]	CARE Dose4D	FAST KV IQ level	FAST KV optimized for	kV	Eff. mAs	CTDIvol [mGy]	DLP [mGy*cm]	Shaped Filter	Feed [mm]	Pitch	Rot. Time [s]	Scan Time [s]	Absolute Start Delay [s]	Contrast Protocol	Slice Thickness [mm]	Increment [mm]	Recon Strategies			Inline Results	Multi Recons	Image Type	Kernel	Iterative Recon	IR Strength	Recon Type	Recon Phase Step	Recon Phase(s)	Recon Phase Unit
Pelvis [factory] (Adult)	>Indications: masses, tumor, lymphom of the urinary bladder, rectum, uterus, prostate, degenerative changes, trauma, osteolysis of the pelvis, w/wo CM >Scan modes: RoutineSpiralAdultPelvis >Decision Tree: Metal implant [factory]; Pelvis soft tissue or bone [factory]	Topogram			ScanplanningTopoAdult		6 x 0.70				110	15	SPR (32 cm) 0.04	0.98	Standard					3.0																
			Topogram 0.70 Tr20 cor																												Tr20					
		Pelvis		Pelvis	RoutineSpiralAdultPelvis	Metal implant [factory]; Pelvis soft tissue or bone [factory]	32 x 0.70	On	170	Non-Contrast	130	107	(32 cm) 11.5	263	Standard		0.80	1.00		3.0																
			Pelvis 5.00 Br40 I3 sag																			5.00	5.00		Metal implant [factory].wo iMAR; Metal implant [factory].w iMAR; Pelvis soft tissue or bone [factory].standard dose	None	Orientation Multi Recons (3)	MPR	Br40	On	3					
			Pelvis 5.00 Br40 I3 cor																																	
			Pelvis 5.00 Br40 I3 ax																																	
			Pelvis 5.00 Br60 I3 ax																			5.00	5.00		Metal implant [factory].wo iMAR; Metal implant [factory].w iMAR; Pelvis soft tissue or bone [factory].standard dose	None	None	MPR	Br60	On	3					
			Pelvis 5.00 Br60 I3 sag																			5.00	5.00		Metal implant [factory].wo iMAR; Metal implant [factory].w iMAR; Pelvis soft tissue or bone [factory].standard dose	None	Orientation Multi Recons (3)	MPR	Br60	On	3					
			Pelvis 5.00 Br60 I3 cor																																	
			Pelvis 5.00 Br60 I3 ax																																	
			Pelvis 5.00 Br40 I3 ax																			5.00	5.00		Metal implant [factory].wo iMAR; Metal implant [factory].w iMAR; Pelvis soft tissue or bone [factory].low dose	None	None	MPR	Br40	On	3					
			Pelvis Br60 I3 cor Radial																						Metal implant [factory].wo iMAR; Metal implant [factory].w iMAR; Pelvis soft tissue or bone [factory].low dose	Radial Ranges	None	VRT	Br60	On	3					
			Pelvis 5.00 Br40 I3 sag																			5.00	5.00		Metal implant [factory].w iMAR; Pelvis soft tissue or bone [factory].standard dose; Pelvis soft tissue or bone [factory].low dose	None	Orientation Multi Recons (3)	MPR	Br40	On	3					
			Pelvis 5.00 Br40 I3 cor																																	
	Pelvis 5.00 Br40 I3 ax																																			
Hip [factory] (Adult)	>Indications: evaluation of joint cavity, trauma, dysplasia, necrosis of the head of the hip joints, congruence evaluations, orthopedic indications, wo CM >Scan modes: RoutineSpiralAdultPelvis >Decision Tree: Side [factory]; Metal implant [factory]	Topogram			ScanplanningTopoAdult		6 x 0.70				110	15	SPR (32 cm) 0.04	0.98	Standard					3.0																
			Topogram 0.70 Tr20 cor																											Tr20						
		Hip		Pelvis	RoutineSpiralAdultPelvis	Side [factory]; Metal implant [factory]	32 x 0.70	On	170	Bone/Calcium	130	109	(32 cm) 11.7	151	Standard		0.80	1.00		3.0																
			Hip 3.00 Br60 I3 ax																			3.00	3.00		Side [factory].recon both; Side [factory].recon right; Side [factory].recon left; Metal implant [factory].wo iMAR; Metal implant [factory].w iMAR	None	Impression Multi Recons (2)	MPR	Br60	On	3					
			Hip 3.00 Br40 I3 ax																										Br40	On	3					
			Hip 3.00 Br60 I3 cor R																			3.00	3.00		Side [factory].recon both; Side [factory].recon right; Metal implant [factory].wo iMAR; Metal implant [factory].w iMAR	Parallel Ranges	None	MPR	Br60	On	3					
			Hip Br40 I3 cor Radial R																						Side [factory].recon both; Side [factory].recon right; Metal implant [factory].wo iMAR; Metal implant [factory].w iMAR	Radial Ranges	None	VRT	Br40	On	3					
			Hip 3.00 Br60 I3 cor L																			3.00	3.00		Side [factory].recon both; Side [factory].recon left; Metal implant [factory].wo iMAR; Metal implant [factory].w iMAR	Parallel Ranges	None	MPR	Br60	On	3					
			Hip Br40 I3 cor Radial L																						Side [factory].recon both; Side [factory].recon left; Metal implant [factory].wo iMAR; Metal implant [factory].w iMAR	Radial Ranges	None	VRT	Br40	On	3					
			Hip 3.00 Br40 I3 ax																			3.00	2.00		Side [factory].recon both; Side [factory].recon right; Side [factory].recon left; Metal implant [factory].w iMAR	None	None	MPR	Br40	On	3					

Protocol Name	Protocol Description	Range Name	Recon Name	Organ Characteristic	Scan Mode	Decision Trees	Acquisition [mm]	CARE Dose4D	FAST KV IQ level	FAST KV optimized for	kV	Eff. mAs	mA	CTDIvol [mGy]	DLP [mGy*cm]	Shaped Filter	Feed [mm]	Pitch	Rot. Time [s]	Scan Time [s]	Absolute Start Delay [s]	Contrast Protocol	Slice Thickness [mm]	Increment [mm]	Recon Strategies	Inline Results	Multi Recons	Image Type	Kernel	Iterative Recon	IR Strength	Recon Type	Recon Phase Step	Recon Phase(s)	Recon Phase Unit	
C-Spine [factory] (Adult)	>Indications: prolapse/hernia, tumors, degenerative changes, bone tumor, trauma, osteolysis of the cervical spine, wo CM >Scan modes: RoutineSpiralAdultNeck >Decision Tree: C>Spine disc or bone [factory]	Topogram			ScanplanningTopoAdult		6 x 0.70				130		30	SPR (32 cm) 0.12	3.01	Standard					3.0															
		Topogram 0.70 Tr20 sag																												Tr20						
		C-Spine		Neck	RoutineSpiralAdultNeck	C-Spine disc or bone [factory]	32 x 0.70	On	195	Non-Contrast	130	129		(32 cm) 13.8	261	Standard		0.80	1.00		3.0															
		C-Spine 2.00 Br40 I3 ax																				2.00	2.00		C-Spine disc or bone [factory].standard dose	None	Impression Multi Recons (2)	MPR	Br40	On	3					
		C-Spine 2.00 Br60 I3 ax																											Br60	On	3					
		C-Spine 2.00 Br40 I3 ax Spine																				2.00	2.00		C-Spine disc or bone [factory].standard dose	Spine Ranges	None	MPR	Br40	On	3					
		C-Spine 2.00 Br40 I3 sag Spine																				2.00	2.00		C-Spine disc or bone [factory].standard dose	Spine Ranges	None	MPR	Br40	On	3					
		C-Spine 2.00 Br60 I3 ax																				2.00	2.00		C-Spine disc or bone [factory].low dose	None	Impression Multi Recons (2)	MPR	Br60	On	3					
		C-Spine 2.00 Br40 I3 ax																											Br40	On	3					
		C-Spine 2.00 Br60 I3 ax Spine																				2.00	2.00		C-Spine disc or bone [factory].low dose	Spine Ranges	None	MPR	Br60	On	3					
		C-Spine 2.00 Br60 I3 sag Spine																											Br60	On	3					
		C-Spine Br40 I3 sag Radial																											Radial Ranges	None	VRT	Br40	On	3		
T-Spine [factory] (Adult)	>Indications: prolapse/hernia, degenerative changes, osteolysis, trauma, tumors of the thoracic spine, wo CM >Scan modes: RoutineSpiralAdultSpine >Decision Tree: T> and L>Spine disc or bone [factory]	Topogram			ScanplanningTopoAdult		6 x 0.70				130		30	SPR (32 cm) 0.10	4.90	Standard					3.0															
		Topogram 0.70 Tr20 sag																												Tr20						
		T-Spine		Spine	RoutineSpiralAdultSpine	T- and L-Spine disc or bone [factory]	32 x 0.70	On	285	Non-Contrast	130	182		(32 cm) 19.5	369	Standard		0.80	1.00		3.0															
		T-Spine 3.00 Br40 I3 ax																				3.00	3.00		T- and L-Spine disc or bone [factory].standard dose	None	Impression Multi Recons (2)	MPR	Br40	On	3					
		T-Spine 3.00 Br60 I3 ax																											Br60	On	3					
		T-Spine 3.00 Br40 I3 ax Spine																				3.00	3.00		T- and L-Spine disc or bone [factory].standard dose	Spine Ranges	None	MPR	Br40	On	3					
		T-Spine 3.00 Br40 I3 sag Spine																				3.00	3.00		T- and L-Spine disc or bone [factory].standard dose	Spine Ranges	None	MPR	Br40	On	3					
		T-Spine 3.00 Br60 I3 ax																											Br40	On	3					
		T-Spine 3.00 Br40 I3 ax																				3.00	3.00		T- and L-Spine disc or bone [factory].low dose	None	Impression Multi Recons (2)	MPR	Br60	On	3					
		T-Spine 3.00 Br60 I3 ax Spine																											Br40	On	3					
		T-Spine 3.00 Br60 I3 sag Spine																				3.00	3.00		T- and L-Spine disc or bone [factory].low dose	Spine Ranges	None	MPR	Br60	On	3					
		T-Spine Br40 I3 sag Radial																											Radial Ranges	None	VRT	Br40	On	3		

Protocol Name	Protocol Description	Range Name	Recon Name	Organ Characteristic	Scan Mode	Decision Trees	Acquisition [mm]	CARE Dose4D	FAST kV IQ level	FAST kV optimized for	kV	Eff. mAs	mAs	CTDIvol [mGy]	DLP [mGy*cm]	Shaped Filter	Feed [mm]	Pitch	Rot. Time [s]	Scan Time [s]	Absolute Start Delay [s]	Contrast Protocol	Slice Thickness [mm]	Increment [mm]	Recon Strategies			Inline Results	Multi Recons	Image Type	Kernel	Iterative Recon	IR Strength	Recon Type	Recon Phase Step	Recon Phase(s)	Recon Phase Unit		
L-Spine [factory] (Adult)	>Indications: prolapse/hernia, tumor, degenerative changes, trauma, osteolysis of the lumbar spine, wo CM >Scan modes: RoutineSpiralAdultSpine >Decision Tree: T- and L-Spine disc or bone [factory]	Topogram			ScanplanningTopoAdult		6 x 0.70				130		30	SPR (32 cm) 0.10	4.90	Standard					3.0																		
			Topogram 0.70 Tr20 sag																													Tr20							
		L-Spine		Spine	RoutineSpiralAdultSpine	T- and L-Spine disc or bone [factory]	32 x 0.70	On	285	Non-Contrast	130	182		(32 cm) 19.5	369	Standard		0.80	1.00		3.0																		
			L-Spine 3.00 Br40 I3 ax																				3.00	3.00		T- and L-Spine disc or bone [factory].standard dose	None	Impression Multi Recons (2)	MPR	Br40	On	3							
			L-Spine 3.00 Br60 I3 ax																							T- and L-Spine disc or bone [factory].standard dose	Spine Ranges	None	MPR	Br40	On	3							
			L-Spine 3.00 Br40 I3 sag Spine																				3.00	3.00		T- and L-Spine disc or bone [factory].standard dose	Spine Ranges	None	MPR	Br40	On	3							
			L-Spine 3.00 Br60 I3 ax																				3.00	3.00		T- and L-Spine disc or bone [factory].low dose	None	Impression Multi Recons (2)	MPR	Br60	On	3							
			L-Spine 3.00 Br40 I3 ax																				3.00	3.00		T- and L-Spine disc or bone [factory].low dose	Spine Ranges	None	MPR	Br60	On	3							
			L-Spine 3.00 Br60 I3 sag Spine																				3.00	3.00		T- and L-Spine disc or bone [factory].low dose	Spine Ranges	None	MPR	Br60	On	3							
			L-Spine 3.00 Br60 I3 sag Spine																				3.00	3.00		T- and L-Spine disc or bone [factory].low dose	Spine Ranges	None	MPR	Br60	On	3							
			L-Spine Br40 I3 sag Radial																								T- and L-Spine disc or bone [factory].low dose	Radial Ranges	None	VRT	Br40	On	3						
Osteo [factory] (Adult)	>Indications: quantitative assessment of vertebral bone mineral density (BMD) in the diagnosis and follow-up of osteopenia and osteoporosis, wo CM >Scan modes: RoutineSequenceAdultOsteo2d	Topogram			ScanplanningTopoAdult		6 x 0.70				130		30	SPR (32 cm) 0.10	4.90	Standard					3.0																		
			Topogram 0.70 Tr20 sag																													Tr20							
		Osteo		Osteo	RoutineSequenceAdultOsteo2d		1 x 10.00	Off			80	219		(32 cm) 5.37	5.37	Standard	0.0		1.00	Full 1.00	3.0																		
			Osteo 10.00 Or40 Evaluation																				10.00						None		Or40								
		Osteo		Osteo	RoutineSequenceAdultOsteo2d		1 x 10.00	Off			80	219		(32 cm) 5.37	5.37	Standard	0.0		1.00	Full 1.00	15.0																		
			Osteo 10.00 Or40 Evaluation																				10.00						None		Or40								
		Osteo		Osteo	RoutineSequenceAdultOsteo2d		1 x 10.00	Off			80	219		(32 cm) 5.37	5.37	Standard	0.0		1.00	Full 1.00	33.0																		
	Osteo 10.00 Or40 Evaluation																				10.00						None		Or40										
Hand [factory] (Adult)	>Indications: trauma, disorders of the hand bones, wo CM >Scan modes: RoutineSpiralAdultExtremitiesSn >Decision Tree: Side [factory]	Topogram			ScanplanningTopoAdult		6 x 0.70				110		15	SPR (32 cm) 0.04	0.98	Standard					3.0																		
			Topogram 0.70 Tr20 cor																													Tr20							
		Hand		Extremities	RoutineSpiralAdultExtremitiesSn	Side [factory]	32 x 0.70	On	75	Auto (Bone/Calcium)	Sn 110	366		(32 cm) 5.22	112	Standard		0.40	1.00		3.0																		
			Hand 3.00 Br60 I3 sag R																				3.00	3.00		Side [factory].recon both; Side [factory].recon right	None	Orientation Multi Recons (3)	MPR	Br60	On	3							
			Hand 3.00 Br60 I3 cor R																																				
			Hand 3.00 Br60 I3 ax R																				3.00	3.00		Side [factory].recon both; Side [factory].recon right	None	None	MPR	Br40	On	3							
			Hand 3.00 Br60 I3 sag L																				3.00	3.00		Side [factory].recon both; Side [factory].recon left	None	Orientation Multi Recons (3)	MPR	Br60	On	3							
			Hand 3.00 Br60 I3 cor L																																				
			Hand 3.00 Br60 I3 ax L																				3.00	3.00		Side [factory].recon both; Side [factory].recon left	None	None	MPR	Br40	On	3							

Protocol Name	Protocol Description	Range Name	Recon Name	Organ Characteristic	Scan Mode	Decision Trees	Acquisition [mm]	CARE Dose4D	FAST KV IQ level	FAST KV optimized for	kV	Eff. mAs	CTDIvol [mGy]	DLP [mGy*cm]	Shaped Filter	Feed [mm]	Pitch	Rot. Time [s]	Scan Time [s]	Absolute Start Delay [s]	Contrast Protocol	Slice Thickness [mm]	Increment [mm]	Recon Strategies			Inline Results	Multi Recons	Image Type	Kernel	Iterative Recon	IR Strength	Recon Type	Recon Phase Step	Recon Phase(s)	Recon Phase Unit
Hand DE Gout [factory] (Adult)	>Indications: visualization, characterization, quantification of monosodium urate crystal deposits, gout, gouty arthritis of the hand bones, wo CM >Scan modes: RoutineSpiralAdultExtremitiesTwinspiral (max. 100 cm) >Decision Tree: Side [factory]	Topogram			ScanplanningTopoAdult		6 x 0.70				110		15	SPR (32 cm) 0.04	0.98	Standard				3.0																
			Topogram 0.70 Tr20 cor																																	
		Hand		Extremities	RoutineSpiralAdultExtremitiesTwinspiral (max. 100 cm)	Side [factory]	32 x 0.70	On	75	Non-Contrast				(32 cm) 4.91	112	Standard		0.80	0.80		3.0															
			Hand 3.00 Br60 I3 sag R																			3.00	3.00	Side [factory].recon both; Side [factory].recon right		None	Orientation Multi Recons (3)	MPR	Br60	On	3					
			Hand 3.00 Br60 I3 cor R																																	
			Hand 3.00 Br60 I3 ax R																																	
			Hand 3.00 Br40 I3 ax R																			3.00	3.00	Side [factory].recon both; Side [factory].recon right		None	None	MPR	Br40	On	3					
			Hand 0.80 Qr40 I3 cor R																			0.80	0.50	Side [factory].recon both; Side [factory].recon right		DE Parallel Ranges	None	MPR	Qr40	On	3					
			Hand 3.00 Br60 I3 sag L																			3.00	3.00	Side [factory].recon both; Side [factory].recon left		None	Orientation Multi Recons (3)	MPR	Br60	On	3					
			Hand 3.00 Br60 I3 cor L																																	
			Hand 3.00 Br60 I3 ax L																																	
			Hand 3.00 Br40 I3 ax L																			3.00	3.00	Side [factory].recon both; Side [factory].recon left		None	None	MPR	Br40	On	3					
			Hand 0.80 Qr40 I3 cor L																			0.80	0.50	Side [factory].recon both; Side [factory].recon left		DE Parallel Ranges	None	MPR	Qr40	On	3					
		Arm or Wrist [factory] (Adult)	>Indications: trauma, disorders of the humerus, radius, ulna bones, wrist joint, wo CM >Scan modes: RoutineSpiralAdultExtremitiesSn >Decision Tree: Side [factory]	Topogram			ScanplanningTopoAdult		6 x 0.70				110		15	SPR (32 cm) 0.04	0.98	Standard				3.0														
	Topogram 0.70 Tr20 cor																																			
Arm or Wrist				Extremities	RoutineSpiralAdultExtremitiesSn	Side [factory]	32 x 0.70	On	75	Auto (Bone/Calcium)	Sn 110	366	(32 cm) 5.22	39.3	Standard		0.40	1.00		3.0																
	Arm or Wrist 2.00 Br60 I3 sag R																					2.00	2.00	Side [factory].recon both; Side [factory].recon right		None	Orientation Multi Recons (3)	MPR	Br60	On	3					
	Arm or Wrist 2.00 Br60 I3 cor R																																			
	Arm or Wrist 2.00 Br60 I3 ax R																																			
	Arm or Wrist 2.00 Br40 I3 ax R																					2.00	2.00	Side [factory].recon both; Side [factory].recon right		None	None	MPR	Br40	On	3					
	Arm or Wrist 2.00 Br60 I3 sag L																					2.00	2.00	Side [factory].recon both; Side [factory].recon left		None	Orientation Multi Recons (3)	MPR	Br60	On	3					
	Arm or Wrist 2.00 Br60 I3 cor L																																			
	Arm or Wrist 2.00 Br60 I3 ax L																																			
	Arm or Wrist 2.00 Br40 I3 ax L																					2.00	2.00	Side [factory].recon both; Side [factory].recon left		None	None	MPR	Br40	On	3					

Protocol Name	Protocol Description	Range Name	Recon Name	Organ Characteristic	Scan Mode	Decision Trees	Acquisition [mm]	CARE Dose4D	FAST KV IQ level	FAST KV optimized for	KV	Eff. mAs	CTDIvol [mGy]	DLP [mGy*cm]	Shaped Filter	Feed [mm]	Pitch	Rot. Time [s]	Scan Time [s]	Absolute Start Delay [s]	Contrast Protocol	Slice Thickness [mm]	Increment [mm]	Recon Strategies			Inline Results	Multi Recons	Image Type	Kernel	Iterative Recon	IR Strength	Recon Type	Recon Phase Step	Recon Phase(s)	Recon Phase Unit	
Elbow [factory] (Adult)	>Indications: trauma, disorders of the elbow joint, wo CM >Scan modes: RoutineSpiralAdultExtremitiesSn >Decision Tree: Side [factory]	Topogram			ScanplanningTopoAdult		6 x 0.70				110		15	SPR (32 cm) 0.04	0.98	Standard					3.0																
		Topogram 0.70 Tr20 cor																													Tr20						
		Elbow		Extremities	RoutineSpiralAdultExtremitiesSn	Side [factory]	32 x 0.70	On	75	Auto (Bone/Calcium)	Sn 110	366	(32 cm) 5.22	60.2	Standard		0.40	1.00		3.0			2.00	2.00	Side [factory].recon both; Side [factory].recon right	None	Orientation Multi Recons (3)	MPR	Br60	On	3						
		Elbow 2.00 Br60 I3 sag R																				2.00	2.00	Side [factory].recon both; Side [factory].recon right	None	None	MPR	Br40	On	3							
		Elbow 2.00 Br60 I3 cor R																				2.00	2.00	Side [factory].recon both; Side [factory].recon right	None	Orientation Multi Recons (3)	MPR	Br60	On	3							
		Elbow 2.00 Br60 I3 ax R																																			
		Elbow 2.00 Br60 I3 sag L																				2.00	2.00	Side [factory].recon both; Side [factory].recon left	None	None	MPR	Br40	On	3							
		Elbow 2.00 Br60 I3 cor L																																			
		Elbow 2.00 Br60 I3 ax L																				2.00	2.00	Side [factory].recon both; Side [factory].recon left	None	None	MPR	Br40	On	3							
		Elbow 2.00 Br40 I3 ax L																																			
Elbow DE Gout [factory] (Adult)	>Indications: visualization, characterization, quantification of monosodium urate crystal deposits, gout, gouty arthritis of the elbow joint, wo CM >Scan modes: RoutineSpiralAdultExtremitiesTwinspiral (max. 100 cm) >Decision Tree: Side [factory]	Topogram			ScanplanningTopoAdult		6 x 0.70				110		15	SPR (32 cm) 0.04	0.98	Standard					3.0																
		Topogram 0.70 Tr20 cor																													Tr20						
		Elbow		Extremities	RoutineSpiralAdultExtremitiesTwinspiral (max. 100 cm)	Side [factory]	32 x 0.70	On	75	Non-Contrast			(32 cm) 4.91	63.2	Standard		0.80	0.80		3.0			2.00	2.00	Side [factory].recon both; Side [factory].recon right	None	Orientation Multi Recons (3)	MPR	Br60	On	3						
		Elbow 2.00 Br60 I3 sag R																				2.00	2.00	Side [factory].recon both; Side [factory].recon right	None	Orientation Multi Recons (3)	MPR	Br60	On	3							
		Elbow 2.00 Br60 I3 cor R																																			
		Elbow 2.00 Br60 I3 ax R																				2.00	2.00	Side [factory].recon both; Side [factory].recon right	None	None	MPR	Br40	On	3							
		Elbow 2.00 Br40 I3 ax R																				0.80	0.50	Side [factory].recon both; Side [factory].recon right	DE Parallel Ranges	None	MPR	Qr40	On	3							
		Elbow 0.80 Qr40 I3 cor R																				2.00	2.00	Side [factory].recon both; Side [factory].recon left	None	Orientation Multi Recons (3)	MPR	Br60	On	3							
		Elbow 2.00 Br60 I3 sag L																																			
		Elbow 2.00 Br60 I3 cor L																																			
Foot [factory] (Adult)	>Indications: trauma, disorders of the foot bones, wo CM >Scan modes: RoutineSpiralAdultExtremitiesSn >Decision Tree: Side [factory]	Topogram			ScanplanningTopoAdult		6 x 0.70				110		15	SPR (32 cm) 0.04	0.98	Standard					3.0																
		Topogram 0.70 Tr20 cor																													Tr20						
		Foot		Extremities	RoutineSpiralAdultExtremitiesSn	Side [factory]	32 x 0.70	On	75	Auto (Bone/Calcium)	Sn 110	366	(32 cm) 5.22	60.2	Standard		0.40	1.00		3.0																	
		Foot 3.00 Br60 I3 sag R																				3.00	3.00	Side [factory].recon both; Side [factory].recon right	None	Orientation Multi Recons (3)	MPR	Br60	On	3							
		Foot 3.00 Br60 I3 cor R																																			
		Foot 3.00 Br60 I3 ax R																																			
		Foot 3.00 Br40 I3 ax R																				3.00	3.00	Side [factory].recon both; Side [factory].recon right	None	None	MPR	Br40	On	3							
		Foot 3.00 Br60 I3 sag L																				3.00	3.00	Side [factory].recon both; Side [factory].recon left	None	Orientation Multi Recons (3)	MPR	Br60	On	3							
		Foot 3.00 Br60 I3 cor L																																			
		Foot 3.00 Br60 I3 ax L																																			
Foot 3.00 Br40 I3 ax L																				3.00	3.00	Side [factory].recon both; Side [factory].recon left	None	None	MPR	Br40	On	3									

Protocol Name	Protocol Description	Range Name	Recon Name	Organ Characteristic	Scan Mode	Decision Trees	Acquisition [mm]	CARE Dose4D	FAST KV IQ level	FAST KV optimized for	kV	Eff. mAs	mA	CTDIvol [mGy]	DLP [mGy*cm]	Shaped Filter	Feed [mm]	Pitch	Rot. Time [s]	Scan Time [s]	Absolute Start Delay [s]	Contrast Protocol	Slice Thickness [mm]	Increment [mm]	Recon Strategies	Inline Results	Multi Recons	Image Type	Kernel	Iterative Recon	IR Strength	Recon Type	Recon Phase Step	Recon Phase(s)	Recon Phase Unit			
Foot DE Gout [factory] (Adult)	>Indications: visualization, characterization, quantification of monosodium urate crystal deposits, gout, gouty arthritis of the foot bones, wo CM >Scan modes: RoutineSpiralAdultExtremitiesTwinspiral (max. 100 cm) >Decision Tree: Side [factory]	Topogram			ScanplanningTopoAdult		6 x 0.70				110		15	SPR (32 cm) 0.04	0.98	Standard						3.0																
			Topogram 0.70 Tr20 cor																											Tr20								
		Foot		Extremities	RoutineSpiralAdultExtremitiesTwinspiral (max. 100 cm)	Side [factory]	32 x 0.70	On	75	Non-Contrast					(32 cm) 4.91	63.2	Standard		0.80	0.80		3.0																
			Foot 3.00 Br60 I3 sag R																					3.00	3.00		Side [factory].recon both; Side [factory].recon right	None	Orientation Multi Recons (3)	MPR	Br60	On	3					
			Foot 3.00 Br60 I3 cor R																																			
			Foot 3.00 Br60 I3 ax R																																			
			Foot 3.00 Br40 I3 ax R																					3.00	3.00		Side [factory].recon both; Side [factory].recon right	None	None	MPR	Br40	On	3					
			Foot 0.80 Qr40 I3 sag R																					0.80	0.50		Side [factory].recon both; Side [factory].recon right	DE Parallel Ranges	None	MPR	Qr40	On	3					
			Foot 3.00 Br60 I3 sag L																					3.00	3.00		Side [factory].recon both; Side [factory].recon left	None	Orientation Multi Recons (3)	MPR	Br60	On	3					
			Foot 3.00 Br60 I3 cor L																																			
			Foot 3.00 Br60 I3 ax L																																			
			Foot 3.00 Br40 I3 ax L																					3.00	3.00		Side [factory].recon both; Side [factory].recon left	None	None	MPR	Br40	On	3					
			Foot 0.80 Qr40 I3 sag L																					0.80	0.50		Side [factory].recon both; Side [factory].recon left	DE Parallel Ranges	None	MPR	Qr40	On	3					
		Leg or Ankle [factory] (Adult)	>Indications: trauma, disorders of the femur, fibula, tibia bones, ankle joint, wo CM >Scan modes: RoutineSpiralAdultExtremitiesSn >Decision Tree: Side [factory]	Topogram			ScanplanningTopoAdult		6 x 0.70				130		30	SPR (32 cm) 0.12	3.01	Standard						3.0														
	Topogram 0.70 Tr20 sag																													Tr20								
Leg or Ankle				Extremities	RoutineSpiralAdultExtremitiesSn	Side [factory]	32 x 0.70	On	75	Auto (Bone/Calcium)	Sn 110	366			(32 cm) 5.22	60.2	Standard		0.40	1.00		3.0																
	Leg or Ankle 2.00 Br60 I3 sag R																							2.00	2.00		Side [factory].recon both; Side [factory].recon right	None	Orientation Multi Recons (3)	MPR	Br60	On	3					
	Leg or Ankle 2.00 Br60 I3 cor R																																					
	Leg or Ankle 2.00 Br60 I3 ax R																																					
	Leg or Ankle 2.00 Br40 I3 ax R																							2.00	2.00		Side [factory].recon both; Side [factory].recon right	None	None	MPR	Br40	On	3					
	Leg or Ankle 2.00 Br60 I3 sag L																							2.00	2.00		Side [factory].recon both; Side [factory].recon left	None	Orientation Multi Recons (3)	MPR	Br60	On	3					
	Leg or Ankle 2.00 Br60 I3 cor L																																					
	Leg or Ankle 2.00 Br60 I3 ax																																					
	Leg or Ankle 2.00 Br40 I3 ax L																								2.00	2.00		Side [factory].recon both; Side [factory].recon left	None	None	MPR	Br40	On	3				

Protocol Name	Protocol Description	Range Name	Recon Name	Organ Characteristic	Scan Mode	Decision Trees	Acquisition [mm]	CARE Dose4D	FAST KV IQ level	FAST KV optimized for	kV	Eff. mAs	mA	CTDIvol [mGy]	DLP [mGy*cm]	Shaped Filter	Feed [mm]	Pitch	Rot. Time [s]	Scan Time [s]	Absolute Start Delay [s]	Contrast Protocol	Slice Thickness [mm]	Increment [mm]	Recon Strategies			Inline Results	Multi Recons	Image Type	Kernel	Iterative Recon	IR Strength	Recon Type	Recon Phase Step	Recon Phase(s)	Recon Phase Unit			
Brain DSA [factory] (Adult)	>Indications: stenosis, aneurysm, occlusion, plaques, abnormalities of the vertebral arteries, COW, w and wo CM >Scan modes: RoutineSpiralAdultHead, PremonitoringSequenceAdult2d, MonitoringSequenceAdult2d, RoutineSpiralAdultAngioHead >Contrast Protocol: CTA [factory]																																							
		Topogram			ScanplanningTopoAdult		6 x 0.70				130		30	SPR (16 cm) 0.25	6.23	Standard						3.0																		
			Topogram 0.70 Tr20 sag																														Tr20							
		Non Contrast Brain		Brain	RoutineSpiralAdultHead		32 x 0.70	On	36	Auto (Non-Contrast)	80	72	(16 cm) 4.41	62.0	Standard		0.55	1.00		3.0					1.50	1.00		None	None		Hv36	On	3							
		Premonitoring		Non Contrast Brain 1.50 Hv36 I3																																				
				Vascular Body			1 x 10.00	On	25	Soft Tissue with Contrast	110	22	(32 cm) 1.41	1.41	Standard	0.0		0.80	Full 0.80	3.0																				
																																			</					

Protocol Name	Protocol Description	Range Name	Recon Name	Organ Characteristic	Scan Mode	Decision Trees	Acquisition [mm]	CARE Dose4D	FAST kV IQ level	FAST kV optimized for	kV	Eff. mAs	mA	CTDIvol [mGy]	DLP [mGy*cm]	Shaped Filter	Feed [mm]	Pitch	Rot. Time [s]	Scan Time [s]	Absolute Start Delay [s]	Contrast Protocol	Slice Thickness [mm]	Increment [mm]	Recon Strategies			Inline Results	Multi Recons	Image Type	Kernel	Iterative Recon	IR Strength	Recon Type	Recon Phase Step	Recon Phase(s)	Recon Phase Unit	
Brain and Neck DSA [factory] (Adult)	>Indications: stenosis, aneurysm, occlusion, plaques, abnormalities of the carotid and vertebral arteries, COW, w and wo CM >Scan modes: RoutineSpiralAdultAngioNeck, PremonitoringSequenceAdult2d, MonitoringSequenceAdult2d, RoutineSpiralAdultAngioNeck >Contrast Protocol: Brain and Neck CTA [factory]	Topogram			ScanplanningTopoAdult		6 x 0.70				110		15	SPR (32 cm) 0.03	1.60	Standard					7.0																	
			Topogram 0.70 Tr20 cor																													Tr20						
		Non Contrast Brain and Neck		Neck	RoutineSpiralAdultAngioNeck		32 x 0.70	On	46	Auto (Non-Contrast)	110	42		(32 cm) 2.94	129	Standard		1.20	0.80		7.0																	
			Non Contrast Brain and Neck 1.50 Bv36 I3																				1.50	1.00			None	None		Bv36	On	3						
		Premonitoring		Vascular Body			1 x 10.00	On	25	Soft Tissue with Contrast	110	22		(32 cm) 1.41	1.41	Standard	0.0		0.80	Full 0.80	3.0																	
			Premonitoring 10.00 Br36																					10.00					None		Br36							
		Monitoring		Vascular Body			1 x 10.00	On	25	Soft Tissue with Contrast	110	22		(32 cm) 42.3	42.3	Standard	0.0		0.80	Full 0.80	10.0	Brain and Neck CTA [factory]																
			Monitoring 10.00 Br36																					10.00						None		Br36						
		CTA Brain and Neck		Neck	RoutineSpiralAdultAngioNeck		32 x 0.70	On	145	Auto (Vascular)	110	111		(32 cm) 7.78	341	Standard		1.20	0.80			Brain and Neck CTA [factory]																
	CTA Brain and Neck 1.50 Bv36 I3																					1.50	1.00			None	None		Bv36	On	3							
Neck CTA [factory] (Adult)	>Indications: stenosis, occlusion, plaques, abnormalities of the carotid arteries, w CM >Scan modes: PremonitoringSequenceAdult2d, MonitoringSequenceAdult2d, RoutineSpiralAdultAngioNeck >Contrast Protocol: CTA [factory]	Topogram			ScanplanningTopoAdult		6 x 0.70				110		15	SPR (32 cm) 0.03	1.60	Standard					7.0																	
			Topogram 0.70 Tr20 cor																													Tr20						
		Premonitoring		Vascular Body			1 x 10.00	On	25	Soft Tissue with Contrast	110	22		(32 cm) 1.41	1.41	Standard	0.0		0.80	Full 0.80	3.0																	
			Premonitoring 10.00 Br36																					10.00						None		Br36						
		Monitoring		Vascular Body			1 x 10.00	On	25	Soft Tissue with Contrast	110	22		(32 cm) 42.3	42.3	Standard	0.0		0.80	Full 0.80	10.0	CTA [factory]																
			Monitoring 10.00 Br36																					10.00						None		Br36						
		CTA Neck		Neck	RoutineSpiralAdultAngioNeck		32 x 0.70	On	145	Auto (Vascular)	110	111		(32 cm) 7.78	263	Standard		1.20	0.80			CTA [factory]																
			CTA Neck 3.00 Bv36 I3 sag																					3.00	3.00			None	Orientation Multi Recons (3)	MIP Thin	Bv36	On	3					
			CTA Neck 3.00 Bv36 I3 cor																																			
			CTA Neck 3.00 Bv36 I3 ax																																			
			CTA Neck 3.00 Bv36 I3 cor																					3.00	3.00			Parallel Ranges	None	MIP Thin	Bv36	On	3					
			CTA Neck 3.00 Br40 I3 Parallel Vascular																					3.00	3.00			Parallel Vascular Ranges	None	MPR	Br40	On	3					
			CTA Neck 3.00 Br40 I3 Parallel Vascular																					3.00	3.00			Parallel Vascular Ranges	None	MPR	Br40	On	3					

Protocol Name	Protocol Description	Range Name	Recon Name	Organ Characteristic	Scan Mode	Decision Trees	Acquisition [mm]	CARE Dose4D	FAST KV IQ level	FAST KV optimized for	kV	Eff. mAs	CTDIvol [mGy]	DLP [mGy*cm]	Shaped Filter	Feed [mm]	Pitch	Rot. Time [s]	Scan Time [s]	Absolute Start Delay [s]	Contrast Protocol	Slice Thickness [mm]	Increment [mm]	Recon Strategies			Inline Results	Multi Recons	Image Type	Kernel	Iterative Recon	IR Strength	Recon Type	Recon Phase Step	Recon Phase(s)	Recon Phase Unit	
Neck DSA [factory] (Adult)	>Indications: stenosis, occlusion, plaques, abnormalities of the carotid arteries, w and wo CM >Scan modes: RoutineSpiralAdultAngioNeck, PremonitoringSequenceAdult2d, MonitoringSequenceAdult2d, RoutineSpiralAdultAngioNeck >Contrast Protocol: CTA [factory]	Topogram			ScanplanningTopoAdult		6 x 0.70				110	15	SPR (32 cm) 0.04	0.98	Standard					7.0																	
			Topogram 0.70 Tr20 cor																												Tr20						
		Non Contrast Neck	Neck	RoutineSpiralAdultAngioNeck		32 x 0.70	On	46	Auto (Non-Contrast)	110	42	(32 cm) 2.94	99.6	Standard		1.20	0.80		7.0																		
		Non Contrast Neck 1.50 Bv36 I3																				1.50	1.00			None	None			Bv36	On	3					
		Premonitoring		Vascular Body		1 x 10.00	On	25	Soft Tissue with Contrast	110	22	(32 cm) 1.41	1.41	Standard	0.0		0.80	Full 0.80	3.0																		
		Premonitoring 10.00 Br36																				10.00					None			Br36							
		Monitoring		Vascular Body		1 x 10.00	On	25	Soft Tissue with Contrast	110	22	(32 cm) 42.3	42.3	Standard	0.0		0.80	Full 0.80	10.0	CTA [factory]																	
		Monitoring 10.00 Br36																				10.00					None			Br36							
		CTA Neck	Neck	RoutineSpiralAdultAngioNeck		32 x 0.70	On	145	Auto (Vascular)	110	111	(32 cm) 7.78	263	Standard		1.20	0.80			CTA [factory]																	
CTA Neck 1.50 Bv36 I3																				1.50	1.00			None	None			Bv36	On	3							
Thorax CTA [factory] (Adult)	>Indications: vascular anomalies, aneurysms, dissection, ruptures of the thoracic aorta, w CM >Scan modes: PremonitoringSequenceAdult2d, MonitoringSequenceAdult2d, RoutineSpiralAdultAngioBody >Contrast Protocol: CTA [factory]	Topogram			ScanplanningTopoAdult		6 x 0.70				110	15	SPR (32 cm) 0.03	1.60	Standard					6.0																	
			Topogram 0.70 Tr20 cor																												Tr20						
		Premonitoring		Vascular Body		1 x 10.00	On	25	Soft Tissue with Contrast	110	22	(32 cm) 1.41	1.41	Standard	0.0		0.80	Full 0.80	3.0																		
		Premonitoring 10.00 Br36																				10.00					None			Br36							
		Monitoring		Vascular Body		1 x 10.00	On	25	Soft Tissue with Contrast	110	22	(32 cm) 42.3	42.3	Standard	0.0		0.80	Full 0.80	10.0	CTA [factory]																	
		Monitoring 10.00 Br36																				10.00					None			Br36							
		CTA Thorax	Vascular Body	RoutineSpiralAdultAngioBody		32 x 0.70	On	70	Auto (Vascular)	110	56	(32 cm) 3.92	133	Standard		1.50	0.80			CTA [factory]																	
			CTA Thorax 5.00 Bv36 I3 sag																			5.00	5.00				None	Orientation Multi Recons (3)	MIP Thin	Bv36	On	3					
			CTA Thorax 5.00 Bv36 I3 cor																																		
			CTA Thorax 5.00 Bv36 I3 ax																																		
			CTA Thorax 5.00 BI60 I3 ax																			5.00	5.00				None	Impression Multi Recons (2)	MIP Thin	BI60	On	3					
			CTA Thorax 5.00 Br40 I3 ax																												Br40	On	3				
			CTA Thorax 3.00 Bv36 I3 Radial Vascular																			3.00					Radial Vascular Ranges	None	MIP Thin	Bv36	On	3					

Protocol Name	Protocol Description	Range Name	Recon Name	Organ Characteristic	Scan Mode	Decision Trees	Acquisition [mm]	CARE Dose4D	FAST kV IQ level	FAST kV optimized for	kV	Eff. mAs	CTDIvol [mGy]	DLP [mGy*cm]	Shaped Filter	Feed [mm]	Pitch	Rot. Time [s]	Scan Time [s]	Absolute Start Delay [s]	Contrast Protocol	Slice Thickness [mm]	Increment [mm]	Recon Strategies			Inline Results	Multi Recons	Image Type	Kernel	Iterative Recon	IR Strength	Recon Type	Recon Phase Step	Recon Phase(s)	Recon Phase Unit	
Pulmonary CTA [factory] (Adult)	>Indications: pulmonary embolism, w CM >Scan modes: PremonitoringSequenceAdult2d, MonitoringSequenceAdult2d, RoutineSpiralAdultAngioBody >Contrast Protocol: Pulmonary CTA [factory]	Topogram			ScanplanningTopoAdult		6 x 0.70				110	15	SPR (32 cm) 0.03	1.60	Standard					4.0																	
			Topogram 0.70 Tr20 cor																												Tr20						
		Premonitoring		Vascular Body			1 x 10.00	On	25	Soft Tissue with Contrast	110	22	(32 cm) 1.41	1.41	Standard	0.0		0.80	Full 0.80	3.0																	
			Premonitoring 10.00 Br36																				10.00					None			Br36						
		Monitoring		Vascular Body			1 x 10.00	On	25	Soft Tissue with Contrast	110	22	(32 cm) 42.3	42.3	Standard	0.0		0.80	Full 0.80	10.0	Pulmonary CTA [factory]																
			Monitoring 10.00 Br36																				10.00						None			Br36					
		CTA Pulmonary		Vascular Body	RoutineSpiralAdultAngioBody		32 x 0.70	On	70	Auto (Vascular)	110	56	(32 cm) 3.92	133	Standard			1.50	0.80			Pulmonary CTA [factory]															
			CTA Pulmonary 3.00 Bv36 I3 sag																				3.00	3.00			None	Orientation Multi Recons (3)	MIP Thin	Bv36	On	3					
			CTA Pulmonary 3.00 Bv36 I3 cor																																		
			CTA Pulmonary 3.00 Bv36 I3 ax																																		
			CTA Pulmonary 3.00 BI60 I3 ax																				3.00	3.00			None	Impression Multi Recons (2)	MIP Thin	BI60	On	3					
			CTA Pulmonary 3.00 Br40 I3 ax																											Br40	On	3					
Coronary Calcium Scoring [factory] (Adult)	>Indications: identification, quantification of calcified lesions in the coronary arteries, wo CM >Possible application: CT CaScoring >Scan modes: CardiacSpiralAdult	Topogram			ScanplanningTopoAdult		6 x 0.70				110	15	SPR (32 cm) 0.03	1.60	Standard					10.0																	
			Topogram 0.70 Tr20 cor																												Tr20						
		Calcium Scoring		Cardio	CardiacSpiralAdult		32 x 0.70	On	19	Auto (Bone/Calcium)	130	12	(32 cm) 8.19	147	Standard			0.31	0.80	10.0																	
			Calcium Scoring 3.00 Sa36 BestDiast ax																				3.00	1.50			None	None	MPR	Sa36	Off	3	BestPhase (diastolic)		75	%	
Abdomen CTA [factory] (Adult)	>Indications: vascular anomalies, aneurysms, dissection, ruptures of the abdominal aorta, renal arteries, w CM >Scan modes: PremonitoringSequenceAdult2d, MonitoringSequenceAdult2d, RoutineSpiralAdultAngioBody >Contrast Protocol: CTA [factory]	Topogram			ScanplanningTopoAdult		6 x 0.70				110	15	SPR (32 cm) 0.03	1.60	Standard					6.0																	
			Topogram 0.70 Tr20 cor																												Tr20						
		Premonitoring		Vascular Body			1 x 10.00	On	25	Soft Tissue with Contrast	110	22	(32 cm) 1.41	1.41	Standard	0.0		0.80	Full 0.80	3.0																	
			Premonitoring 10.00 Br36																				10.00					None			Br36						
		Monitoring		Vascular Body			1 x 10.00	On	25	Soft Tissue with Contrast	110	22	(32 cm) 42.3	42.3	Standard	0.0		0.80	Full 0.80	10.0	CTA [factory]																
			Monitoring 10.00 Br36																				10.00					None			Br36						
		CTA Abdomen		Vascular Body	RoutineSpiralAdultAngioBody		32 x 0.70	On	115	Auto (Vascular)	110	92	(32 cm) 6.44	283	Standard			1.50	0.80			CTA [factory]															
			CTA Abdomen 5.00 Bv36 I3 sag																				5.00	5.00			None	Orientation Multi Recons (3)	MIP Thin	Bv36	On	3					
			CTA Abdomen 5.00 Bv36 I3 cor																																		
			CTA Abdomen 5.00 Bv36 I3 ax																																		

Protocol Name	Protocol Description	Range Name	Recon Name	Organ Characteristic	Scan Mode	Decision Trees	Acquisition [mm]	CARE Dose4D	FAST kV IQ level	FAST kV optimized for	kV	Eff. mAs	mA	CTDIvol [mGy]	DLP [mGy*cm]	Shaped Filter	Feed [mm]	Pitch	Rot. Time [s]	Scan Time [s]	Absolute Start Delay [s]	Contrast Protocol	Slice Thickness [mm]	Increment [mm]	Recon Strategies			Inline Results	Multi Recons	Image Type	Kernel	Iterative Recon	IR Strength	Recon Type	Recon Phase Step	Recon Phase(s)	Recon Phase Unit			
Whole Body CTA [factory] (Adult)	>Indications: vascular anomalies, aneurysms, dissection, ruptures of the whole aorta, w CM >Scan modes: PremonitoringSequenceAdult2d, MonitoringSequenceAdult2d, RoutineSpiralAdultAngioBody >Contrast Protocol: CTA [factory]	Topogram			ScanplanningTopoAdult		6 x 0.70				110		15	SPR (32 cm) 0.03	2.83	Standard					6.0																			
			Topogram 0.70 Tr20 cor																																					
		Premonitoring		Vascular Body			1 x 10.00	On	25	Soft Tissue with Contrast	110	22		(32 cm) 1.41	1.41	Standard	0.0		0.80	Full 0.80	3.0																			
			Premonitoring 10.00 Br36												(32 cm) 42.3	42.3	Standard	0.0		0.80	Full 0.80	10.0	CTA [factory]	10.00					None		Br36									
		Monitoring		Vascular Body			1 x 10.00	On	25	Soft Tissue with Contrast	110	22		(32 cm) 42.3	42.3	Standard	0.0		0.80	Full 0.80	10.0	CTA [factory]																		
			Monitoring 10.00 Br36																					10.00							None		Br36							
		CTA Whole Body		Vascular Body	RoutineSpiralAdultAngioBody		32 x 0.70	On	115	Auto (Vascular)	110	92		(32 cm) 6.44	412	Standard		1.50	0.80				CTA [factory]																	
			CTA Whole Body 5.00 Bv36 I3 sag																						5.00	5.00				None	Orientation Multi Recons (3)	MIP Thin	Bv36	On	3					
			CTA Whole Body 5.00 Bv36 I3 cor																																					
			CTA Whole Body 5.00 Bv36 I3 ax																																					
			CTA Whole Body 3.00 Bv36 I3 Radial Vascular																						3.00						Radial Vascular Ranges	None	MIP Thin	Bv36	On	3				
		Extremity CTA [factory] (Adult)	>Indications: vascular anomalies, aneurysms of the femoral arteries, w CM >Scan modes: PremonitoringSequenceAdult2d, MonitoringSequenceAdult2d, RoutineSpiralAdultAngioExtremities >Contrast Protocol: CTA [factory]	Topogram			ScanplanningTopoAdult		6 x 0.70				110		15	SPR (32 cm) 0.03	2.83	Standard					3.0																	
	Topogram 0.70 Tr20 cor																																							
Premonitoring				Vascular Body			1 x 10.00	On	25	Soft Tissue with Contrast	110	22		(32 cm) 1.41	1.41	Standard	0.0		0.80	Full 0.80	3.0																			
	Premonitoring 10.00 Br36														(32 cm) 42.3	42.3	Standard	0.0		0.80	Full 0.80	10.0	CTA [factory]	10.00					None		Br36									
Monitoring				Vascular Body			1 x 10.00	On	25	Soft Tissue with Contrast	110	22		(32 cm) 42.3	42.3	Standard	0.0		0.80	Full 0.80	10.0	CTA [factory]																		
	Monitoring 10.00 Br36																							10.00							None		Br36							
CTA Extremity				Runoff	RoutineSpiralAdultAngioExtremities		32 x 0.70	On	155	Auto (Vascular)	110	124		(32 cm) 8.69	726	Standard		1.00	0.80				CTA [factory]																	
	CTA Extremity 5.00 Bv36 I3 sag																								5.00	5.00				None	Orientation Multi Recons (3)	MIP Thin	Bv36	On	3					
	CTA Extremity 5.00 Bv36 I3 cor																																							
	CTA Extremity 5.00 Bv36 I3 ax																																							
	CTA Extremity 3.00 Bv36 I3 Radial Vascular R																								3.00					Radial Vascular Ranges	None	MIP Thin	Bv36	On	3					
	CTA Extremity 3.00 Bv36 I3 Radial Vascular L																									3.00				Radial Vascular Ranges	None	MIP Thin	Bv36	On	3					

Protocol Name	Protocol Description	Range Name	Recon Name	Organ Characteristic	Scan Mode	Decision Trees	Acquisition [mm]	CARE Dose4D	FAST kV IQ level	FAST kV optimized for	kV	Eff. mAs	mAs	CTDIvol [mGy]	DLP [mGy*cm]	Shaped Filter	Feed [mm]	Pitch	Rot. Time [s]	Scan Time [s]	Absolute Start Delay [s]	Contrast Protocol	Slice Thickness [mm]	Increment [mm]	Recon Strategies			Inline Results	Multi Recons	Image Type	Kernel	Iterative Recon	IR Strength	Recon Type	Recon Phase Step	Recon Phase(s)	Recon Phase Unit		
Head Trauma [factory] (Adult)	>Indications: emergency head trauma, hemorrhage, wo CM >Scan modes: RoutineSpiralAdultHead																																						
		Topogram			ScanplanningTopoAdult		6 x 0.70				110		15	SPR (16 cm) 0.08	2.10	Standard						3.0																	
			Topogram 0.70 Tr20 cor																														Tr20						
		Head		Brain	RoutineSpiralAdultHead		32 x 0.70	Off			130	329		(16 cm) 73.0	1025	Standard		0.55	1.00		3.0																		
			Head 5.00 Hr40 ax																				5.00	5.00				None	Orientation Multi Recons (3)	MPR	Hr40	Off	2						
			Head 5.00 Hr40 sag																																				
			Head 5.00 Hr40 cor																																				
			Head 4.00 Hr40 BrainHem																				4.00	4.00				Brain Hemorrhage	None			Off	2						
			Head 1.50 Hr40 SkullUnf																				1.50	0.60				Skull Unfolding	None			Hr40	Off	3					
			Head 0.80 Hr60 ax																				0.80	0.60				None	Orientation Multi Recons (3)	MPR	Hr60	Off	2						
	Head 0.80 Hr60 sag																																						
	Head 0.80 Hr60 cor																																						
	Head 0.80 Hr60 SkullUnf																						0.80	0.60				Skull Unfolding	None			Hr60	Off	2					
Body Trauma [factory] (Adult)	>Indications: fast screening for emergency whole body trauma, hemorrhage, w/wo CM >Scan modes: RoutineSpiralAdultAbdomen																																						
		Topogram			ScanplanningTopoAdult		6 x 0.70				110		15	SPR (32 cm) 0.03	2.83	Standard						6.0																	
			Topogram 0.70 Tr20 cor																														Tr20						
		Body		Abdomen	RoutineSpiralAdultAbdomen		32 x 0.70	On	207	Auto (Non-Contrast)	130	131		(32 cm) 14.1	874	Standard		0.60	0.80		6.0																		
			Body 1.50 Br40 ax																				1.50	1.00				None	Impression Multi Recons (3)	MPR	Br40	Off	3						
			Body 1.50 Br60 ax																												Br60	Off	3						
			Body 1.50 BI60 ax																												BI60	Off	3						
			Body Br60 Rib																									Radial Rib Ranges	None	MPR	Br60	Off	3						

Protocol Name	Protocol Description	Range Name	Recon Name	Organ Characteristic	Scan Mode	Decision Trees	Acquisition [mm]	CARE Dose4D	FAST KV IQ level	FAST KV optimized for	kV	Eff. mAs	CTDIvol [mGy]	DLP [mGy*cm]	Shaped Filter	Feed [mm]	Pitch	Rot. Time [s]	Scan Time [s]	Absolute Start Delay [s]	Contrast Protocol	Slice Thickness [mm]	Increment [mm]	Recon Strategies			Inline Results	Multi Recons	Image Type	Kernel	Iterative Recon	IR Strength	Recon Type	Recon Phase Step	Recon Phase(s)	Recon Phase Unit		
Poly Trauma [factory] (Adult)	>Indications: emergency whole body trauma, hemorrhage, w/wo CMT >Scan modes: RoutineSpiralAdultHead, RoutineSpiralAdultNeck, RoutineSpiralAdultAbdomen	Topogram			ScanplanningTopoAdult		6 x 0.70				130	30	SPR (16 cm) 0.20	10.1	Standard					3.0																		
		Topogram 0.70 Tr20 sag																																				
		Head		Brain	RoutineSpiralAdultHead		32 x 0.70	On	429	Auto (Non-Contrast)	130	287		(16 cm) 63.7	767	Standard		0.55	1.00		3.0																	
			Head 5.00 Hr40 ax																				5.00	5.00			None	Orientation Multi Recons (3)	MPR	Hr40	Off	3						
			Head 5.00 Hr40 sag																																			
			Head 5.00 Hr40 cor																																			
			Head 4.00 Hr40 BrainHem																					4.00	4.00			Brain Hemorrhage	None			Off	3					
			Head 1.50 Hr40 SkullUnf																					1.50	0.60			Skull Unfolding	None		Hr40	Off	3					
			Head 5.00 Hr60																					5.00	5.00			None	None		Hr60	Off	3					
			Head 0.80 Hr60 SkullUnf																					0.80	0.60			Skull Unfolding	None		Hr60	Off	3					
		Neck		Neck	RoutineSpiralAdultNeck		32 x 0.70	On	222	Auto (Non-Contrast)	130	147		(32 cm) 15.8	298	Standard		0.80	1.00		17.0																	
			Neck 5.00 Br40 ax																					5.00	5.00			None	Orientation Multi Recons (3)	MPR	Br40	Off	3					
			Neck 5.00 Br40 sag																																			
			Neck 5.00 Br40 cor																																			
			Neck 5.00 Br60 sag																					5.00	5.00			None	None	MPR	Br60	Off	3					
		Topogram			ScanplanningTopoAdult		6 x 0.70					110	15		SPR (32 cm) 0.03	2.83	Standard					6.0																
			Topogram 0.70 Tr20 cor																																			
		Thorax and Abdomen		Abdomen	RoutineSpiralAdultAbdomen		32 x 0.70	On	207	Auto (Non-Contrast)	130	131		(32 cm) 14.1	875	Standard		0.60	0.80		6.0																	
			Thorax and Abdomen 5.00 Br40 ax																					5.00	5.00			None	Orientation Multi Recons (2)	MPR	Br40	Off	3					
			Thorax and Abdomen 5.00 Br40 cor																																			
			Thorax and Abdomen 5.00 BI60 sag																					5.00	5.00			None	Orientation Multi Recons (3)	MIP Thin	BI60	Off	3					
			Thorax and Abdomen 5.00 BI60 cor																																			
			Thorax and Abdomen 5.00 BI60 ax																																			
			Thorax and Abdomen 5.00 Br60 sag																					5.00	5.00			None	Orientation Multi Recons (2)	MPR	Br60	Off	3					
			Thorax and Abdomen 5.00 Br60 cor																																			
			Thorax and Abdomen Br60 Rib																																			

Protocol Name	Protocol Description	Range Name	Recon Name	Organ Characteristic	Scan Mode	Decision Trees	Acquisition [mm]	CARE Dose4D	FAST kV IQ level	FAST kV optimized for	kV	Eff. mAs	CTDIvol [mGy]	DLP [mGy*cm]	Shaped Filter	Feed [mm]	Pitch	Rot. Time [s]	Scan Time [s]	Absolute Start Delay [s]	Contrast Protocol	Slice Thickness [mm]	Increment [mm]	Recon Strategies			Inline Results	Multi Recons	Image Type	Kernel	Iterative Recon	IR Strength	Recon Type	Recon Phase Step	Recon Phase(s)	Recon Phase Unit	
RTP Head [factory] (Adult)	>Indications: radiotherapy planning of the head region, w/wvo CM >Scan modes: RoutineSpiralAdultHead >Decision Tree: RTP Brain [factory]	Topogram			ScanplanningTopoAdult		6 x 0.70				130	30	SPR (16 cm) 0.20	10.1	Standard					3.0																	
			Topogram 0.70 Tr20 sag																												Tr20						
		Head		Brain	RoutineSpiralAdultHead	RTP Brain [factory]	32 x 0.70	Off			130	230	(16 cm) 51.0	870	Standard		0.55	1.00		3.0																	
			Head 2.00 Qr40 I2 Head 1.00 Qr40 I2																			2.00	2.00	RTP Brain [factory].standard			None	None		Qr40	On	2					
																						1.00	1.00	RTP Brain [factory].stereotactic			None	None		Qr40	On	2					
RTP Head and Neck [factory] (Adult)	>Indications: radiotherapy planning of the head, neck ,shoulder region, w/wvo CM >Scan modes: RoutineSpiralAdultNeck	Topogram			ScanplanningTopoAdult		6 x 0.70				110	15	SPR (32 cm) 0.03	1.60	Standard					3.0																	
			Topogram 0.70 Tr20 cor																												Tr20						
		Head and Neck		Neck	RoutineSpiralAdultNeck		32 x 0.70	On	155	Auto (Non-Contrast)	130	102	(32 cm) 10.9	360	Standard		0.80	1.00		3.0																	
			Head and Neck 1.50 Qr40 I3																			1.50	1.50					None	None		Qr40	On	3				
RTP Thorax [factory] (Adult)	>Indications: radiotherapy planning of the thorax region, w/wvo CM >Scan modes: RoutineSpiralAdultThorax	Topogram			ScanplanningTopoAdult		6 x 0.70				110	15	SPR (32 cm) 0.03	1.60	Standard					3.0																	
			Topogram 0.70 Tr20 cor																												Tr20						
		Thorax		Thorax	RoutineSpiralAdultThorax		32 x 0.70	On	80	Auto (Non-Contrast)	130	51	(32 cm) 5.47	180	Standard		0.80	0.80		3.0																	
			Thorax 1.50 Qr40 I3																			1.50	1.50					None	None		Qr40	On	3				
RTP Respiratory Motion Management [factory] (Adult)	>Indications: motion management for radiotherapy planning for moving tumors, w/wvo CM >Scan modes: RespiratorySpiralAdult	Topogram			ScanplanningTopoAdult		6 x 0.70				110	15	SPR (32 cm) 0.03	1.60	Standard					3.0																	
			Topogram 0.70 Tr20 cor																												Tr20						
		Respiratory Motion		Respiratory	RespiratorySpiralAdult		32 x 0.70	On	26	Auto (Non-Contrast)	130	17	(32 cm) 30.0	978	Standard		0.12	0.80		3.0																	
			Respiratory Motion 1.50 Qr40 I3 Average CT																			1.50	1.50					None	None		Qr40	On	3	Average CT			
			Respiratory Motion 1.50 Qr40 I3 100% Ex - 80% In																			1.50	1.50					None	None		Qr40	On	3	Multi Phase	10	100% Ex; 80% Ex; 60% Ex; 40% Ex; 20% Ex; 0% Ex; 20% In; 40% In; 60% In; 80% In	
RTP Abdomen [factory] (Adult)	>Indications: radiotherapy planning of the abdomen region, w/wvo CM >Scan modes: RoutineSpiralAdultAbdomen	Topogram			ScanplanningTopoAdult		6 x 0.70				110	15	SPR (32 cm) 0.03	1.60	Standard					3.0																	
			Topogram 0.70 Tr20 cor																												Tr20						
		Abdomen		Abdomen	RoutineSpiralAdultAbdomen		32 x 0.70	On	145	Auto (Non-Contrast)	130	92	(32 cm) 9.87	325	Standard		0.80	0.80		3.0																	
			Abdomen 1.50 Qr40 I3																			1.50	1.50					None	None		Qr40	On	3				
RTP Pelvis [factory] (Adult)	>Indications: radiotherapy planning of the pelvis region, w/wvo CM >Scan modes: RoutineSpiralAdultPelvis	Topogram			ScanplanningTopoAdult		6 x 0.70				110	15	SPR (32 cm) 0.04	0.98	Standard					3.0																	
			Topogram 0.70 Tr20 cor																												Tr20						
		Pelvis		Pelvis	RoutineSpiralAdultPelvis		32 x 0.70	On	170	Auto (Non-Contrast)	130	107	(32 cm) 11.5	263	Standard		0.80	1.00		3.0																	
			Pelvis 1.50 Qr40 I3																			1.50	1.50					None	None		Qr40	On	3				

Protocol Name	Protocol Description	Range Name	Recon Name	Organ Characteristic	Scan Mode	Decision Trees	Acquisition [mm]	CARE Dose4D	FAST kV level	FAST kV optimized for	kV	Eff. mAs	CTDIvol [mGy]	DLP [mGy*cm]	Shaped Filter	Feed [mm]	Pitch	Rot. Time [s]	Scan Time [s]	Absolute Start Delay [s]	Contrast Protocol	Slice Thickness [mm]	Increment [mm]	Recon Strategies			Inline Results	Multi Recons	Image Type	Kernel	Iterative Recon	IR Strength	Recon Type	Recon Phase Step	Recon Phase(s)	Recon Phase Unit
RTP Extremity [factory] (Adult)	>Indications: radiotherapy planning of the extremities region, w/wo CM >Scan modes: RoutineSpiralAdultExtremities	Topogram			ScanplanningTopoAdult		6 x 0.70				110	15	SPR (32 cm) 0.04	0.98	Standard					3.0																
		Topogram 0.70 Tr20 cor																													Tr20					
		Extremity		Extremities	RoutineSpiralAdultExtremities		32 x 0.70	On	75	Auto (Bone/Calcium)	130	52	(32 cm) 5.58	71.9	Standard		0.80	1.00		3.0																
		Extremity 1.50 Qr40 I3																				1.50	1.50			None	None			Qr40	On	3				
RTP Direct Density commissioning [factory] (Adult)	>Indications: commissioning of Direct Density for radiotherapy planning scan protocols >Scan modes: RoutineSpiralAdultAbdomen	Topogram			ScanplanningTopoAdult		6 x 0.70				110	15	SPR (32 cm) 0.03	1.60	Standard					3.0																
		Topogram 0.70 Tr20 cor																													Tr20					
		80 kV		Abdomen	RoutineSpiralAdultAbdomen		32 x 0.70	Off			80	559	(32 cm) 15.0	469	Standard		0.35	0.80		3.0																
		80 kV 1.50 Qr40 I3																				1.50	1.00			None	None			Qr40	On	3				
		80 kV 1.50 Sd40 I3																				1.50	1.00			None	None			Sd40	On	3				
		80 kV 1.50 Sm40 I3																				1.50	1.00			None	None			Sm40	On	3				
		110 kV		Abdomen	RoutineSpiralAdultAbdomen		32 x 0.70	Off			110	214	(32 cm) 15.0	470	Standard		0.35	0.80		3.0																
		110 kV 1.50 Qr40 I3																				1.50	1.00			None	None			Qr40	On	3				
		110 kV 1.50 Sd40 I3																				1.50	1.00			None	None			Sd40	On	3				
		110 kV 1.50 Sm40 I3																				1.50	1.00			None	None			Sm40	On	3				
		130 kV		Abdomen	RoutineSpiralAdultAbdomen		32 x 0.70	Off			130	140	(32 cm) 15.0	471	Standard		0.35	0.80		3.0																
		130 kV 1.50 Qr40 I3																				1.50	1.00			None	None			Qr40	On	3				
		130 kV 1.50 Sd40 I3																				1.50	1.00			None	None			Sd40	On	3				
		130 kV 1.50 Sm40 I3																				1.50	1.00			None	None			Sm40	On	3				
CTDI measurement Head [factory] (Adult)	>Indications: special head protocol for ACR accreditation only; don't use for clinical cases! >Scan modes: RoutineSequenceAdultHead	Topogram			ScanplanningTopoAdult		6 x 0.70				110	15	SPR (16 cm) 0.08	2.10	Standard					3.0																
		Topogram 0.70 Tr20 cor																													Tr20					
		Head		Brain	RoutineSequenceAdultHead		32 x 0.70	Off			130	87	(16 cm) 21.6	303	Standard	20.0		1.00	Full 1.00	3.0																
		Head 5.00 Hr40																				5.00	5.00			None	None			Hr40	Off	3				
CTDI measurement Body [factory] (Adult)	>Indications: special body protocol for ACR accreditation only; don't use for clinical cases! >Scan modes: RoutineSequenceAdultAbdomen	Topogram			ScanplanningTopoAdult		6 x 0.70				110	15	SPR (32 cm) 0.04	0.98	Standard					3.0																
		Topogram 0.70 Tr20 cor																													Tr20					
		Body		Abdomen	RoutineSequenceAdultAbdomen		32 x 0.70	Off			130	82	(32 cm) 9.86	138	Standard	20.0		1.00	Full 1.00	3.0																
		Body 5.00 Br40																				5.00	5.00			None	None			Br40	Off	3				
Thorax Intervention 3D [factory] (Adult)	>Indications: biopsy, abscess in the thorax w/wo CM >Scan modes: InterventionSpiralAdultThorax, InterventionSequenceAdultThorax, InterventionFluoroAdultThorax2d	Topogram			ScanplanningTopoAdult		6 x 0.70				110	15	SPR (32 cm) 0.03	1.60	Standard					6.0																
		Topogram 0.70 Tr20 cor																													Tr20					
		i-Spiral Thorax		Thorax	InterventionSpiralAdultThorax		32 x 0.70	On	80	Auto (Non-Contrast)	130	51	(32 cm) 5.47	39.5	Standard		0.80	0.80		6.0																
		i-Spiral Thorax 5.00 BI60																				5.00	4.80			None	None			BI60	Off	3				
		i-Sequence Thorax		Thorax	InterventionSequenceAdultThorax		32 x 0.70	Off			130	17	(32 cm) 1.82	4.09	Standard	0.0		0.80	Full 0.80																	
		i-Sequence Thorax 4.80 Br40																				4.80	4.40			None	None			Br40	Off	3				
		i-Fluoro Thorax		Thorax	InterventionFluoroAdultThorax2d		3 x 5.00	Off			130	17			Standard			0.80																		
		i-Fluoro Thorax 5.00 Br40																				5.00	5.00				None			Br40						

Protocol Name	Protocol Description	Range Name	Recon Name	Organ Characteristic	Scan Mode	Decision Trees	Acquisition [mm]	CARE Dose4D	FAST kV IQ level	FAST kV optimized for	kV	Eff. mAs	CTDIvol [mGy]	DLP [mGy*cm]	Shaped Filter	Feed [mm]	Pitch	Rot. Time [s]	Scan Time [s]	Absolute Start Delay [s]	Contrast Protocol	Slice Thickness [mm]	Increment [mm]	Recon Strategies			Inline Results	Multi Recons	Image Type	Kernel	Iterative Recon	IR Strength	Recon Type	Recon Phase Step	Recon Phase(s)	Recon Phase Unit		
Thorax Intervention [factory] (Adult)	>Indications: biopsy, abscess in the thorax w/wo CM >Scan modes: RoutineSpiralAdultThorax, InterventionSequenceAdultThorax2d, InterventionFluoroAdultThorax2d, RoutineSpiralAdultThorax >Decision Tree: Thorax FoV [factory]																																					
		Topogram			ScanplanningTopoAdult		6 x 0.70				110	15	SPR (32 cm) 0.03	1.60	Standard						6.0																	
			Topogram 0.70 Tr20 cor																												Tr20							
		Planning Thorax		Thorax	RoutineSpiralAdultThorax	Thorax FoV [factory]	32 x 0.70	On	80	Auto (Non-Contrast)	130	51	(32 cm) 5.47	180	Standard		0.80	0.80		6.0																		
			Planning Thorax 5.00 BI60 I3																			5.00	5.00		Thorax FoV [factory].fast acquisition	None	None		BI60	On	3							
			Planning Thorax 5.00 Br40 I3																			5.00	5.00		Thorax FoV [factory].fast acquisition	None	None		Br40	On	3							
			Planning Thorax 5.00 BI60 I3																			5.00	5.00		Thorax FoV [factory].standard acquisition	None	None		BI60	On	3							
			Planning Thorax 5.00 Br40 I3																			5.00	5.00		Thorax FoV [factory].standard acquisition	None	None		Br40	On	3							
		i-Sequence Thorax		Thorax	InterventionSequenceAdultThorax2d		3 x 5.00	Off			130	17	(32 cm) 2.11	3.16	Standard	0.0		0.80	Full 0.80																			
			i-Sequence Thorax 5.00 Br40																			5.00						None		Br40								
		i-Fluoro Thorax		Thorax	InterventionFluoroAdultThorax2d		3 x 5.00	Off			130	17				Standard			0.80																			
			i-Fluoro Thorax 5.00 Br40																			5.00	5.00					None		Br40								
		Control Thorax		Thorax	RoutineSpiralAdultThorax	Thorax FoV [factory]	32 x 0.70	On	80	Auto (Non-Contrast)	130	51	(32 cm) 5.47	180	Standard		0.80	0.80		6.0																		
			Control Thorax 5.00 BI60 I3																			5.00	5.00		Thorax FoV [factory].fast acquisition	None	None		BI60	On	3							
			Control Thorax 5.00 Br40 I3																			5.00	5.00		Thorax FoV [factory].fast acquisition	None	None		Br40	On	3							
			Control Thorax 5.00 BI60 I3																			5.00	5.00		Thorax FoV [factory].standard acquisition	None	None		BI60	On	3							
			Control Thorax 5.00 Br40 I3																			5.00	5.00		Thorax FoV [factory].standard acquisition	None	None		Br40	On	3							
Abdomen Intervention 3D [factory] (Adult)	>Indications: biopsy, abscess in the abdomen, w/wo CM >Scan modes: InterventionSpiralAdultAbdomen, InterventionSequenceAdultAbdomen, InterventionFluoroAdultAbdomen2d	Topogram			ScanplanningTopoAdult		6 x 0.70				110	15	SPR (32 cm) 0.03	1.60	Standard					6.0																		
			Topogram 0.70 Tr20 cor																												Tr20							
		i-Spiral Abdomen		Abdomen	InterventionSpiralAdultAbdomen		32 x 0.70	On	145	Auto (Non-Contrast)	130	92	(32 cm) 9.87	71.3	Standard		0.80	0.80		6.0																		
			i-Spiral Abdomen 5.00 Br40																			5.00	4.80					None	None		Br40	Off	3					
		i-Sequence Abdomen		Abdomen	InterventionSequenceAdultAbdomen		32 x 0.70	Off			130	31	(32 cm) 3.33	7.45	Standard	0.0		0.80	Full 0.80																			
			i-Sequence Abdomen 4.80 Br40																			4.80	4.40					None	None		Br40	Off	3					
		i-Fluoro Abdomen		Abdomen	InterventionFluoroAdultAbdomen2d		3 x 5.00	Off			130	31			Standard			0.80																				
			i-Fluoro Abdomen 5.00 Br40																			5.00	5.00					None		Br40								

Protocol Name	Protocol Description	Range Name	Recon Name	Organ Characteristic	Scan Mode	Decision Trees	Acquisition [mm]	CARE Dose4D	FAST KV IQ level	FAST KV optimized for	kV	Eff. mAs	mA	CTDIvol [mGy]	DLP [mGy*cm]	Shaped Filter	Feed [mm]	Pitch	Rot. Time [s]	Scan Time [s]	Absolute Start Delay [s]	Contrast Protocol	Slice Thickness [mm]	Increment [mm]	Recon Strategies		Inline Results	Multi Recons	Image Type	Kernel	Iterative Recon	IR Strength	Recon Type	Recon Phase Step	Recon Phase(s)	Recon Phase Unit		
Abdomen Intervention [factory] (Adult)	>Indications: biopsy, abscess in the abdomen, w/wvo CM >Scan modes: RoutineSpiralAdultAbdomen, InterventionSequenceAdultAbdomen2d, InterventionFluoroAdultAbdomen2d	Topogram			ScanplanningTopoAdult		6 x 0.70				110		15	SPR (32 cm) 0.03	1.60	Standard						6.0																
		Topogram 0.70 Tr20 cor																												Tr20								
		Planning Abdomen		Abdomen	RoutineSpiralAdultAbdomen		32 x 0.70	On	145	Auto (Non-Contrast)	130	92		(32 cm) 9.87	325	Standard		0.80	0.80		6.0				5.00	5.00			None	None	Br40	On	3					
		Planning Abdomen 5.00 Br40 I3																											None		Br40							
		i-Sequence Abdomen		Abdomen	InterventionSequenceAdultAbdomen2d		3 x 5.00	Off			130	31		(32 cm) 3.85	5.77	Standard	0.0		0.80	Full 0.80																		
		i-Sequence Abdomen 5.00 Br40																											None		Br40							
		i-Fluoro Abdomen		Abdomen	InterventionFluoroAdultAbdomen2d		3 x 5.00	Off			130	31				Standard			0.80																			
		i-Fluoro Abdomen 5.00 Br40																								5.00	5.00				None		Br40					
Spine Pain Therapy 3D [factory] (Adult)	>Indications: periradicular therapy, nerve root block, facet joint injection, arthrogram in the spine region, w/wvo CM >Scan modes: InterventionSpiralAdultSpineSn, InterventionSequenceAdultSpineSn, InterventionFluoroAdultSpine2d	Topogram			ScanplanningTopoAdult		6 x 0.70				130		30	SPR (32 cm) 0.10	4.90	Standard						3.0																
		Topogram 0.70 Tr20 sag																												Tr20								
		i-Spiral Spine		Spine	InterventionSpiralAdultSpineSn		32 x 0.70	On	70	Auto (Bone/Calcium)	Sn 130	160		(32 cm) 4.59	29.5	Standard		0.60	0.80		6.0				3.00	2.80			None	None	Br60	Off	3					
		i-Spiral Spine 3.00 Br60																											None									
		i-Sequence Spine		Spine	InterventionSequenceAdultSpineSn		32 x 0.70	Off			Sn 130	55		(32 cm) 1.58	3.53	Standard	0.0		0.80	Full 0.80																		
		i-Sequence Spine 2.60 Br60																											None	None	Br60	Off	3					
		i-Fluoro Spine		Spine	InterventionFluoroAdultSpine2d		3 x 3.00	Off			130	15				Standard			0.80																			
		i-Fluoro Spine 3.00 Br60																								3.00	3.00				None		Br60					
Spine Pain Therapy [factory] (Adult)	>Indications: periradicular therapy, nerve root block, facet joint injection, arthrogram in the spine region, w/wvo CM >Scan modes: RoutineSpiralAdultSpineSn, InterventionSequenceAdultSpine2d, InterventionFluoroAdultSpine2d	Topogram			ScanplanningTopoAdult		6 x 0.70				130		30	SPR (32 cm) 0.10	4.90	Standard						3.0																
		Topogram 0.70 Tr20 sag																												Tr20								
		Planning Spine		Spine	RoutineSpiralAdultSpineSn		32 x 0.70	On	70	Auto (Bone/Calcium)	Sn 130	160		(32 cm) 4.59	83.5	Standard		0.60	0.80		6.0				3.00	3.00			None	None	Br60	On	3					
		Planning Spine 3.00 Br60 I3																											None	None								
		Planning Spine 3.00 Br40 I3																											None	None	Br40	On	3					
		i-Sequence Spine		Spine	InterventionSequenceAdultSpine2d		3 x 5.00	Off			130	15		(32 cm) 1.86	2.79	Standard	0.0		0.80	Full 0.80																		
		i-Sequence Spine 5.00 Br60																											None		Br60							
		i-Fluoro Spine		Spine	InterventionFluoroAdultSpine2d		3 x 3.00	Off			130	15				Standard			0.80																			
i-Fluoro Spine 3.00 Br60																								3.00	3.00				None		Br60							

Protocol Name	Protocol Description	Range Name	Recon Name	Organ Characteristic	Scan Mode	Decision Trees	Acquisition [mm]	CARE Dose4D	FAST kV IQ level	FAST kV optimized for	kV	Eff. mAs	mAs	CTDIvol [mGy]	DLP [mGy*cm]	Shaped Filter	Feed [mm]	Pitch	Rot. Time [s]	Scan Time [s]	Absolute Start Delay [s]	Contrast Protocol	Slice Thickness [mm]	Increment [mm]	Recon Strategies	Inline Results	Multi Recons	Image Type	Kernel	Iterative Recon	IR Strength	Recon Type	Recon Phase Step	Recon Phase(s)	Recon Phase Unit
Brain [factory] (Child)	>Indications: head trauma, hemorrhage, hydrocephalus, abnormalities, inflammation, wo CM >Scan modes: RoutineSpiralChildHeadXcare	Topogram			ScanplanningTopoChild		6 x 0.70				130		15	SPR (16 cm) 0.12	3.11	Standard					3.0														
			Topogram 0.70 Tr20 sag																										Tr20						
		Brain		Brain	RoutineSpiralChildHeadXcare		32 x 0.70	On	345	Auto (Non-Contrast)	110	314		(16 cm) 46.9	565	Standard		0.55	1.00		3.0														
			Brain 5.00 Hr40 I2 Brain 5.00 Hr60 I2																				5.00	5.00		None	Impression Multi Recons (2)		Hr40 Hr60	On	2				
Brain Multi Phase [factory] (Child)	>Indications: tumors, hydrocephalus, abnormalities, inflammation, w and wo CM >Scan modes: RoutineSpiralChildHeadXcare	Topogram			ScanplanningTopoChild		6 x 0.70				130		15	SPR (16 cm) 0.12	3.11	Standard					3.0														
			Topogram 0.70 Tr20 sag																										Tr20						
		Non Contrast Brain		Brain	RoutineSpiralChildHeadXcare		32 x 0.70	On	345	Auto (Non-Contrast)	110	314		(16 cm) 46.9	565	Standard		0.55	1.00		3.0														
			Non Contrast Brain 5.00 Hr40 I2 Non Contrast Brain 5.00 Hr60 I2																				5.00	5.00		None	Impression Multi Recons (2)		Hr40 Hr60	On	2				
		Venous Phase Brain		Brain	RoutineSpiralChildHeadXcare		32 x 0.70	On	345	Soft Tissue with Contrast	110	280		(16 cm) 41.8	504	Standard		0.55	1.00		100.0														
			Venous Phase Brain 5.00 Hr40 I2 Venous Phase Brain 5.00 Hr60 I2																				5.00	5.00		None	Impression Multi Recons (2)		Hr40 Hr60	On	2				
Inner Ear [factory] (Child)	>Indications: inflammation, tumorous processes of pyramids, acoustic neuroma, pathologies of the mastoid process, cerebellopontine angle tumors, cochlea implants, trauma, wo CM >Scan modes: RoutineSpiralChildHeadSn	Topogram			ScanplanningTopoChild		6 x 0.70				130		15	SPR (16 cm) 0.12	3.11	Standard					3.0														
			Topogram 0.70 Tr20 sag																										Tr20						
		Inner Ear		Brain	RoutineSpiralChildHeadSn		32 x 0.70	On	110	Auto (Bone/Calcium)	Sn 110	475		(16 cm) 12.9	91.1	Standard		0.55	1.00		3.0														
			Inner Ear 1.50 Hr60 I3 ax Inner Ear 1.50 Hr40 I3 ax Inner Ear 0.80 Hr60 I3 ax																				1.50	1.50		None	Impression Multi Recons (2)	MPR	Hr60 Hr40	On	3				
Sinus [factory] (Child)	>Indications: sinusitis, mucocoele, pneumatization, polypoid, tumor, corrections, wo CM >Scan modes: RoutineSpiralChildHeadSnXcare	Topogram			ScanplanningTopoChild		6 x 0.70				130		15	SPR (16 cm) 0.12	3.11	Standard					3.0														
			Topogram 0.70 Tr20 sag																										Tr20						
		Sinus		Brain	RoutineSpiralChildHeadSnXcare		32 x 0.70	On	50	Auto (Bone/Calcium)	Sn 110	236		(16 cm) 6.43	64.6	Standard		0.55	1.00		3.0														
			Sinus 3.00 Hr60 I3 ax Sinus 3.00 Hr40 I3 ax Sinus 3.00 Hr60 I3 ax																				3.00	3.00		None	Impression Multi Recons (2)	MPR	Hr60 Hr40	On	3				
Orbit [factory] (Child)	>Indications: fracture, tumor, trauma, wo CM >Scan modes: RoutineSpiralChildHeadSnXcare	Topogram			ScanplanningTopoChild		6 x 0.70				130		15	SPR (16 cm) 0.12	3.11	Standard					3.0														
			Topogram 0.70 Tr20 sag Topogram Recon 1																										Tr20						
		Orbit		Brain	RoutineSpiralChildHeadSnXcare		32 x 0.70	On	110	Auto (Bone/Calcium)	Sn 110	475		(16 cm) 12.9	91.1	Standard		0.55	1.00		3.0														
			Orbit 2.00 Hr60 I3 ax Orbit 2.00 Hr40 I3 ax Orbit 2.00 Hr60 I3 ax																				2.00	2.00		None	Impression Multi Recons (2)	MPR	Hr60 Hr40	On	3				
			Orbit 2.00 Hr60 I3 ax																				2.00	2.00		None	None	MPR	Hr60	On	3				
			Orbit Hr40 I3 cor Radial																							Radial Ranges	None	VRT	Hr40	On	3				

Protocol Name	Protocol Description	Range Name	Recon Name	Organ Characteristic	Scan Mode	Decision Trees	Acquisition [mm]	CARE Dose4D	FAST kV IQ level	FAST kV optimized for	kV	Eff. mAs	mA	CTDIvol [mGy]	DLP [mGy*cm]	Shaped Filter	Feed [mm]	Pitch	Rot. Time [s]	Scan Time [s]	Absolute Start Delay [s]	Contrast Protocol	Slice Thickness [mm]	Increment [mm]	Recon Strategies	Inline Results	Multi Recons	Image Type	Kernel	Iterative Recon	IR Strength	Recon Type	Recon Phase Step	Recon Phase(s)	Recon Phase Unit
Neck [factory] (Child)	>Indications: masses, tumor, metastases, lymphoma, lymph nodes in the larynx or pharynx region, w/w CM >Scan modes: RoutineSpiralChildNeck	Topogram			ScanplanningTopoChild		6 x 0.70				130	15		SPR (32 cm) 0.06	1.51	Standard					7.0														
		Topogram 0.70 Tr20 sag																											Tr20						
		Neck	Neck	RoutineSpiralChildNeck		32 x 0.70	On	155	Auto (Non-Contrast)	110	143		(32 cm) 10.0	189	Standard		1.20	1.00		7.0															
		Neck 5.00 Br40 I3 sag Neck 5.00 Br40 I3 cor Neck 5.00 Br40 I3 ax																					5.00	5.00		None	Orientation Multi Recons (3)	MPR	Br40	On	3				
Shoulder [factory] (Child)	>Indications: evaluation of joint cavities, trauma, dislocations, orthopedic indications, shoulder joint, wo CM >Scan modes: RoutineSpiralChildShoulder	Topogram			ScanplanningTopoChild		6 x 0.70				80	20		SPR (32 cm) 0.02	0.50	Standard					3.0														
		Topogram 0.70 Tr20 cor																											Tr20						
		Shoulder	Shoulder	RoutineSpiralChildShoulder		32 x 0.70	On	130	Auto (Bone/Calcium)	110	121		(32 cm) 8.48	144	Standard		1.50	1.00		3.0															
		Shoulder 5.00 Br60 I3 ax Shoulder 3.00 Br60 I3 sag R Shoulder 3.00 Br60 I3 cor R Shoulder 3.00 Br60 I3 ax R																					5.00	5.00		None	None	MPR	Br60	On	3				
		Shoulder 3.00 Br60 I3 sag R Shoulder 3.00 Br60 I3 cor R Shoulder 3.00 Br60 I3 ax R																					3.00	3.00		Parallel Ranges	Orientation Multi Recons (3)	MPR	Br60	On	3				
		Shoulder 3.00 Br40 I3 sag R Shoulder 3.00 Br40 I3 cor R Shoulder 3.00 Br40 I3 ax R																					3.00	3.00		Parallel Ranges	Orientation Multi Recons (3)	MPR	Br40	On	3				
		Shoulder 3.00 Br60 I3 sag L Shoulder 3.00 Br60 I3 cor L Shoulder 3.00 Br60 I3 ax L																					3.00	3.00		Parallel Ranges	Orientation Multi Recons (3)	MPR	Br60	On	3				
		Shoulder 3.00 Br40 I3 sag L Shoulder 3.00 Br40 I3 cor L Shoulder 3.00 Br40 I3 ax L																					3.00	3.00		Parallel Ranges	Orientation Multi Recons (3)	MPR	Br40	On	3				
Thorax [factory] (Child)	>Indications: masses, tumor, metastases, lymphoma, lymph nodes, tuberculosis, empyema, infection, early visualization of pulmonary nodules, screening of the lung, w/w CM >Scan modes: RoutineSpiralChildThorax	Topogram			ScanplanningTopoChild		6 x 0.70				80	20		SPR (32 cm) 0.02	0.81	Standard					6.0														
		Topogram 0.70 Tr20 cor																											Tr20						
		Thorax	Thorax	RoutineSpiralChildThorax		32 x 0.70	On	80	Auto (Non-Contrast)	110	76		(32 cm) 5.32	181	Standard		1.50	0.80		6.0															
		Thorax 5.00 Br40 I3 cor Thorax 5.00 Br40 I3 ax Thorax 5.00 BI60 I3 sag Thorax 5.00 BI60 I3 cor Thorax 5.00 BI60 I3 ax																					5.00	5.00		None	Orientation Multi Recons (2)	MPR	Br40	On	3				
		Thorax 5.00 BI60 I3 sag Thorax 5.00 BI60 I3 cor Thorax 5.00 BI60 I3 ax																					5.00	5.00		None	Orientation Multi Recons (3)	MIP Thin	BI60	On	3				
		Thorax 10.00 BI60 I3 cor Thorax 1.00 BI60 I3																					10.00	10.00		None	None	MIP Thin	BI60	On	3				
																							1.00	0.70		Lung CAD	None		BI60	On	3				

Protocol Name	Protocol Description	Range Name	Recon Name	Organ Characteristic	Scan Mode	Decision Trees	Acquisition [mm]	CARE Dose4D	FAST kV IQ level	FAST kV optimized for	kV	Eff. mAs	CTDIvol [mGy]	DLP [mGy*cm]	Shaped Filter	Feed [mm]	Pitch	Rot. Time [s]	Scan Time [s]	Absolute Start Delay [s]	Contrast Protocol	Slice Thickness [mm]	Increment [mm]	Recon Strategies	Inline Results	Multi Recons	Image Type	Kernel	Iterative Recon	IR Strength	Recon Type	Recon Phase Step	Recon Phase(s)	Recon Phase Unit
Thorax and Abdomen [factory] (Child)	>Indications: masses, tumor, lymphoma, lymph nodes of the lung, liver, pancreas, gall bladder, spleen, kidneys, w CM, 2 phases preconfigured >Scan modes: RoutineSpiralChildAbdomenFlexdose	Topogram			ScanplanningTopoChild		6 x 0.70				80	20	SPR (32 cm) 0.02	0.81	Standard					6.0														
		Topogram 0.70 Tr20 cor																										Tr20						
		Thorax and Abdomen		Abdomen	RoutineSpiralChildAbdomenFlexdose		32 x 0.70	On	145	Auto (Non-Contrast)	110	113	(32 cm) 7.90	425	Standard		1.20	0.80		6.0														
		Thorax and Abdomen 5.00 Br40 I3 cor																				5.00	5.00		None	Orientation Multi Recons (2)	MPR	Br40	On	3				
		Thorax and Abdomen 5.00 Br40 I3 ax																																
		Thorax and Abdomen 5.00 BI60 I3 sag Lung																				5.00	5.00		None	Orientation Multi Recons (3)	MIP Thin	BI60	On	3				
		Thorax and Abdomen 5.00 BI60 I3 cor Lung																																
		Thorax and Abdomen 5.00 BI60 I3 ax Lung																																
Abdomen [factory] (Child)	>Indications: masses, tumor, lymphoma, lymph nodes of the liver, pancreas, gall bladder, spleen, kidneys, intestine, w/w CM >Scan modes: RoutineSpiralChildAbdomen	Topogram			ScanplanningTopoChild		6 x 0.70				80	20	SPR (32 cm) 0.02	0.81	Standard					6.0														
		Topogram 0.70 Tr20 cor																										Tr20						
		Abdomen		Abdomen	RoutineSpiralChildAbdomen		32 x 0.70	On	145	Auto (Non-Contrast)	110	140	(32 cm) 9.81	332	Standard		1.20	0.80		6.0														
		Abdomen 5.00 Br40 I3 cor Abdomen 5.00 Br40 I3 ax																				5.00	5.00		None	Orientation Multi Recons (2)	MPR	Br40	On	3				
Abdomen Multi Phase [factory] (Child)	>Indications: masses, tumor, lymphoma, lymph nodes of the liver, pancreas, gall bladder, spleen, kidneys, intestine, w and w CM, 3 phases preconfigured >Scan modes: RoutineSpiralChildAbdomen, PremonitoringSequenceChild2d, MonitoringSequenceChild2d, RoutineSpiralChildAbdomen	Topogram			ScanplanningTopoChild		6 x 0.70				80	20	SPR (32 cm) 0.02	0.81	Standard					6.0														
		Topogram 0.70 Tr20 cor																										Tr20						
		Non Contrast Liver		Abdomen	RoutineSpiralChildAbdomen		32 x 0.70	On	145	Auto (Non-Contrast)	110	140	(32 cm) 9.81	234	Standard		1.20	0.80		6.0														
		Non Contrast Liver 5.00 Br40 I3 cor Non Contrast Liver 5.00 Br40 I3 ax																				5.00	5.00		None	Orientation Multi Recons (2)	MPR	Br40	On	3				
		Premonitoring		Vascular Body			1 x 5.00	On	25	Soft Tissue with Contrast	80	47	(32 cm) 1.15	0.58	Standard	0.0		0.80	Full 0.80	3.0														
		Premonitoring 5.00 Br36																				5.00				None		Br36						
		Monitoring		Vascular Body			1 x 5.00	On	25	Soft Tissue with Contrast	80	47	(32 cm) 34.6	17.3	Standard	0.0		0.80	Full 0.80	10.0														
		Monitoring 5.00 Br36																				5.00				None		Br36						
		Arterial Phase Abdomen		Abdomen	RoutineSpiralChildAbdomen		32 x 0.70	On	145	Auto (Soft Tissue with Contrast)	110	125	(32 cm) 8.76	297	Standard		1.20	0.80																
		Arterial Phase Abdomen 5.00 Br40 I3 cor Arterial Phase Abdomen 5.00 Br40 I3 ax																				5.00	5.00		None	Orientation Multi Recons (2)	MPR	Br40	On	3				
		Venous Phase Abdomen		Abdomen	RoutineSpiralChildAbdomen		32 x 0.70	On	145	Auto (Soft Tissue with Contrast)	110	125	(32 cm) 8.76	209	Standard		1.20	0.80		71.0														
		Venous Phase Abdomen 5.00 Br40 I3 cor Venous Phase Abdomen 5.00 Br40 I3 ax																				5.00	5.00		None	Orientation Multi Recons (2)	MPR	Br40	On	3				

Protocol Name	Protocol Description	Range Name	Recon Name	Organ Characteristic	Scan Mode	Decision Trees	Acquisition [mm]	CARE Dose4D	FAST kV IQ level	FAST kV optimized for	kV	Eff. mAs	CTDIvol [mGy]	DLP [mGy*cm]	Shaped Filter	Feed [mm]	Pitch	Rot. Time [s]	Scan Time [s]	Absolute Start Delay [s]	Contrast Protocol	Slice Thickness [mm]	Increment [mm]	Recon Strategies	Inline Results	Multi Recons	Image Type	Kernel	Iterative Recon	IR Strength	Recon Type	Recon Phase Step	Recon Phase(s)	Recon Phase Unit	
Whole Body [factory] (Child)	>Indications: masses, tumor, lymphoma, lymph nodes of the lung, liver, pancreas, gall bladder, spleen, kidneys, w/w CM >Scan modes: RoutineSpiralChildAbdomen																																		
		Topogram			ScanplanningTopoChild		6 x 0.70				80	20	SPR (32 cm) 0.02	0.81	Standard						6.0														
		Topogram 0.70 Tr20 cor																										Tr20							
		Whole Body		Abdomen	RoutineSpiralChildAbdomen		32 x 0.70	On	145	Auto (Non-Contrast)	110	140	(32 cm) 9.81	528	Standard		1.20	0.80		6.0															
			Whole Body 5.00 Br40 I3 cor																				5.00	5.00		None	Orientation Multi Recons (2)	MPR	Br40	On	3				
			Whole Body 5.00 Br40 I3 ax																																
			Whole Body 5.00 BI60 I3 sag Lung																				5.00	5.00		None	Orientation Multi Recons (3)	MIP Thin	BI60	On	3				
			Whole Body 5.00 BI60 I3 cor Lung																																
			Whole Body 5.00 BI60 I3 ax Lung																																
			Whole Body 1.00 BI60 I3																				1.00	0.70		Lung CAD	None		BI60	On	3				
			Whole Body 5.00 Br60 I3 sag																				5.00	5.00		None	Orientation Multi Recons (3)	MPR	Br60	On	3				
			Whole Body 5.00 Br60 I3 cor																																
			Whole Body 5.00 Br60 I3 ax																																
	Whole Body 5.00 Br40 I3 ax																				5.00	5.00		None	None	MPR	Br40	On	3						
Pelvis [factory] (Child)	>Indications: masses, tumor, lymphom of the urinary bladder, rectum, bone tumor, trauma, osteolysis of the pelvis, w/w CM >Scan modes: RoutineSpiralChildPelvis	Topogram			ScanplanningTopoChild		6 x 0.70				80	20	SPR (32 cm) 0.02	0.50	Standard					3.0															
		Topogram 0.70 Tr20 cor																										Tr20							
		Pelvis		Pelvis	RoutineSpiralChildPelvis		32 x 0.70	On	170	Auto (Non-Contrast)	110	165	(32 cm) 11.6	277	Standard		1.50	1.00		3.0															
			Pelvis 5.00 Br40 I3 sag																				5.00	5.00		None	Orientation Multi Recons (3)	MPR	Br40	On	3				
			Pelvis 5.00 Br40 I3 cor																																
			Pelvis 5.00 Br40 I3 ax																				5.00	5.00		None	None	MPR	Br60	On	3				
			Pelvis 5.00 Br60 I3 ax																																
			Pelvis Br40 I3 cor Radial																							Radial Ranges	None	VRT	Br40	On	3				
Hip [factory] (Child)	>Indications: evaluation of joint cavity, trauma, congruence evaluations, orthopedic indications, wo CM >Scan modes: RoutineSpiralChildPelvis	Topogram			ScanplanningTopoChild		6 x 0.70				80	20	SPR (32 cm) 0.02	0.50	Standard					3.0															
		Topogram 0.70 Tr20 cor																										Tr20							
		Hip		Pelvis	RoutineSpiralChildPelvis		32 x 0.70	On	170	Auto (Non-Contrast)	110	165	(32 cm) 11.6	161	Standard		1.50	1.00		3.0															
			Hip 3.00 Br60 I3 ax																				3.00	3.00		None	Impression Multi Recons (2)	MPR	Br60	On	3				
			Hip 3.00 Br40 I3 ax																										Br40	On	3				
			Hip 3.00 Br60 I3 cor R																				3.00	3.00		Parallel Ranges	None	MPR	Br60	On	3				
			Hip 3.00 Br60 I3 cor L																				3.00	3.00		Parallel Ranges	None	MPR	Br60	On	3				
C-Spine [factory] (Child)	>Indications: tumors, malformations, bone tumor, trauma, osteolysis of the cervical spine, wo CM >Scan modes: RoutineSpiralChildNeck	Topogram			ScanplanningTopoChild		6 x 0.70				130	15	SPR (32 cm) 0.06	1.51	Standard					7.0															
		Topogram 0.70 Tr20 sag																										Tr20							
		C-Spine		Neck	RoutineSpiralChildNeck		32 x 0.70	On	195	Auto (Non-Contrast)	110	174	(32 cm) 12.2	207	Standard		1.50	1.00		7.0															
			C-Spine 2.00 Br40 I3 ax																				2.00	2.00		None	Impression Multi Recons (2)	MPR	Br40	On	3				
			C-Spine 2.00 Br60 I3 ax																										Br60	On	3				
			C-Spine 2.00 Br40 I3 ax																				2.00	2.00		None	None	MPR	Br40	On	3				
			C-Spine 2.00 Br60 I3 ax																				2.00	2.00		None	Impression Multi Recons (2)	MPR	Br60	On	3				
			C-Spine 2.00 Br40 I3 ax																																
			C-Spine 0.80 Br60 I3 ax																				0.80	0.60		None	None	MPR	Br60	On	3				
			C-Spine Br40 I3 sag Radial																							Radial Ranges	None	VRT	Br40	On	3				

Protocol Name	Protocol Description	Range Name	Recon Name	Organ Characteristic	Scan Mode	Decision Trees	Acquisition [mm]	CARE Dose4D	FAST kV IQ level	FAST kV optimized for	kV	Eff. mAs	mA	CTDIvol [mGy]	DLP [mGy*cm]	Shaped Filter	Feed [mm]	Pitch	Rot. Time [s]	Scan Time [s]	Absolute Start Delay [s]	Contrast Protocol	Slice Thickness [mm]	Increment [mm]	Recon Strategies	Inline Results	Multi Recons	Image Type	Kernel	Iterative Recon	IR Strength	Recon Type	Recon Phase Step	Recon Phase(s)	Recon Phase Unit
T- or L-Spine [factory] (Child)	>Indications: tumors, malformations, bone tumor, trauma, osteolysis of the lumbar spine, wo CM >Scan modes: RoutineSpiralChildSpine	Topogram			ScanplanningTopoChild		6 x 0.70				130	15		SPR (32 cm) 0.06	1.51	Standard					3.0														
		Topogram 0.70 Tr20 sag																											Tr20						
		T- or L-Spine		Spine	RoutineSpiralChildSpine		32 x 0.70	On	285	Auto (Non-Contrast)	110	174		(32 cm) 12.2	207	Standard		1.50	1.00		3.0														
		T- or L-Spine 3.00 Br40 I3 ax																					3.00	3.00		None	Impression Multi Recons (2)	MPR	Br40	On	3				
		T- or L-Spine 3.00 Br60 I3 ax																								None	Impression Multi Recons (2)	MPR	Br60	On	3				
		T- or L-Spine 3.00 Br40 I3 ax																					3.00	3.00		None	None	MPR	Br40	On	3				
		T- or L-Spine 3.00 Br60 I3 ax																					3.00	3.00		None	Impression Multi Recons (2)	MPR	Br60	On	3				
		T- or L-Spine 3.00 Br40 I3 ax																								None	None	MPR	Br40	On	3				
		T- or L-Spine 3.00 Br60 I3 ax																					3.00	3.00		None	None	MPR	Br60	On	3				
Hand [factory] (Child)	>Indications: trauma, disorders of the hand bones, wo CM >Scan modes: RoutineSpiralChildExtremitiesSn	Topogram			ScanplanningTopoChild		6 x 0.70				80	20		SPR (32 cm) 0.02	0.50	Standard					3.0														
Topogram 0.70 Tr20 cor																												Tr20							
Hand			Extremities	RoutineSpiralChildExtremitiesSn		32 x 0.70	On	75	Auto (Bone/Calcium)	Sn 110	326		(32 cm) 4.65	107	Standard		0.80	1.00		3.0															
Hand 3.00 Br60 I3 sag R																						3.00	3.00		None	Orientation Multi Recons (3)	MPR	Br60	On	3					
Hand 3.00 Br60 I3 cor R																																			
Hand 3.00 Br60 I3 ax R																																			
Hand 3.00 Br40 I3 ax R																						3.00	3.00		None	None	MPR	Br40	On	3					
Hand 3.00 Br60 I3 sag L																						3.00	3.00		None	Orientation Multi Recons (3)	MPR	Br60	On	3					
Hand 3.00 Br60 I3 cor L																																			
Arm or Wrist [factory] (Child)	>Indications: trauma, disorders of the humerus, radius, ulna bones, wrist joint, wo CM >Scan modes: RoutineSpiralChildExtremitiesSn	Topogram			ScanplanningTopoChild		6 x 0.70				80	20		SPR (32 cm) 0.02	0.50	Standard					3.0														
Topogram 0.70 Tr20 cor																												Tr20							
Arm or Wrist			Extremities	RoutineSpiralChildExtremitiesSn		32 x 0.70	On	75	Auto (Bone/Calcium)	Sn 110	326		(32 cm) 4.65	41.4	Standard		0.80	1.00		3.0															
Arm or Wrist 2.00 Br60 I3 sag R																						2.00	2.00		None	Orientation Multi Recons (3)	MPR	Br60	On	3					
Arm or Wrist 2.00 Br60 I3 cor R																																			
Arm or Wrist 2.00 Br60 I3 ax R																																			
Arm or Wrist 2.00 Br40 I3 ax R																						2.00	2.00		None	None	MPR	Br40	On	3					
Arm or Wrist 2.00 Br60 I3 sag L																																			
Arm or Wrist 2.00 Br60 I3 cor L																						2.00	2.00		None	Orientation Multi Recons (3)	MPR	Br60	On	3					
Arm or Wrist 2.00 Br60 I3 ax L																																			
Arm or Wrist 2.00 Br40 I3 ax L																																			

Protocol Name	Protocol Description	Range Name	Recon Name	Organ Characteristic	Scan Mode	Decision Trees	Acquisition [mm]	CARE Dose4D	FAST kV IQ level	FAST kV optimized for	kV	Eff. mAs	mA	CTDIvol [mGy]	DLP [mGy*cm]	Shaped Filter	Feed [mm]	Pitch	Rot. Time [s]	Scan Time [s]	Absolute Start Delay [s]	Contrast Protocol	Slice Thickness [mm]	Increment [mm]	Recon Strategies	Inline Results	Multi Recons	Image Type	Kernel	Iterative Recon	IR Strength	Recon Type	Recon Phase Step	Recon Phase(s)	Recon Phase Unit		
Elbow [factory] (Child)	>Indications: trauma, disorders of the elbow joint, wo CM >Scan modes: RoutineSpiralChildExtremitiesSn	Topogram			ScanplanningTopoChild		6 x 0.70				80	20		SPR (32 cm) 0.02	0.50	Standard						3.0															
		Topogram 0.70 Tr20 cor																												Tr20							
		Elbow		Extremities	RoutineSpiralChildExtremitiesSn		32 x 0.70	On	75	Auto (Bone/Calcium)	Sn 110	326		(32 cm) 4.65	60.0	Standard		0.80	1.00			3.0															
		Elbow 2.00 Br60 I3 sag R																					2.00	2.00		None	Orientation Multi Recons (3)	MPR	Br60	On	3						
		Elbow 2.00 Br60 I3 cor R																																			
		Elbow 2.00 Br60 I3 ax R																																			
		Elbow 2.00 Br40 I3 ax R																					2.00	2.00		None	None	MPR	Br40	On	3						
		Elbow 2.00 Br60 I3 sag L																					2.00	2.00		None	Orientation Multi Recons (3)	MPR	Br60	On	3						
		Elbow 2.00 Br60 I3 cor L																																			
		Elbow 2.00 Br60 I3 ax L																																			
Foot [factory] (Child)	>Indications: trauma, disorders of the foot bones, wo CM >Scan modes: RoutineSpiralChildExtremitiesSn	Topogram			ScanplanningTopoChild		6 x 0.70				80	20		SPR (32 cm) 0.02	0.50	Standard						3.0															
Topogram 0.70 Tr20 cor																													Tr20								
Foot			Extremities	RoutineSpiralChildExtremitiesSn		32 x 0.70	On	75	Auto (Bone/Calcium)	Sn 110	326		(32 cm) 4.65	60.0	Standard		0.80	1.00			3.0																
Foot 3.00 Br60 I3 sag R																						3.00	3.00		None	Orientation Multi Recons (3)	MPR	Br60	On	3							
Foot 3.00 Br60 I3 cor R																																					
Foot 3.00 Br60 I3 ax R																																					
Foot 3.00 Br40 I3 ax R																						3.00	3.00		None	None	MPR	Br40	On	3							
Foot 3.00 Br60 I3 sag L																																					
Foot 3.00 Br60 I3 cor L																																					
Foot 3.00 Br60 I3 ax L																																					
Leg or Ankle [factory] (Child)	>Indications: trauma, disorders of the femur, fibula, tibia bones, ankle joint, wo CM >Scan modes: RoutineSpiralChildExtremitiesSn	Topogram			ScanplanningTopoChild		6 x 0.70				130	15		SPR (32 cm) 0.06	1.51	Standard						3.0															
Topogram 0.70 Tr20 sag																													Tr20								
Leg or Ankle			Extremities	RoutineSpiralChildExtremitiesSn		32 x 0.70	On	75	Auto (Bone/Calcium)	Sn 110	326		(32 cm) 4.65	60.0	Standard		0.80	1.00			3.0																
Leg or Ankle 2.00 Br60 I3 sag R																						2.00	2.00		None	Orientation Multi Recons (3)	MPR	Br60	On	3							
Leg or Ankle 2.00 Br60 I3 cor R																																					
Leg or Ankle 2.00 Br60 I3 ax R																																					
Leg or Ankle 2.00 Br40 I3 ax R																																					
Leg or Ankle 2.00 Br60 I3 sag L																																					
Leg or Ankle 2.00 Br60 I3 cor L																																					
Leg or Ankle 2.00 Br60 I3 ax L																																					
Leg or Ankle 2.00 Br40 I3 ax L																																					

Protocol Name	Protocol Description	Range Name	Recon Name	Organ Characteristic	Scan Mode	Decision Trees	Acquisition [mm]	CARE Dose4D	FAST kV IQ level	FAST kV optimized for	kV	Eff. mAs	mA	CTDIvol [mGy]	DLP [mGy*cm]	Shaped Filter	Feed [mm]	Pitch	Rot. Time [s]	Scan Time [s]	Absolute Start Delay [s]	Contrast Protocol	Slice Thickness [mm]	Increment [mm]	Recon Strategies	Inline Results	Multi Recons	Image Type	Kernel	Iterative Recon	IR Strength	Recon Type	Recon Phase Step	Recon Phase(s)	Recon Phase Unit			
Knee [factory] (Child)	>Indications: trauma, disorders of the knee joint, w/o CM >Scan modes: RoutineSpiralChildExtremitiesSn	Topogram			ScanplanningTopoChild		6 x 0.70				80	20		SPR (32 cm) 0.02	0.50	Standard					3.0																	
			Topogram 0.70 Tr20 cor																										Tr20									
		Knee		Extremities	RoutineSpiralChildExtremitiesSn		32 x 0.70	On	75	Auto (Bone/Calcium)	Sn 110	326		(32 cm) 4.65	73.9	Standard		0.80	1.00		3.0																	
			Knee 3.00 Br60 I3 sag R Knee 3.00 Br60 I3 cor R Knee 3.00 Br60 I3 ax R																					3.00	3.00		None	Orientation Multi Recons (3)	MPR	Br60	On	3						
			Knee 3.00 Br40 I3 ax R																					3.00	3.00		None	None	MPR	Br40	On	3						
			Knee 3.00 Br60 I3 cor R																					3.00	3.00		Parallel Ranges	None	MPR	Br60	On	3						
			Knee 3.00 Br60 I3 sag L Knee 3.00 Br60 I3 cor L																					3.00	3.00		None	Orientation Multi Recons (3)	MPR	Br60	On	3						
			Knee 3.00 Br60 I3 ax L																																			
			Knee 3.00 Br40 I3 ax L																						3.00	3.00		None	None	MPR	Br40	On	3					
			Knee 3.00 Br60 I3 cor L																						3.00	3.00		Parallel Ranges	None	MPR	Br60	On	3					
		Brain CTA [factory] (Child)	>Indications: stenosis, aneurysm, occlusion, plaques, abnormalities of the vertebral arteries, COW, w CM >Scan modes: PremonitoringSequenceChild2d, MonitoringSequenceChild2d, RoutineSpiralChildAngio Head	Topogram			ScanplanningTopoChild		6 x 0.70				130	15		SPR (32 cm) 0.06	1.51	Standard					3.0															
					Topogram 0.70 Tr20 sag																										Tr20							
Premonitoring				Vascular Body			1 x 5.00	On	25	Soft Tissue with Contrast	80	47		(32 cm) 1.15	0.58	Standard	0.0		0.80	Full 0.80	3.0																	
	Premonitoring 5.00 Br36																							5.00				None		Br36								
Monitoring				Vascular Body			1 x 5.00	On	25	Soft Tissue with Contrast	80	47		(32 cm) 34.6	17.3	Standard	0.0		0.80	Full 0.80	10.0																	
	Monitoring 5.00 Br36																							5.00				None		Br36								
CTA Brain				Head/Vascular Head	RoutineSpiralChildAngio Head		32 x 0.70	On	105	Auto (Vascular)	80	94		(16 cm) 5.76	80.5	Standard		1.50	0.80																			
	CTA Brain 3.00 Hv36 I3 cor CTA Brain 3.00 Hv36 I3 sag CTA Brain 3.00 Hv36 I3 ax																							3.00	3.00		None	Orientation Multi Recons (3)	MIP Thin	Hv36	On	3						
Neck CTA [factory] (Child)	>Indications: stenosis, occlusion, plaques, abnormalities of the carotid arteries, w CM >Scan modes: PremonitoringSequenceChild2d, MonitoringSequenceChild2d, RoutineSpiralChildAngio Neck	Topogram			ScanplanningTopoChild		6 x 0.70				80	20		SPR (32 cm) 0.02	0.50	Standard					7.0																	
			Topogram 0.70 Tr20 cor																										Tr20									
		Premonitoring		Vascular Body			1 x 5.00	On	25	Soft Tissue with Contrast	80	47		(32 cm) 1.15	0.58	Standard	0.0		0.80	Full 0.80	3.0																	
			Premonitoring 5.00 Br36																					5.00				None		Br36								
		Monitoring		Vascular Body			1 x 5.00	On	25	Soft Tissue with Contrast	80	47		(32 cm) 34.6	17.3	Standard	0.0		0.80	Full 0.80	10.0																	
			Monitoring 5.00 Br36																					5.00				None		Br36								
		CTA Neck		Neck	RoutineSpiralChildAngio Neck		32 x 0.70	On	145	Auto (Vascular)	80	140		(32 cm) 3.74	89.7	Standard		1.50	0.80																			
			CTA Neck 3.00 Bv36 I3 sag CTA Neck 3.00 Bv36 I3 cor CTA Neck 3.00 Bv36 I3 ax																					3.00	3.00		None	Orientation Multi Recons (3)	MIP Thin	Bv36	On	3						
			CTA Neck 2.00 Bv36 I3 cor																					2.00	2.00		Parallel Ranges	None	MIP Thin	Bv36	On	3						
			CTA Neck 3.00 Br40 I3 Parallel Vascular																					3.00	3.00		Parallel Vascular Ranges	None	MPR	Br40	On	3						
			CTA Neck 3.00 Br40 I3 Parallel Vascular																					3.00	3.00		Parallel Vascular Ranges	None	MPR	Br40	On	3						

Protocol Name	Protocol Description	Range Name	Recon Name	Organ Characteristic	Scan Mode	Decision Trees	Acquisition [mm]	CARE Dose4D	FAST kV IQ level	FAST kV optimized for	kV	Eff. mAs	CTDIvol [mGy]	DLP [mGy*cm]	Shaped Filter	Feed [mm]	Pitch	Rot. Time [s]	Scan Time [s]	Absolute Start Delay [s]	Contrast Protocol	Slice Thickness [mm]	Increment [mm]	Recon Strategies	Inline Results	Multi Recons	Image Type	Kernel	Iterative Recon	IR Strength	Recon Type	Recon Phase Step	Recon Phase(s)	Recon Phase Unit			
Thorax CTA [factory] (Child)	>Indications: vascular anomalies, aneurysms, dissection, ruptures of the thoracic aorta, w CM >Scan modes: PremonitoringSequenceChild2d, MonitoringSequenceChild2d, RoutineSpiralChildAngioBody	Topogram			ScanplanningTopoChild		6 x 0.70				80	20	SPR (32 cm) 0.02	0.81	Standard					6.0																	
		Topogram 0.70 Tr20 cor																										Tr20									
		Premonitoring		Vascular Body			1 x 5.00	On	25	Soft Tissue with Contrast	80	47	(32 cm) 1.15	0.58	Standard	0.0		0.80	Full 0.80	3.0																	
		Premonitoring 5.00 Br36																				5.00			None		Br36										
		Monitoring		Vascular Body			1 x 5.00	On	25	Soft Tissue with Contrast	80	47	(32 cm) 34.6	17.3	Standard	0.0		0.80	Full 0.80	10.0																	
		Monitoring 5.00 Br36																				5.00				None		Br36									
		CTA Thorax		Vascular Body	RoutineSpiralChildAngioBody		32 x 0.70	On	70	Auto (Vascular)	80	92	(32 cm) 2.46	59.0	Standard		1.50	0.80																			
		CTA Thorax 5.00 Bv36 I3 sag																					5.00	5.00		None	Orientation Multi Recons (3)	MIP Thin	Bv36	On	3						
		CTA Thorax 5.00 Bv36 I3 cor																																			
		CTA Thorax 5.00 Bv36 I3 ax																																			
		CTA Thorax 5.00 BI60 I3 ax																					5.00	5.00		None	Impression Multi Recons (2)	MIP Thin	BI60	On	3						
		CTA Thorax 5.00 Br40 I3 ax																											Br40	On	3						
Pulmonary CTA [factory] (Child)	>Indications: pulmonary embolism, w CM >Scan modes: PremonitoringSequenceChild2d, MonitoringSequenceChild2d, RoutineSpiralChildAngioBody	Topogram			ScanplanningTopoChild		6 x 0.70				80	20	SPR (32 cm) 0.02	0.81	Standard					4.0																	
		Topogram 0.70 Tr20 cor																										Tr20									
		Premonitoring		Vascular Body			1 x 5.00	On	25	Soft Tissue with Contrast	80	47	(32 cm) 1.15	0.58	Standard	0.0		0.80	Full 0.80	3.0																	
		Premonitoring 5.00 Br36																				5.00			None		Br36										
		Monitoring		Vascular Body			1 x 5.00	On	25	Soft Tissue with Contrast	80	47	(32 cm) 34.6	17.3	Standard	0.0		0.80	Full 0.80	10.0																	
		Monitoring 5.00 Br36																				5.00				None		Br36									
		CTA Pulmonary		Vascular Body	RoutineSpiralChildAngioBody		32 x 0.70	On	70	Auto (Vascular)	80	92	(32 cm) 2.46	59.0	Standard		1.50	0.80																			
		CTA Pulmonary 3.00 Bv36 I3 sag																					3.00	3.00		None	Orientation Multi Recons (3)	MIP Thin	Bv36	On	3						
		CTA Pulmonary 3.00 Bv36 I3 cor																																			
CTA Pulmonary 3.00 Bv36 I3 ax																																					
CTA Pulmonary 3.00 BI60 I3 ax																						3.00	3.00		None	Impression Multi Recons (2)	MIP Thin	BI60	On	3							
CTA Pulmonary 3.00 Br40 I3 ax																											Br40	On	3								

Protocol Name	Protocol Description	Range Name	Recon Name	Organ Characteristic	Scan Mode	Decision Trees	Acquisition [mm]	CARE Dose4D	FAST kV IQ level	FAST kV optimized for	kV	Eff. mAs	CTDIvol [mGy]	DLP [mGy*cm]	Shaped Filter	Feed [mm]	Pitch	Rot. Time [s]	Scan Time [s]	Absolute Start Delay [s]	Contrast Protocol	Slice Thickness [mm]	Increment [mm]	Recon Strategies	Inline Results	Multi Recons	Image Type	Kernel	Iterative Recon	IR Strength	Recon Type	Recon Phase Step	Recon Phase(s)	Recon Phase Unit
Abdomen CTA [factory] (Child)	>Indications: vascular anomalies, aneurysms, dissection, ruptures of the abdominal aorta, renal arteries, w CM >Scan modes: PremonitoringSequenceChild2d, MonitoringSequenceChild2d, RoutineSpiralChildAngio Body	Topogram			ScanplanningTopoChild		6 x 0.70				80	20	SPR (32 cm) 0.02	0.81	Standard					6.0														
			Topogram 0.70 Tr20 cor																									Tr20						
		Premonitoring		Vascular Body			1 x 5.00	On	25	Soft Tissue with Contrast	80	47	(32 cm) 1.15	0.58	Standard	0.0		0.80	Full 0.80	3.0														
			Premonitoring 5.00 Br36																			5.00				None		Br36						
		Monitoring		Vascular Body			1 x 5.00	On	25	Soft Tissue with Contrast	80	47	(32 cm) 34.6	17.3	Standard	0.0		0.80	Full 0.80	10.0														
			Monitoring 5.00 Br36																			5.00				None		Br36						
		CTA Abdomen		Vascular Body	RoutineSpiralChildAngio Body		32 x 0.70	On	115	Auto (Vascular)	80	148	(32 cm) 3.96	134	Standard		1.50	0.80																
			CTA Abdomen 5.00 Bv36 I3 sag																			5.00	5.00		None	Orientation Multi Recons (3)	MIP Thin	Bv36	On	3				
			CTA Abdomen 5.00 Bv36 I3 cor																															
			CTA Abdomen 5.00 Bv36 I3 ax																															
Whole Body CTA [factory] (Child)	>Indications: vascular anomalies, aneurysms, dissection, ruptures of the whole aorta, w CM >Scan modes: PremonitoringSequenceChild2d, MonitoringSequenceChild2d, RoutineSpiralChildAngio Body	Topogram			ScanplanningTopoChild		6 x 0.70				80	20	SPR (32 cm) 0.02	0.81	Standard					6.0														
			Topogram 0.70 Tr20 cor																									Tr20						
		Premonitoring		Vascular Body			1 x 5.00	On	25	Soft Tissue with Contrast	80	47	(32 cm) 1.15	0.58	Standard	0.0		0.80	Full 0.80	3.0														
			Premonitoring 5.00 Br36																			5.00				None		Br36						
		Monitoring		Vascular Body			1 x 5.00	On	25	Soft Tissue with Contrast	80	47	(32 cm) 34.6	17.3	Standard	0.0		0.80	Full 0.80	10.0														
			Monitoring 5.00 Br36																			5.00				None		Br36						
		CTA Whole Body		Vascular Body	RoutineSpiralChildAngio Body		32 x 0.70	On	115	Auto (Vascular)	80	148	(32 cm) 3.96	214	Standard		1.50	0.80																
			CTA Whole Body 5.00 Bv36 I3 sag																			5.00	5.00		None	Orientation Multi Recons (3)	MIP Thin	Bv36	On	3				
			CTA Whole Body 5.00 Bv36 I3 cor																															
			CTA Whole Body 5.00 Bv36 I3 ax																															
Extremity CTA [factory] (Child)	>Indications: vascular anomalies, aneurysms of the femoral arteries, w CM >Scan modes: PremonitoringSequenceChild2d, MonitoringSequenceChild2d, RoutineSpiralChildAngio Extremities	Topogram			ScanplanningTopoChild		6 x 0.70				80	20	SPR (32 cm) 0.02	0.81	Standard					3.0														
			Topogram 0.70 Tr20 cor																									Tr20						
		Premonitoring		Vascular Body			1 x 5.00	On	25	Soft Tissue with Contrast	80	47	(32 cm) 1.15	0.58	Standard	0.0		0.80	Full 0.80	3.0														
			Premonitoring 5.00 Br36																			5.00				None		Br36						
		Monitoring		Vascular Body			1 x 5.00	On	25	Soft Tissue with Contrast	80	47	(32 cm) 34.6	17.3	Standard	0.0		0.80	Full 0.80	10.0														
			Monitoring 5.00 Br36																			5.00				None		Br36						
		CTA Extremity		Runoff	RoutineSpiralChildAngio Extremities		32 x 0.70	On	155	Auto (Vascular)	80	197	(32 cm) 5.27	284	Standard		1.20	0.80																
			CTA Extremity 5.00 Bv36 I3 sag																			5.00	5.00		None	Orientation Multi Recons (3)	MIP Thin	Bv36	On	3				
			CTA Extremity 5.00 Bv36 I3 cor																															
			CTA Extremity 5.00 Bv36 I3 ax																															

Protocol Name	Protocol Description	Range Name	Recon Name	Organ Characteristic	Scan Mode	Decision Trees	Acquisition [mm]	CARE Dose4D	FAST kV IQ level	FAST kV optimized for	kV	Eff. mAs	CTDIvol [mGy]	DLP [mGy*cm]	Shaped Filter	Feed [mm]	Pitch	Rot. Time [s]	Scan Time [s]	Absolute Start Delay [s]	Contrast Protocol	Slice Thickness [mm]	Increment [mm]	Recon Strategies	Inline Results	Multi Recons	Image Type	Kernel	Iterative Recon	IR Strength	Recon Type	Recon Phase Step	Recon Phase(s)	Recon Phase Unit	
Head Trauma [factory] (Child)	>Indications: emergency head trauma, hemorrhage, w/o CM >Scan modes: RoutineSpiralChildHead	Topogram			ScanplanningTopoChild		6 x 0.70				80	20	SPR (16 cm) 0.05	1.14	Standard					3.0															
		Topogram 0.70 Tr20 cor																									Tr20								
		Head	Brain	RoutineSpiralChildHead		32 x 0.70	On	493	Auto (Non-Contrast)	110	220	(16 cm) 32.8	505	Standard		0.95	0.80		3.0																
		Head 5.00 Hr40 ax Head 5.00 Hr40 sag Head 5.00 Hr40 cor																					5.00	5.00		None	Orientation Multi Recons (3)	MPR	Hr40	Off	2				
			Head 0.80 Hr60 ax Head 0.80 Hr60 sag Head 0.80 Hr60 cor																					0.80	0.60		None	Orientation Multi Recons (3)	MPR	Hr60	Off	2			
Body Trauma [factory] (Child)	>Indications: fast screening for emergency whole body trauma, hemorrhage, w/o CM >Scan modes: RoutineSpiralChildAbdomen	Topogram			ScanplanningTopoChild		6 x 0.70				80	20	SPR (32 cm) 0.02	0.81	Standard					6.0															
		Topogram 0.70 Tr20 cor																									Tr20								
		Body	Abdomen	RoutineSpiralChildAbdomen		32 x 0.70	On	207	Auto (Non-Contrast)	110	139	(32 cm) 9.75	526	Standard		1.50	0.80		6.0																
		Body 1.50 Br40 ax Body 1.50 Br60 ax Body 1.50 Bl60 ax																					1.50	1.00		None	Impression Multi Recons (3)	MPR	Br40	Off	3				
																												Br60	Off	3					
																												Bl60	Off	3					
RTP Head [factory] (Child)	>Indications: radiotherapy planning of the head region, w/o CM >Scan modes: RoutineSpiralChildHead	Topogram			ScanplanningTopoChild		6 x 0.70				130	15	SPR (16 cm) 0.10	5.07	Standard					3.0															
		Topogram 0.70 Tr20 sag																									Tr20								
		Head	Brain	RoutineSpiralChildHead		32 x 0.70	On	345	Auto (Non-Contrast)	130	230	(16 cm) 51.0	614	Standard		0.55	1.00		3.0																
		Head 1.50 Qr40 I2																					1.50	1.50		None	None		Qr40	On	2				
RTP Head and Neck [factory] (Child)	>Indications: radiotherapy planning of the head, neck and shoulder region, w/o CM >Scan modes: RoutineSpiralChildNeck	Topogram			ScanplanningTopoChild		6 x 0.70				80	20	SPR (32 cm) 0.02	0.81	Standard					3.0															
		Topogram 0.70 Tr20 cor																									Tr20								
		Head and Neck	Neck	RoutineSpiralChildNeck		32 x 0.70	On	155	Auto (Non-Contrast)	130	102	(32 cm) 10.9	305	Standard		0.80	1.00		3.0																
		Head and Neck 1.50 Qr40 I3																					1.50	1.50		None	None		Qr40	On	3				
RTP Thorax [factory] (Child)	>Indications: radiotherapy planning of the thorax region, w/o CM >Scan modes: RoutineSpiralChildThorax	Topogram			ScanplanningTopoChild		6 x 0.70				80	20	SPR (32 cm) 0.02	0.81	Standard					3.0															
		Topogram 0.70 Tr20 cor																									Tr20								
		Thorax	Thorax	RoutineSpiralChildThorax		32 x 0.70	On	80	Auto (Non-Contrast)	130	51	(32 cm) 5.47	125	Standard		0.80	0.80		3.0																
		Thorax 1.50 Qr40 I3																					1.50	1.50		None	None		Qr40	On	3				
RTP Respiratory Motion Management [factory] (Child)	>Indications: motion management for radiotherapy planning for moving tumors, w/o CM >Scan modes: RespiratorySpiralChild	Topogram			ScanplanningTopoChild		6 x 0.70				80	20	SPR (32 cm) 0.02	0.81	Standard					3.0															
		Topogram 0.70 Tr20 cor																									Tr20								
		Respiratory Motion	Respiratory	RespiratorySpiralChild		32 x 0.70	On	26	Auto (Non-Contrast)	130	17	(32 cm) 30.0	678	Standard		0.12	0.80		3.0																
		Respiratory Motion 1.50 Qr40 I3 Average CT																					1.50	1.50		None	None		Qr40	On	3	Average CT			
		Respiratory Motion 1.50 Qr40 I3 0% - 90%																					1.50	1.50		None	None		Qr40	On	3	Multi Phase	10	0%; 10%; 20%; 30%; 40%; 50%; 60%; 70%; 80%; 90%	Phase (%)

Protocol Name	Protocol Description	Range Name	Recon Name	Organ Characteristic	Scan Mode	Decision Trees	Acquisition [mm]	CARE Dose4D	FAST kV IQ level	FAST kV optimized for	kV	Eff. mAs	CTDIvol [mGy]	DLP [mGy*cm]	Shaped Filter	Feed [mm]	Pitch	Rot. Time [s]	Scan Time [s]	Absolute Start Delay [s]	Contrast Protocol	Slice Thickness [mm]	Increment [mm]	Recon Strategies	Inline Results	Multi Recons	Image Type	Kernel	Iterative Recon	IR Strength	Recon Type	Recon Phase Step	Recon Phase(s)	Recon Phase Unit	
RTP Abdomen [factory] (Child)	>Indications: radiotherapy planning of the abdomen region, w/wo CM >Scan modes: RoutineSpiralChildAbdomen	Topogram			ScanplanningTopoChild		6 x 0.70				80	20	SPR (32 cm) 0.02	0.81	Standard					3.0															
		Topogram 0.70 Tr20 cor																										Tr20							
		Abdomen		Abdomen	RoutineSpiralChildAbdomen		32 x 0.70	On	145	Auto (Non-Contrast)	130	92	(32 cm) 9.87	325	Standard		0.80	0.80		3.0															
		Abdomen 1.50 Qr40 I3																				1.50	1.50		None	None		Qr40	On	3					
RTP Extremity [factory] (Child)	>Indications: radiotherapy planning of the extremities region, w/wo CM >Scan modes: RoutineSpiralChildExtremities	Topogram			ScanplanningTopoChild		6 x 0.70				80	20	SPR (32 cm) 0.02	0.50	Standard					3.0															
		Topogram 0.70 Tr20 cor																										Tr20							
		Extremity		Extremities	RoutineSpiralChildExtremities		32 x 0.70	On	77	Non-Contrast	130	52	(32 cm) 5.58	71.9	Standard		0.80	1.00		3.0															
		Extremity 1.50 Qr40 I3 R Extremity 1.50 Qr40 I3 L																				1.50	1.50		None	None		Qr40	On	3					
CTDI measurement Head [factory] (Child)	>Indications: special head protocol for ACR accreditation only; don't use for clinical cases! >Scan modes: RoutineSequenceChildHead	Topogram			ScanplanningTopoChild		6 x 0.70				80	20	SPR (16 cm) 0.05	1.14	Standard					3.0															
		Topogram 0.70 Tr20 cor																										Tr20							
		Head		Brain	RoutineSequenceChildHead		32 x 0.70	Off			130	43	(16 cm) 10.7	150	Standard	20.0		1.00	Full 1.00	3.0															
		Head 5.00 Hr40																				5.00	5.00		None	None		Hr40	Off	3					
CTDI measurement Body [factory] (Child)	>Indications: special body protocol for ACR accreditation only; don't use for clinical cases! >Scan modes: RoutineSequenceChildAbdomen	Topogram			ScanplanningTopoChild		6 x 0.70				80	20	SPR (32 cm) 0.02	0.50	Standard					3.0															
		Topogram 0.70 Tr20 cor																										Tr20							
		Body		Abdomen	RoutineSequenceChildAbdomen		32 x 0.70	Off			130	41	(32 cm) 4.93	69.0	Standard	20.0		1.00	Full 1.00	3.0															
Brain Baby [factory] (Child)	>Indications: head trauma, hemorrhage, hydrocephalus, abnormalities, inflammation, w/wo CM >Scan modes: RoutineSpiralChildHead	Topogram			ScanplanningTopoChild		6 x 0.70				130	15	SPR (16 cm) 0.12	3.11	Standard					3.0															
		Topogram 0.70 Tr20 sag																										Tr20							
		Brain		Brain	RoutineSpiralChildHead		32 x 0.70	On	338	Auto (Non-Contrast)	110	139	(16 cm) 20.8	249	Standard		1.50	0.80		3.0															
		Brain 3.00 Hr40 I2 Brain 3.00 Hr60 I2																				3.00	3.00		None	Impression Multi Recons (2)		Hr40 Hr60	On	2 2					
Body Baby [factory] (Child)	>Indications: masses, tumor, lymphoma of the thorax,liver, pancreas, kidneys, gastrointestinal studies, w/wo CM >Scan modes: RoutineSpiralChildAbdomen	Topogram			ScanplanningTopoChild		6 x 0.70				80	20	SPR (32 cm) 0.02	0.81	Standard					3.0															
		Topogram 0.70 Tr20 cor																										Tr20							
		Body		Abdomen	RoutineSpiralChildAbdomen		32 x 0.70	On	145	Auto (Non-Contrast)	110	139	(32 cm) 9.75	234	Standard		1.50	0.80		3.0															
		Body 3.00 Br40 I3 cor Body 3.00 Br40 I3 ax																				3.00	3.00		None	Orientation Multi Recons (2)	MPR	Br40	On	3					
		Body 3.00 Br60 I3 sag Body 3.00 Br60 I3 cor																				3.00	3.00		None	Orientation Multi Recons (3)	MPR	Br60	On	3					
		Body 3.00 Br60 I3 ax																																	
Thorax Intervention 3D [factory] (Child)	>Indications: biopsy, abscess in the thorax w/wo CM >Scan modes: InterventionSpiralChildThorax, InterventionSequenceChildThorax, InterventionFluoroChildThorax2d	Topogram			ScanplanningTopoChild		6 x 0.70				80	20	SPR (32 cm) 0.02	0.81	Standard					3.0															
		Topogram 0.70 Tr20 cor																										Tr20							
		i-Spiral Thorax		Thorax	InterventionSpiralChildThorax		32 x 0.70	On	80	Auto (Non-Contrast)	110	76	(32 cm) 5.32	42.0	Standard		1.50	0.80		6.0															
		i-Spiral Thorax 3.20 BI60																				3.20	3.00		None	None		BI60	Off	3					
		i-Sequence Thorax		Thorax	InterventionSequenceChildThorax		32 x 0.70	Off			110	13	(32 cm) 0.91	2.04	Standard	0.0		0.80	Full 0.80																
		i-Sequence Thorax 2.60 Br40																				2.60	2.20		None	None		Br40	Off	3					
		i-Fluoro Thorax		Thorax	InterventionFluoroChildThorax2d		3 x 3.00	Off			110	13			Standard			0.80																	

Protocol Name	Protocol Description	Range Name	Recon Name	Organ Characteristic	Scan Mode	Decision Trees	Acquisition [mm]	CARE Dose4D	FAST kV IQ level	FAST kV optimized for	kV	Eff. mAs	mA	CTDIvol [mGy]	DLP [mGy*cm]	Shaped Filter	Feed [mm]	Pitch	Rot. Time [s]	Scan Time [s]	Absolute Start Delay [s]	Contrast Protocol	Slice Thickness [mm]	Increment [mm]	Recon Strategies	Inline Results	Multi Recons	Image Type	Kernel	Iterative Recon	IR Strength	Recon Type	Recon Phase Step	Recon Phase(s)	Recon Phase Unit		
Thorax Intervention [factory] (Child)	>Indications: biopsy, abscess in the thorax w/wo CM >Scan modes: RoutineSpiralChildThorax2d, InterventionSequenceChildThorax2d, InterventionFluoroChildThorax2d, RoutineSpiralChildThorax	Topogram			ScanplanningTopoChild		6 x 0.70				80	20		SPR (32 cm) 0.02	0.81	Standard					3.0																
			Topogram 0.70 Tr20 cor																										Tr20								
		Planning Thorax		Thorax	RoutineSpiralChildThorax		32 x 0.70	On	80	Auto (Non-Contrast)	110	76		(32 cm) 5.32	128	Standard		1.50	0.80		6.0							None	None		BI60	On	3				
			Planning Thorax 3.00 BI60 I3																					3.00	3.00			None	None		BI60	On	3				
			Planning Thorax 3.00 Br40 I3																					3.00	3.00			None	None		Br40	On	3				
		i-Sequence Thorax		Thorax	InterventionSequenceChildThorax2d		3 x 3.00	Off				110	13		(32 cm) 0.93	0.83	Standard	0.0		0.80	Full 0.80																
			i-Sequence Thorax 3.00 Br40																					3.00				None		Br40							
		i-Fluoro Thorax		Thorax	InterventionFluoroChildThorax2d		3 x 3.00	Off				110	13				Standard			0.80																	
			i-Fluoro Thorax 3.00 BI60																						3.00	3.00			None		BI60						
		Control Thorax		Thorax	RoutineSpiralChildThorax		32 x 0.70	On	80	Auto (Non-Contrast)	110	76		(32 cm) 5.32	128	Standard		1.50	0.80		6.0																
			Control Thorax 3.00 BI60 I3																						3.00	3.00			None	None		BI60	On	3			
			Control Thorax 3.00 Br40 I3																						3.00	3.00			None	None		Br40	On	3			
Abdomen Intervention 3D [factory] (Child)	>Indications: biopsy, abscess in the abdomen, w/wo CM >Scan modes: InterventionSpiralChildAbdomen, InterventionSequenceChildAbdomen, InterventionFluoroChildAbdomen2d	Topogram			ScanplanningTopoChild		6 x 0.70				80	20		SPR (32 cm) 0.02	0.81	Standard					3.0																
			Topogram 0.70 Tr20 cor																										Tr20								
		i-Spiral Abdomen		Abdomen	InterventionSpiralChildAbdomen		32 x 0.70	On	145	Auto (Non-Contrast)	110	140		(32 cm) 9.81	76.4	Standard		1.20	0.80		6.0																
			i-Spiral Abdomen 3.20 Br40																					3.20	3.00			None	None		Br40	Off	3				
		i-Sequence Abdomen		Abdomen	InterventionSequenceChildAbdomen		32 x 0.70	Off				110	23		(32 cm) 1.61	3.61	Standard	0.0		0.80	Full 0.80																
			i-Sequence Abdomen 2.60 Br40																					2.60	2.20			None	None		Br40	Off	3				
		i-Fluoro Abdomen		Abdomen	InterventionFluoroChildAbdomen2d		3 x 3.00	Off				110	23				Standard			0.80																	
	i-Fluoro Abdomen 3.00 Br40																						3.00	3.00			None		Br40								
Abdomen Intervention [factory] (Child)	>Indications: biopsy, abscess in the abdomen, w/wo CM >Scan modes: RoutineSpiralChildAbdomen, InterventionSequenceChildAbdomen2d, InterventionFluoroChildAbdomen2d	Topogram			ScanplanningTopoChild		6 x 0.70				80	20		SPR (32 cm) 0.02	0.81	Standard					3.0																
			Topogram 0.70 Tr20 cor																										Tr20								
		Planning Abdomen		Abdomen	RoutineSpiralChildAbdomen		32 x 0.70	On	145	Auto (Non-Contrast)	110	140		(32 cm) 9.81	332	Standard		1.20	0.80		6.0																
			Planning Abdomen 3.00 Br40 I3																						3.00	3.00			None	None		Br40	On	3			
		i-Sequence Abdomen		Abdomen	InterventionSequenceChildAbdomen2d		3 x 3.00	Off				110	23		(32 cm) 1.64	1.48	Standard	0.0		0.80	Full 0.80																
			i-Sequence Abdomen 3.00 Br40																						3.00				None		Br40						
		i-Fluoro Abdomen		Abdomen	InterventionFluoroChildAbdomen2d		3 x 3.00	Off				110	23				Standard			0.80																	
	i-Fluoro Abdomen 3.00 Br40																						3.00	3.00			None		Br40								

Brain contrast medium necessary [factory]				
Scan Strategy				
	wo CM		w CM	
InfoLevel	None		None	
InfoText				
Scan Mode	RoutineSpiralAdultHead		RoutineSpiralAdultHead	
Injector Coupling	Off		On	
Contrast Protocol			Brain [factory]	
Absolute Start Delay	2		180	
Decision Tree				
wo CM				
Contrast medium?				
yes	no			
w CM	wo CM			
Side [factory]				
Scan Strategy				
	recon both		recon right	
InfoLevel	None		None	
InfoText			None	
Decision Tree				
recon both				
Side?				
both	right			
recon both	recon right			
Metal implant [factory]				
Scan Strategy				
	wo iMAR		w iMAR	
InfoLevel	None		None	
InfoText				
Decision Tree				
wo iMAR				
Metal implant?				
yes	no			
w iMAR	wo iMAR			
Sinus sinusitis or soft tissue [factory]				
Scan Strategy				
	low dose		standard dose	
InfoLevel	None		None	
InfoText				
Scan Mode	RoutineSpiralAdultHeadSnXcare		RoutineSpiralAdultHeadXcare	
CARE Dose4D	On		On	
FAST kV optimized for	Bone/Calcium		Non-Contrast	
FAST kV IQ level	50		100	
Decision Tree				
low dose				
Sinus Sinusitis/soft tissue?				
sinusitis	soft tissue			
low dose	standard dose			
Neck FoV [factory]				
Scan Strategy				
	standard acquisition		fast acquisition	
InfoLevel	Information		None	
InfoText	Please check scan time!			
Scan Mode	RoutineSpiralAdultNeck		RoutineSpiralAdultNeck	
Pitch	0.8		1.5	
Decision Tree				
standard acquisition				
FoV?				
yes	no			
fast acquisition	standard acquisition			

Neck contrast medium necessary [factory]									
Scan Strategy									
	wo CM		w CM						
InfoLevel	None		None						
InfoText									
Scan Mode	RoutineSpiralAdultNeck		RoutineSpiralAdultNeck						
CARE Dose4D	On		On						
FAST kV optimized for	Non-Contrast		Soft Tissue with Contrast						
Injector Coupling	Off		On						
Contrast Protocol			Neck [factory]						
Absolute Start Delay	2		15						
Decision Tree									
wo CM									
Contrast medium?									
yes	no								
w CM	wo CM								
Lung screening [factory]									
Scan Strategy									
	low dose wo X-CARE		low dose w X-CARE		high resolution wo X-CARE		high resolution w X-CARE		
InfoLevel	None		None		None		None		
InfoText									
Scan Mode	RoutineSpiralAdultThoraxSn		RoutineSpiralAdultThoraxSnXcare		RoutineSpiralAdultThorax		RoutineSpiralAdultThoraxSnXcare		
CARE Dose4D	On		On		On		On		
FAST kV IQ level	15		15		100		100		
Pitch	0.8		0.6		0.8		0.6		
Decision Tree									
low dose wo X-CARE									
Lung low dose/high resolution?									
lung cancer screening				high resolution					
low dose wo X-CARE				high resolution wo X-CARE					
Gender				Gender					
Male	Female	Other		Male	Female	Other			
low dose wo X-CARE	low dose w X-CARE	low dose w X-CARE		high resolution wo X-CARE	high resolution w X-CARE	high resolution w X-CARE			
Thorax contrast medium necessary [factory]									
Scan Strategy									
	wo CM		w CM						
InfoLevel	None		None						
InfoText									
Scan Mode	RoutineSpiralAdultThorax		RoutineSpiralAdultThorax						
CARE Dose4D	On		On						
FAST kV optimized for	Non-Contrast		Soft Tissue with Contrast						
Injector Coupling	Off		On						
Contrast Protocol			Thorax [factory]						
Absolute Start Delay	2		25						
Decision Tree									
wo CM									
Contrast medium?									
yes	no								
w CM	wo CM								
Thorax FoV [factory]									
Scan Strategy									
	standard acquisition		fast acquisition						
InfoLevel	Information		None						
InfoText	Please check scan time!								
Scan Mode	RoutineSpiralAdultThorax		RoutineSpiralAdultThorax						
Pitch	0.8		1.5						
Decision Tree									
standard acquisition									
FoV?									
yes	no								
fast acquisition	standard acquisition								

Thorax and Abdomen FoV [factory]			
Scan Strategy			
	standard acquisition	fast acquisition	
InfoLevel	Information	None	
InfoText	Please check scan time!		
Scan Mode	RoutineSpiralAdultAbdomenFlexdose	RoutineSpiralAdultAbdomenFlexdose	
Pitch	0.8	1.2	
Decision Tree			
standard acquisition			
FoV?			
yes	no		
fast acquisition	standard acquisition		
Thorax and Abdomen contrast medium necessary [factory]			
Scan Strategy			
	wo CM	w CM	
InfoLevel	None	None	
InfoText			
Scan Mode	RoutineSpiralAdultAbdomenFlexdose	RoutineSpiralAdultAbdomenFlexdose	
CARE Dose4D	On	On	
FAST kV optimized for	Non-Contrast	Soft Tissue with Contrast	
Injector Coupling	Off	On	
Contrast Protocol		Abdomen [factory]	
Absolute Start Delay	2	20	
Decision Tree			
wo CM			
Contrast medium?			
yes	no		
w CM	wo CM		
Abdomen contrast medium necessary [factory]			
Scan Strategy			
	wo CM	w CM	
InfoLevel	None	None	
InfoText			
Scan Mode	RoutineSpiralAdultAbdomen	RoutineSpiralAdultAbdomen	
CARE Dose4D	On	On	
FAST kV optimized for	Non-Contrast	Soft Tissue with Contrast	
Injector Coupling	Off	On	
Contrast Protocol		Abdomen [factory]	
Absolute Start Delay	2	75	
Decision Tree			
wo CM			
Contrast medium?			
yes	no		
w CM	wo CM		
Abdomen FoV [factory]			
Scan Strategy			
	standard acquisition	fast acquisition	
InfoLevel	Information	None	
InfoText	Please check scan time!		
Scan Mode	RoutineSpiralAdultAbdomen	RoutineSpiralAdultAbdomen	
Pitch	0.8	1.2	
Decision Tree			
standard acquisition			
FoV?			
yes	no		
fast acquisition	standard acquisition		

Kidney TwinSpiral DE possible [factory]			
Scan Strategy			
	TwinSpiral DE	Single Energy	
InfoLevel	None	Important Information	
InfoText		No Dual Energy information. Please check patient instruction!	
Scan Mode	RoutineSpiralAdultAbdomenTwinspiral	RoutineSpiralAdultAbdomenSn	
Decision Tree			
	TwinSpiral DE		
Dual Energy Recommendation			
Dual Energy	Single Energy		
TwinSpiral DE	Single Energy		
Whole Body soft tissue or bone [factory]			
Scan Strategy			
	standard dose	low dose	
InfoLevel	None	None	
InfoText			
Scan Mode	RoutineSpiralAdultAbdomen	RoutineSpiralAdultAbdomenSn	
CARE Dose4D	On	On	
FAST kV optimized for	Non-Contrast	Bone/Calcium	
FAST kV IQ level	145	72	
Decision Tree			
	standard dose		
Soft tissue/bone?			
bone	soft tissue		
low dose	standard dose		
Pelvis soft tissue or bone [factory]			
Scan Strategy			
	standard dose	low dose	
InfoLevel	None	None	
InfoText			
Scan Mode	RoutineSpiralAdultPelvis	RoutineSpiralAdultPelvisSn	
CARE Dose4D	On	On	
FAST kV optimized for	Non-Contrast	Bone/Calcium	
FAST kV IQ level	170	85	
Decision Tree			
	standard dose		
Soft tissue/bone?			
bone	soft tissue		
low dose	standard dose		
Pelvis contrast medium necessary [factory]			
Scan Strategy			
	wo CM	w CM	
InfoLevel	None	None	
InfoText			
Scan Mode	RoutineSpiralAdultPelvis	RoutineSpiralAdultPelvis	
CARE Dose4D	On	On	
FAST kV optimized for	Non-Contrast	Soft Tissue with Contrast	
Injector Coupling	Off	On	
Contrast Protocol		Abdomen [factory]	
Absolute Start Delay	2	85	
Decision Tree			
	wo CM		
Contrast medium?			
yes	no		
w CM	wo CM		
Pelvis TwinSpiral DE possible [factory]			
Scan Strategy			
	TwinSpiral DE	Single Energy	
InfoLevel	None	Important Information	
InfoText		No Dual Energy information. Please check patient instruction!	
Scan Mode	RoutineSpiralAdultPelvisTwinspiral	RoutineSpiralAdultPelvis	
Decision Tree			
	TwinSpiral DE		
Dual Energy Recommendation			
Dual Energy	Single Energy		
TwinSpiral DE	Single Energy		

C-Spine disc or bone [factory]						
Scan Strategy						
	standard dose		low dose			
InfoLevel	None		None			
InfoText						
Scan Mode	RoutineSpiralAdultNeck		RoutineSpiralAdultNeckSn			
CARE Dose4D	On		On			
FAST kV optimized for	Non-Contrast		Bone/Calcium			
FAST kV IQ level	195		97			
Decision Tree						
	standard dose					
	Spine disc/fracture?					
	spine disc	fracture				
	standard dose	low dose				
T- and L-Spine disc or bone [factory]						
Scan Strategy						
	standard dose		low dose			
InfoLevel	None		None			
InfoText						
Scan Mode	RoutineSpiralAdultSpine		RoutineSpiralAdultSpineSn			
CARE Dose4D	On		On			
FAST kV optimized for	Non-Contrast		Bone/Calcium			
FAST kV IQ level	285		142			
Decision Tree						
	standard dose					
	Spine disc/fracture?					
	spine disc	fracture				
	standard dose	low dose				
Spine TwinSpiral DE possible [factory]						
Scan Strategy						
	TwinSpiral DE		Single Energy			
InfoLevel	None		Important Information			
InfoText			No Dual Energy information. Please check patient instruction!			
Scan Mode	RoutineSpiralAdultSpineTwinspiral		RoutineSpiralAdultSpine			
Decision Tree						
	TwinSpiral DE					
	Dual Energy Recommendation					
	Dual Energy	Single Energy				
	TwinSpiral DE	Single Energy				
Spine TwinSpiral metal DE possible [factory]						
Scan Strategy						
	TwinSpiral DE		Single Energy			
InfoLevel	None		Important Information			
InfoText			No Dual Energy information. Please			
Scan Mode	RoutineSpiralAdultSpineTwinspiralMetal		RoutineSpiralAdultSpine			
Decision Tree						
	TwinSpiral DE					
	Dual Energy Recommendation					
	Dual Energy	Single Energy				
	TwinSpiral DE	Single Energy				
Abdomen CTA contrast medium sensitive [factory]						
Scan Strategy						
	standard iodine		less iodine		no iodine	
InfoLevel	None		Important Information		Important Information	
InfoText			Decreased eGFR!		Severely decreased eGFR!	
Scan Mode	RoutineSpiralAdultAngioBody		RoutineSpiralAdultAngioBody		RoutineSpiralAdultAngioBody	
CARE Dose4D	On		On		On	
FAST kV optimized for	Vascular		Soft Tissue with Contrast		Non-Contrast	
Decision Tree						
	standard iodine					
	(e)GFR					
	Normal	Decreased	Severely decreased			
	standard iodine	less iodine	no iodine			

Coronary Calcium Scoring HR variability [factory]				
Scan Strategy				
	standard	end systolic arrhythmia		
InfoLevel	None	None		
InfoText				
Scan Mode	CardiacSpiralAdult	CardiacSpiralAdult		
Phase Start	75	350		
Phase End	75	350		
Phase Unit	%	ms		
Decision Tree				
standard				
Heart Rate Variability				
Regular	Unstable			Arrhythmic
standard	standard			end systolic arrhythmia
RTP Brain [factory]				
Scan Strategy				
	standard	stereotactic		
InfoLevel	None	None		
InfoText				
Decision Tree				
standard				
Stereotaxic?				
yes	no			
stereotactic	standard			
Test Run [factory]				
Scan Strategy				
	Information for Test Run			
InfoLevel	Important Information			
InfoText	Due to a high table speed a test run is recommended.			
Decision Tree				
Information for Test Run				
TwinSpiral contrast medium [factory]				
Scan Strategy				
	Information for TwinSpiral			
InfoLevel	Important Information			
InfoText	To ensure constant contrast distribution, the delay needs to be larger than 75 seconds.			
Decision Tree				
Information for TwinSpiral				

Info template [factory]			
Scan Strategy			
	template 1		
InfoLevel	None		
InfoText			
Decision Tree			
template 1			
Extremity DE template [factory]			
Scan Strategy			
	template 1		
InfoLevel	None		
InfoText			
Scan Mode	RoutineSpiralAdultExtremities		
Decision Tree			
template 1			
Brain CTA template [factory]			
Scan Strategy			
	template 1		
InfoLevel	None		
InfoText			
Scan Mode	RoutineSpiralAdultAngioHead		
Decision Tree			
template 1			
Neck CTA template [factory]			
Scan Strategy			
	template 1		
InfoLevel	None		
InfoText			
Scan Mode	RoutineSpiralAdultAngioNeck		
Decision Tree			
template 1			
Body CTA template [factory]			
Scan Strategy			
	template 1		
InfoLevel	None		
InfoText			
Scan Mode	RoutineSpiralAdultAngioBody		
Decision Tree			
template 1			
Extremity CTA template [factory]			
Scan Strategy			
	template 1		
InfoLevel	None		
InfoText			
Scan Mode	RoutineSpiralAdultAngioExtremities		
Decision Tree			
template 1			

Contrast Protocol Name	CM Name	Iodine Concentration [mg/ml]	Phase Type	Flow [ml/s]	Volume [ml]	Duration [s]	CM Ratio [%]
Factory Contrast Protocols							
Brain [factory]		300					
Contrast			Contrast	2.5	60	24	100
Saline			Saline	2.5	40	16	0
Brain Perfusion [factory]		370					
Contrast			Contrast	6.0	41	7	100
Saline			Saline	6.0	40	7	0
Neck [factory]		300					
Contrast			Contrast	2.5	70	28	100
Saline			Saline	2.5	40	16	0
Thorax [factory]		300					
Contrast			Contrast	3.0	70	24	100
Saline			Saline	3.0	40	14	0
Abdomen [factory]		300					
Contrast			Contrast	4.0	100	25	100
Saline			Saline	4.0	40	10	0
CTA [factory]		300					
Contrast			Contrast	4.0	80	20	100
Saline			Saline	4.0	40	10	0
Brain and Neck CTA [factory]		370					
Contrast			Contrast	5.0	80	16	100
Saline			Saline	5.0	40	8	0
Pulmonary CTA [factory]		370					
Contrast			Contrast	5.0	65	13	100
Saline			Saline	5.0	40	8	0



Manufacturer's note:

The CE marking applies only to Medical Devices which have been put on the market according to the EC Regulation.

Unauthorized changes to this product are not covered by the CE mark and the related Declaration of Conformity.

Caution: Federal law restricts this device to sale by or on the order of a physician, dentist, or veterinarian (21 CFR 801.109(b)(1)).

The original language of this document is English.

Siemens Healthineers Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens-healthineers.com

Manufacturer

Siemens Shanghai Medical Equipment Ltd.
278 Zhou Zhu Road
201318 Shanghai
P.R. China

EU Authorized Representative

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany