
CRYOABLATION CASES: TANK ORDERING INFORMATION

To ensure smooth preparation for cryoablation procedures, please follow the guidelines below when ordering gas tanks:

TANK SPECIFICATIONS

- **Type:** Ultra High Purity Argon and Helium
 - **Pressure:** 6000 psi
 - **Connector:** Male adapters
 - **Quantity:** 2 tanks per case
-

ORDERING INSTRUCTIONS

- **Contact:** Josh Masalko to place tank orders
 - **Timing:** Order and **pick up tanks 1–2 days before each case**
-

CRYO EQUIPMENT


- **Handled by:** Fortec Medical
 - **What's Included:** All cryo equipment, **including regulators, except the gas tanks**
-

FORTEC CONTACT

Matt Warr

Territory Manager

 9407 Tradeport Drive, Orlando, FL 32827

 800-963-7101

 mwarr@fortecmedical.com