

# Implant Integrity

Updated

11/18/2023

Indications: implant rupture, implant deformity, assess implant deformity.

IV contrast is not given for this protocol.

<b>Pulse Sequence</b>	<b>PACS Name</b>	<b>plane</b>	<b>fat sat</b>	<b>slice (mm)</b>	<b>gap (mm)</b>	<b>first slice</b>
STIR	STIR AX BILAT	ax	no	4	1	top
T1	T1 AX BILAT	ax	no	5	5	top
STIR	STIR RIGHT	sag	no	4	1	right
T1	T1 RIGHT	sag	no	5	5	right
STIR	STIR LEFT	sag	no	4	1	right
T1	T1 LEFT	sag	no	5	5	right