

# **MRA Neck (Carotids)**

Updated 11/15/22

Reviewed 05/14/25

## **WITHOUT ONLY:**

Pulse Sequence	PACS Name	plane	fat sat	slice (mm)	gap (mm)	first slice
2D TOF	2D TOF AX	ax	no	3	0.9	base

## **RECONS:**

horizontal MIP spinners of both carotids together and each separately

## **WITHOUT & WITH:**

Pulse Sequence	PACS Name	plane	fat sat	slice (mm)	gap (mm)	first slice
2D TOF	2D TOF AX	ax	no	3	0.9	base
3D T1	3D T1 PRE COR	cor	no	1	0.2	back

**CONTRAST** - 2 mL/sec standard dose gadolinium (0.2 mL/kg Clariscan or 0.1 mL/kg Gadavist) followed by 20 mL saline flush.

3D T1	3D T1 PRE COR	cor	no	1	0.2	back
-------	---------------	-----	----	---	-----	------

## **RECONS:**

coronal subtractions

horizontal MIP spinners of both carotids together on the TOF sequence

horizontal MIP spinners of both carotids together and each separately on the post contrast sequence