

# Neck Parotids

Updated

1/15/2021

Remind patient not to swallow during image acquisition to reduce motion artifact.

Go to MRIMaster.com for a guide of proper positioning.

<b>Pulse Sequence</b>	<b>PACS Name</b>	<b>plane</b>	<b>fat sat</b>	<b>slice (mm)</b>	<b>gap (mm)</b>	<b>first slice</b>
T1	T1 AX	ax	no	3	1.2	base
T2	T2 AX	ax	no	3	1.2	base
T2	T2 FS AX	ax	yes	3	1.2	base
T1	T1 COR	cor	no	3	1.2	ant
T2	T2 COR	cor	no	3	1.2	ant
T2	T2 FS COR	cor	yes	3	1.2	ant

**CONTRAST** - 2 mL/sec standard dose gadolinium (0.2 mL/kg Clariscan or 0.1 mL/kg Gadavist) followed by 20 mL saline flush.

T1	T1 POST AX	ax	yes	3	1.2	base
T1	T1 POST COR	cor	yes	3	1.2	ant
T1	T1 POST SAG	sag	yes	3	1.2	left