

# **Bony Pelvis Mass**

Updated

6/15/2017

Use small field-of-view (FOV) centered on the area of interest with borders marked with skin capsule.

Angle imaging planes as per routine protocols.

Go to MRIMaster.com for a guide of proper positioning.

<b>Pulse Sequence</b>	<b>PACS Name</b>	<b>plane</b>	<b>fat sat</b>	<b>slice (mm)</b>	<b>gap (mm)</b>	<b>first slice</b>
T1	T1 AX	ax	no	5	1	top
T1	T1 FS AX	ax	yes	5	1	top
T2	T2 AX	ax	no	5	1	top
T2	T2 FS AX	ax	yes	5	1	top
T1	T1 COR	cor	no	5	1	front
T2	T2 FS COR	cor	yes	6	1	front

**CONTRAST** - 2 mL/sec standard dose gadolinium (0.2 mL/kg Clariscan or 0.1 mL/kg Gadavist) followed by 20 mL saline flush.

T1	T1 FS POST AX	ax	yes	5	1	top
T1	T1 FS POST COR	cor	yes	5	1	front

## **RECONS:**

axial subtractions