
ADULT BONY PELVIS

- **Pelvis**

- Order - XR Pelvis Routine
 - View
 - #1 AP crosswire
-

- **Sacrum/Coccyx**

- Order - XR Sacrum/Coccyx
 - Views
 - #1 AP axial sacrum 15° cephalic angulation
 - #2 AP axial coccyx 10° cephalic angulation
 - #3 Lateral sacrum & coccyx
-

- **Sacroiliac (SI) Joints**

- Order - XR Sacroiliac Joint
 - Views
 - #1 AP axial 30° cephalic angulation
 - #2 RPO 25-30° angulation
 - #3 LPO 25-30° angulation
-