

# Temporomandibular Joints (TMJs)

Updated

11/22/2023

IV contrast is not given for this protocol.

Go to MRIMaster.com for a guide of proper positioning.

<b>Pulse Sequence</b>	<b>PACS Name</b>	<b>plane</b>	<b>fat sat</b>	<b>slice (mm)</b>	<b>gap (mm)</b>	<b>first slice</b>
T1	T1 SAG CLOSED	sag	no	3	0.3	left
T1	T1 SAG OPEN	sag	no	3	0.3	left
T2	T2 FS COR CLOSED	cor	yes	3	0.3	ant
T2	T2 FS AX CLOSED	cor	yes	3	0.3	base
Gradient Cine	GRE SAG RIGHT	sag	no	4	0.3	left
Gradient Cine	GRE SAG LEFT	sag	no	4	0.3	left