
OVERVIEW

- **Primary Synonym:** IR Tunneled Drain Placement
 - > The actual order name viewable in Cerner & Department Order Entry.
 - > This is a new order that will need to be created with details provided below.
- The order requisition prints to **IR Printers**.
 - > Riverside - xmlp_sv rad ip ir re
 - > Southside - sl1040_ir_req
 - > Clay - cc1059_svc_ir_req
 - > St Johns - sj_cathlab_tech_req

ORDER SPECIFICS

- **Order Name:** IR Tunneled Drain Placement
 - CPT - none
 - CDM - 7036
- **Procedure:** (required drop down list)
 - > Right chest/pleural
 - > Left chest/pleural
 - > Larger side chest/pleural
 - > Bilateral chest/pleural
 - > Abdominal/peritoneal
- **Reason for Exam:** (required drop down list)
 - > Pleural effusion
 - > Ascites
- **Fluid Labs:** (required drop down list)
 - > Without fluid labs (therapeutic)
 - > With fluid labs (diagnostic & therapeutic)
- **Reference:**
 - > Ordering provider must order any needed fluid labs. Search for "Pleural Fluid Labs" Power Plan.
- **Priority:** (required drop down list)
 - > Routine (default selection)
 - > Expedite
 - > STAT
- **Transport Mode:** (required drop down list)
 - > Stretcher (default selection)
- **Reference:**
 - > Visit ir.svrads.com for guidelines.